HumanTech

Spine

TRISTAN®

CERVICAL INTERBODY FUSION

0297

HumanTech Spine GmbH

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D-71144 Steinbronn

Tel: +49 (0) 7157/5246-71

Fax: +49 (0) 7157/5246-33

e-mail: info@humantech-spine.de

www.humantech-spine.de

Stand: 25.10.2019

DE

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EN

IMPORTANT NOTICE

This instruction for use must be read carefully and its content must be adhered.

SCOPE

The instructions for use is valid for sterile delivered implants:

- TRISTAN® PEEK...
- TRISTAN® PEEK...-S
- TRISTAN® R-PEEK-Ti...
- TRISTANflex...
- TRISTAN® Ti... sterile,
- and unsterile delivered implants:
- TRISTAN® Ti...

BASIC DESIGN

The TRISTAN® Cervical Interbody Fusion, used as an implant system for the fusion of vertebral bodies in spinal surgery, is comprised of intervertebral implants in different sizes, angles and heights, enabling adaptations to be made while taking into account the unique pathology of the individual patient.

MATERIAL

There are different implant product lines.

Basic material PEEK Optima®:

TRISTAN® PEEK...

The base body of products in the TRISTAN® PEEK range is exclusively made from PEEK Optima®. To achieve a better radiographic contrast, radiographic markers made out of titaniumalloy Ti6Al4V respectively tantalum have been incorporated into the implant.

TRISTAN® R-PEEK-Ti...

The base body of products in the TRISTAN® R-PEEK-Ti range also consists of PEEK Optima®. In addition, the surfaces of the implant has a titanium coating. To achieve a better radiographic contrast, radiographic markers made out of titaniumalloy Ti6Al4V respectively tantalum have been incorporated into the implant. In addition, there are other components made out of titanium alloy Ti6Al4V included (e.g. spikes).

TRISTANflex and TRISTAN® PEEK...-S

The base body of products in the TRISTANflex and TRISTAN® PEEK...-S ranges also consists of PEEK Optima®. To achieve a better radiographic contrast, radiographic markers made out of titaniumalloy Ti6Al4V respectively tantalum have been incorporated into the implant. In addition, there are other components made out of titanium alloy Ti6Al4V included (e.g. spikes).

Basic material Titanium alloy:

TRISTAN® Ti and TRISTAN® Ti... sterile

Products of the TRISTAN® Ti and TRISTAN® Ti...sterile range consist of the titanium alloy Ti6Al4V.

INDICATIONS FOR USE

TRISTAN® is intended for the operative treatment of diseases and injuries of the cervical spine in patients whose overall skeletal growth is fully developed, in particular for such indications as:

- disk herniation
- hard disk herniation
- mechanical instabilities
- calcification of the posterior longitudinal ligament
- osteochondrosis (for TRISTANflex reduced to Modic I and II)
- ppinal canal stenosis

AREA OF APPLICATION

The area of application is the cervical spine from C3 to C7.

GENERAL CONDITIONS FOR USE

- Only surgeons with the required professional training in the field of spinal surgery may insert the implants. The use of the implants must be decided in keeping with the surgical and medical indications, the potential risks and the limitations relating to this type of surgical intervention, as well as with the knowledge of

contraindications, side effects and known precautionary measures and the knowledge of the condition and the physical, chemical and biological properties of the implant.

- The activities and physical activity of the patient have a significant influence on the implant's useful life. Patients must be informed, that every activity increases the risk of loss, loosening, dislocation, migration, bending or breakage of implant or implant components. Informing patients about limitations to their activities in the postoperative phase and postoperatively monitoring patients are crucial factors in assessing the development of the surgical result and the condition of the implant. The above-mentioned effects may occur even if the implant is well integrated and the activity restrictions are complied with. The patient must be informed of this. If such an effect occurs, the surgeon must decide whether a revision of the implant or other measures should be taken, taking into account the well-being of the patient and the risks involved.
- This instruction for use provides essential information, but is not sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a clearly defined post-operative treatment regimen, and to conduct scheduled follow-up examinations.
- It is recommended not to replace components of TRISTAN® system with components or products of other systems from another source or another manufacturer, or make a direct connection to components of other systems. HumanTech Spine GmbH is not liable if this recommendation is not followed.
- Following the notes in the operating instructions (surgical technique) is essential. You will find these in the Internet at www.humantech-spine.de or you can obtain them directly from your HumanTech representative. It is advisable to use the instruments described there and provided by the manufacturer. There may be complications if the implant is used with or without the dedicated instrument. For the use of foreign instruments, no liability can be accepted.
- Proceed with extreme caution in the region of the spinal cord and the roots of the nerves, since damage to the nerves can lead to the impairment of neurological functions.
- Breakage, slippage or incorrect use of the instruments or implants can injure the patient or the operating staff.
- Do not use bone cement, as this material makes the removal of the components difficult or impossible. The heat produced by the hardening process can damage or deform the PEEK implants.
- Residues consisting of implant material and/ or not from implant material should be removed.
- Damaged implants must not be implanted.
- Never reuse the implants. Even if the implant appears to be intact after revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break.
- Never use the implants, if they had contact with body fluids or third-party tissue.
- Handle removed implants in such a way that their reuse is not possible.
- Safety and compatibility of the device in the setting of magnetic resonance (imaging) have not been evaluated. No thermal test or test of migration has been performed on the device in this setting.
- To achieve a good primary fixation, the TRISTANflex and TRISTAN PEEK...-S Cages have spiky and sharp spikes, therefore they must be carefully handled and a sufficient distraction is necessary.
- It is recommended to use an additional anterior plating if the bone quality or the implant site can't be evaluated pre- or interoperatively in question of primary fixation or if the implant site or the bone quality is expected not to ensure a proper primary fixation.
- When using PEEK implants particular care in preparing the intervertebral disc space is needed to keep the placement forces low
- For a better ongrowth, the R-PEEK-Ti implants are coated with a titanium surface. Take care during insertion of these cages as the surface is structured and rough.

CONTRAINDICATIONS

Contraindications may be either relative or absolute. The selection of a particular implant must be weighed carefully against the overall assessment of the patient. The following conditions can have an adverse impact on the chances of successful surgery:

- all destructive and inflammatory diseases of the cervical vertebral body and vertebral motor segments such as tumors, spondylitis and fractures with a high degree of instability, with pronounced static and structural alterations of the vertebral body and segment structure, requiring the replacement of the vertebral body and other stabilizing measures for the reconstruction and stabilization of the segment
- acute infections or significant risks of infections (weakened immune system)
- indications of local inflammation
- fever or leucocytosis
- pathological obesity
- pregnancy
- psychological disorders
- pronounced anatomical deformities resulting from congenital abnormalities
- any other medical or surgical condition hindering possible improvement through the use of the implant, such as the existence of congenital abnormalities, fracture

in the vicinity of the operating point, elevated deposition rate not indicated by other clinical pictures, elevated leucocyte count (white blood cells), or distinct leftward shift in the white blood cell differential blood count.

- Diseases of the joints, bone resorption, osteopenia, osteomalacia and/or osteoporosis are relative contraindications, since these can limit the degree of correction and stabilization that can be achieved.
- Any neuromuscular disease that would place excess strain on the implant during the healing period.
- known hereditary or acquired brittleness of the bones or calcification problems
- for TRISTANflex: Osteochondrosis Modic III
- spondylolisthesis that cannot be reduced to the first degree
- Suspicion of a metal allergy or intolerance, and documented metal allergy or intolerance. Appropriate tests should be carried out.
- all cases in which the use of components of different metals or alloys is necessary
- already existing fusion in the segments treated
- all cases in which fusion is not necessary
- all cases in which the implant component selected for use is too large or too small to achieve a satisfactory result
- all patients with inadequate tissue structure on the operative side or inadequate bone or bone quality
- all patients in which the use of the implant would interfere with anatomical structures or restrict physiological performance
- all patients not willing to follow postoperative instructions
- all cases not described in the indications

SIDE EFFECTS / COMPLICATIONS

The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

- loosening of the implant
- shifting of the implant
- breakage of the implant, respectively detached parts of the implant and wear
- foreign body reactions on the implant, including possible tumor formation, autoimmune disease and/or scar formation
- scar formation with possible neurological impairments, nerve compression or pain
- neurological complications
- paralysis
- soft tissue lesions
- neurological or spinal lesions of the dura mater as a result of surgical trauma
- The migration of the implant in the posterior direction can cause a loss of neurological function, the occurrence of radiculopathy, dural tears and/or pain.
- Neurovascular impairment, including paralysis
- temporary or permanent retrograde ejaculation in men, or other severe injuries
- Discharge of cerebrospinal fluid
- inability to perform daily tasks
- urinary retention or loss of control over the bladder, or other types of urinary system impairment
- difficulty swallowing (dysphagia)
- alteration of the curvature and stiffness of the spine
- Failure to achieve the desired result of surgery, malpositioning of the implant
- partial loss of the degree of correction achieved during the operation
- superficial or deep-seated infections and inflammations, such as discitis or arachnoiditis
- bone fracture above or below the segments treated
- late osseointegration or no visible fusion and pseudarthrosis
- fracture, micro fracture, resorption, damage to or penetration of vertebral bone and/or the bone transplant and/or the bone transplant donor site above and below the segment treated
- adjacent segment disease
- bleeding and/or hematomas
- deep-seated venous thrombosis, thrombophlebitis, pulmonary embolism
- complications at the bone transplant donor site
- herniated disk, disk destruction or disk degeneration at, above or below the segment treated
- loss or enhancement of spinal mobility or functions
- damage to the reproductive system, sterility and sexual functional disturbances
- development of respiratory problems, for example pulmonary embolism, atelectasis, bronchitis, or pneumonia
- alteration of the mental state
- alteration in the growth of the fused spine
- allergic reaction to the materials used
- decrease in bone density
- bursitis
- death

PACKAGING, LABELING AND STORAGE

- Handle and store the implant components carefully. Damage to the implant can significantly reduce the stability and long-term stability of the implant system and can cause cracking and/or higher internal stresses, possibly resulting in the breakage of the implant.
- Storage of the implants and instruments should be at room temperature. Ambient influences, such as salty air, humidity, and chemicals, must not be allowed to act on the implants.

- Thorough inspection is recommended before operating in order to ensure that the instruments or implants have not been damaged during storage or previous procedures.

The implants are delivered in both NON STERILE and STERILE packaged condition:

Unsterile delivered implants:

UNSTERILE delivered implants are labeled as UNSTERILE and must be cleaned and sterilised before use (see REPROCESSING and STERILISATION).

The implants are supplied as a complete set or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilisation in the original pouch packaging is not allowed.

Sterile delivered implants:

STERILE delivered implants are sterilized by a validated ethylene oxide or gamma sterilization procedure and labeled as STERILE. Cleaning, processing and sterilization prior to use need not be performed.

The implants are supplied individually packaged in a blister-pouch- or blister-blister-combination and protected by a covering box. The implants may be used only when the label on the outer packaging and also the inner packaging are intact. If the packaging is damaged or already open, the sterility of the implant is not guaranteed and the implant may not be used.

The implants may not be used when the shelf life indicated has been exceeded.

Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or with damaged sterile packaging is not intended.

REPROCESSING

Unsterile delivered products of the TRISTAN®-System must be processed prior to sterilisation. The necessary steps are described in the reprocessing instructions for Ti implants.

STERILISATION

Unsterile delivered products of the TRISTAN-System must be sterilised or re-sterilised prior to use by means of validated steam sterilisation in an autoclave. The minimum recommended sterilisation parameters are as follows:

Steam	prevakuum	132°C - 134°C	5 minutes
Steam	prevakuum	121°C	30 minutes

PRODUCT COMPLAINTS

Any persons active in health services (for example a customer or user of this product system) who wishes to make complaints of any kind or is not satisfied with the handling of the product with regard to quality, identity, durability, stability, safety, effectiveness, and/or function should inform their HumanTech representative.

In the event that a TRISTAN® implant exhibits a "functional defect" (i.e. it does not fulfill the performance specifications or does not function as intended) or this is suspected, immediately inform your HumanTech representative.

In the event that a HumanTech product exhibits a functional defect, contributing to the death or serious injury of the patient, immediately inform your representative by telephone, fax or in writing.

Should you wish to make a complaint, we ask that you inform us in writing, specifying the name, item number and batch number of the components and your name and address, together with as detailed as possible a description of the fault.

OTHER INFORMATION

In the event of complaints, suggestions or references to the contents of these instructions or the use of the product, please contact us at the address above.


© HUMANTECH Spine GmbH. All rights reserved.

	Manufacturer
	Do not re-use
	Batch number
	Reference number
	Read the documentation
	Keep dry
	Unsterile
	sterilized by irradiation
	Sterilisation with ethylene oxide

	Do not re-sterilize
	Date of expiry
	Do not use when packaging is damaged
	Attention

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