ADONIS®

LUMBAR INTERBODY FUSION

Gebrauchsanweisung	DE	2 - 6
Instruction for use	EN	7 - 10
Manual de instrucciones	ES	11 - 15
Instructions d'utilisation	FR	16 - 20
Istruzioni per l'uso	IT	21 - 25
Instruções de utilização	PT	26 - 30
Инструкция по применению	RU	31 - 35
Gebruiksaanwijzing	NL	36 - 40
Návod k použití	CZ	41 - 45
Lietošanas instrukcija	LV	46 - 49
Οδηγίες χρήσης	EL	50 - 54
Inštrukcie na používanie	SK	55 - 59
Revision		60

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IMPORTANT NOTICE

This instruction for use must be read carefully and its content must be adhered to.

SCOPE

The instructions for use is valid for sterile delivered implants:

- ADONIS® PEEK ...
- ADONIS® R-PEEK-Ti...
- ADONIS® Ti... sterile
- ADONIS® LLIF Ti 3D ...

and non sterile delivered implants:

- ADONIS® Ti...

INTENDED USE

The implants of the ADONIS® Lumbar Interbody Fusion System are used for long-term use to stabilize and support the fusion of lumbosacral vertebral bodies in spinal surgery during the normal healing process after surgical correction of diseases of the spine. The implants are used to replace the intervertebral disc to restore the original disc height, as a support element between two or more lumbar or lumbosacral vertebral bodies, and to correct the profile of the spine.

BASIC STRUCTURE

Different types of implants are available for the various surgical approaches to the spine:

posterior approaches:

- ADONIS® PLIF – posterior lumbar interbody fusion
- ADONIS® TLIF – transforaminal lumbar interbody fusion
- ADONIS® UniLIF – unilateral lumbar interbody fusion

anterior approach:

- ADONIS® ALIF – anterior lumbar interbody fusion

Lateral approach:

- ADONIS® LLIF – lateral lumbar interbody fusion

The implants are available in different dimensions and angulations, so that adaptations taking into account the unique anatomy of individual patients are possible.

MATERIAL

There are different implant product lines.

Basic material PEEK Optima®:

ADONIS® PEEK

For products of the ADONIS® PEEK line, the base body is entirely made of PEEK Optima®. To achieve a better radiographic contrast, radiographic markers made out of titanium alloy Ti6Al4V respectively tantalum have been incorporated into the implant.

ADONIS® R-PEEK-Ti

For products of the ADONIS® R-PEEK-Ti line, the base body is also made of PEEK Optima®. In addition, the surfaces of the implant have a titanium coating. To

achieve a better radiographic contrast, radiographic markers made out of titanium alloy Ti6Al4V respectively tantalum have been incorporated into the implant.

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Products of the ADONIS® Ti line consist of the titanium alloy Ti6Al4V.

INDICATIONS

ADONIS® is intended for the surgical treatment of diseases and injuries of the lumbar and lumbosacral spine in patients whose general skeletal growth is completed, especially for indications such as:

- Prolapsed intervertebral disc
- Hard prolapsed intervertebral disc
- Mechanical instabilities
- Calcification of the posterior longitudinal ligament
- Osteochondrosis
- Lumbal spinal stenosis

The ADONIS® Lumbar Interbody Fusion System is intended for use with an additional dorsal fixation system, e.g. VENUS®. Implants of the ADONIS® system cannot be used for stand-alone use.

Implants of the ADONIS® Lumbar Interbody Fusion System can be used with bone grafts or bone substitute material.

SUITABLE PERSONNEL

The implants may only be implanted by surgeons who have completed the necessary training in spinal surgery.

GENERAL CONDITIONS FOR USE

The use of the implants must be decided in accordance with the surgical and medical indications, the potential risks and restrictions associated with this type of surgical procedure, as well as in knowledge of the contraindications, side effects and defined precautions in awareness of both the nature and the physical, chemical and biological properties of the implant.

- The activities / movement behaviour of the patient has a significant influence on the implant's useful life. Patients must be informed, that any activity increases the risk of loosening, dislocation, migration, bending or fracture of the implant or implant components. It is crucial to inform patients about limitations to their activities in the postoperative phase and to monitor patients postoperatively in order to assess the development of the fusion and the condition of the implant. Even when a solid bone fusion has occurred, implant components may still bend, break or loosen. Therefore, the patient must be informed that implant components can also bend, break or loosen if the restrictions in activities are properly followed.
- The information in this instruction for use is necessary but not sufficient for the use of the system. This information is not a substitute for the

following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for implantation; the surgeon's ability to ensure the patient's willingness to cooperate, to adhere to an appropriately defined post-operative treatment program, and to undergo scheduled follow-up examinations.

- In case of complications, the surgeon must decide whether a revision of the implant or other measures should be carried out, taking into account the patient's condition and the possible risks involved.
- Components of the ADONIS® system may not be replaced by components/ products from other systems, from another source or from a different manufacturer. Furthermore, no direct connection of components/ products of the systems to components of other systems may be established. If this recommendation is not complied with or if the products are otherwise used or used improperly, HumanTech Spine GmbH assumes no responsibility.
- The instruction in the surgical technique must be adhered to. You will find these in the Internet at www.humantech-spine.de or you can obtain them directly from your HumanTech representative. It is advisable to use the instruments described there and provided by the manufacturer. No liability can be assumed for the use of third-party instruments.
- Proceed with extreme care in the area of the spinal cord and the nerve roots, as any damage to nerves can result in loss of neurological functions.
- Breakage, slippage or incorrect use of the instruments or implants can result in injury to the patient or surgical personell.
- bone cement may not be used, as this material makes the removal of components difficult or even impossible. The heat produced by the hardening process can damage or deform the PEEK implants.
- Residues consisting of implant material and/ or not from implant material should be removed.
- Damaged implants must not be implanted.
- Under no circumstance may the implants be re-used. Even if the implant appears to be intact after revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break.
- Implants that have already come into contact with a patient's body fluids or tissues or have been soiled must not be reused.
- Once removed implants shall be treated in such a way that they cannot be reused.
- The safety and compatibility of the implants were not assessed in terms of magnetic resonance (imaging). No thermal or migration test has been performed

under these influences.

- When using X-ray equipment for diagnosis or control, the national limit values for the radiation dose must be observed.
- When using PEEK implants particular care in preparing the intervertebral disc space is needed to keep the placement forces low
- When inserting the PEEK-Ti coated and 3D-printed implants, special care is required, as the coating / implant geometry for the improvement of the ingrowth behavior, consists of a structured, rough surface.

CONTRAINDICATIONS

Contraindications may be either relative or absolute. The selection of a particular implant must be weighed carefully against the overall assessment of the patient. The following conditions can have an adverse impact on the chances of successful surgery:

- acute infections or significant risks of infections (weakened immune system)
- all destructive and inflammatory diseases of the vertebral body or movement segment such as tumours, spondylitis and fractures with a high degree of instability, with a strong static and structural change of the vertebral body and segment structure, which require vertebral body replacement and other stabilising measures for reconstruction and stabilisation of the segment.
- indications of local inflammation
- fever or leucocytosis
- morbid obesity
- pregnancy
- mental disease
- Drug or alcohol abuse
- Wound healing disorders
- Any neuromuscular disease that would place an exceptionally high load on the implant during the healing period
- Severe anatomical deformities caused by congenital abnormalities
- Extreme malalignments which impair the stability of the instrumentation
- Missing intact adjacent segments
- bone resorption, osteopenia, osteomalacia and/ or osteoporosis are relative contraindications, as these can limit the degree of achievable correction/ stabilization
- known hereditary or acquired bone brittleness or calcification problems
- Suspected allergy or intolerance, as well as documented allergy or intolerance, to the material used. Appropriate tests must be carried out
- all cases in which fusion is not required
- all cases where the implant components selected for use are too big or too small to achieve a satisfactory result
- all cases where the use of components of different metals or alloys is necessary.

- any patient with inadequate tissue structure at the surgical site or with inadequate bone stock or bone quality
- all patients in whom the use of the implant would interfere with anatomical structures or restrict physiological performance
- any patient who is unwilling to follow postoperative instructions
- any other medical or surgical condition that prevents possible improvement through the use of the implant, such as the presence of congenital abnormalities, fracture near the surgical site, increase in deposition rate not indicated by other clinical pictures, increase in leucocyte count (white blood cells), or marked left shift in the white blood cell differential blood count.
- spondylolisthesis that cannot be reduced to the first degree
- already existing fusion in the segments treated
- all cases not described in the indications

SIDE EFFECTS / COMPLICATIONS

- The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:
 - implant loosening
 - implant migration
 - implant subsidence
 - implant breakage, respectively detached parts of the implant
 - Signs of abrasion, for example occurrence of microparticles in the area of the implant
 - neurological complications
 - foreign body reactions to the implant, including possible tumor formation, autoimmune disease and/or scar formation
 - Neurological complications with paralysis, functional failures, loss of sensitivity and pain, e.g. due to scarring or durotomy, irritation of the nerve roots (radiculopathy) and peripheral nerves
 - Neurovascular impairment, including paralysis
 - Soft tissue complications like lesion or heterotopic ossification
 - Pressure on surrounding tissue, nerves or organs; injury of surrounding tissue, nerves or organs
 - neurological or spinal lesions of the dura mater as a result of surgical trauma with the risk of loss of cerebrospinal fluid (CSF) or a CSF fistula
 - Injury to the lymphatic vessels / lymph leakage
 - Failure to achieve the desired surgical result, implant malpositioning
 - fracture, micro fracture, resorption, damage to or penetration of vertebral bone and/or bone graft and/or the bone graft extraction site above or below the treated segment(s)
 - Negative impairment of the adjacent segments
 - late osseointegration or no visible fusion and pseudarthrosis
 - superficial or deep-seated infections and inflammations, such as discitis or arachnoiditis
- allergic reaction to the materials used
- bleeding and/or hematomas in particular of the vertebral arteries
- deep vein thrombosis, thrombophlebitis, pulmonary embolism
- complications at the bone transplant donor site
- herniated disk, disk destruction or disk degeneration at, above or below the treated segment(s)
- loss or enhancement of spinal mobility or functions
- Physiological restrictions, such as adjacent segment degeneration
- damage to the reproductive system, sterility and sexual dysfunction
- temporary or permanent retrograde ejaculation in men, or other serious injuries
- inability to perform daily tasks
- urinary retention or loss of bladder control, or other types of impairment of the urological system
- change of mental state
- changed growth of the fused spine
- partial loss of the degree of correction achieved during the operation
- changes in the curvature and stiffness of the spine
- Negative complications of bony structures, e.g. fracture of a vertebral body, necrosis, reduction in bone density
- bone fracture above or below the segments treated
- development of respiratory problems, for example pulmonary embolism, atelectasis, bronchitis, or pneumonia
- bursitis
- death
- Revision surgery
- As well as all general surgical risks, like high blood pressure, coagulation disorder (coagulopathy), blunt trauma, fever, pleural effusion and neck stiffness

PACKAGING, LABELING, TRANSPORT AND STORAGE

- The implant components must be handled, transported and stored with care. Damages to the product packaging or damages to the implant itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increased internal loads, which can result in a fracture of the implant.
- The implants and instruments should be stored at room temperature. Environmental influences such as salty air, humidity, chemicals etc. are not allowed to act on the implants.
- Before the surgery, a careful inspection of the components of the ADONIS® Lumbar Interbody Fusion System to be used must be carried out to rule out damage caused by storage, transport or previous procedures.
- The implants are delivered in both NON STERILE and STERILE packaged condition:

Non sterile delivered implants:

NON STERILE delivered implants are labeled as NON STERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILIZATION). The implants are supplied as an implant system in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

Sterile delivered implants:

STERILE delivered implants are sterilized by a validated gamma sterilization or ethylene oxide procedure and labeled as STERILE. Cleaning, disinfection and sterilization prior to use must not be performed. The implants are delivered individually packaged in a blister and protected by an outer carton. The implants may only be used if the label on the outer packaging and also the inner packaging are intact. If the packaging is damaged or already open, the sterility of the implant is not guaranteed and the implant must not be used. The implants must not be used if the indicated expiry date has been exceeded. Processing, reprocessing, sterilization or re-sterilization of the products after opening the sterile packaging or if the sterile packaging is damaged is not intended. HumanTech Spine GmbH assumes no responsibility for the use of re-sterilized implants regardless of the person who performed the re-sterilization or the method used. The rules of asepsis must be observed when removing the implant from the sterile packaging. The sterile packaging may only be opened immediately before the implant is inserted. It is recommended that a replacement implant is always kept available.

CLEANING, DISINFECTION AND STERILIZATION

NON STERILE packaged delivered products of the ADONIS® Lumbar Interbody Fusion System must be cleaned, disinfected and sterilized prior to use. All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the instruction “Processing unsterile delivered Spine”. The latest version can always be found on our homepage <https://www.humantech-spine.de/378-en-IFUs.html>, as well as on request from HumanTech Spine GmbH.

DISPOSAL

The products must be disposed of in accordance with local regulations and environmental legislation taking into account the level of contamination.

PRODUCT COMPLAINTS

Any healthcare professional (e.g. customer or user of this

product system) who has complaints of any kind or is dissatisfied using the product in terms of quality, identity, strength, durability, safety, effectiveness and/or function should notify the appropriate HumanTech representative. If an ADONIS® implant should ever exhibit a malfunction (i.e. does not meet the performance specifications or does not function as intended), or if the above is suspected, the HumanTech representative should be informed immediately.

If a HumanTech product should ever malfunction causing or contributing to the death or serious injury of a patient, the representative is to be informed immediately by telephone, fax or in written form.





If you have any complaints we ask you to state the name, article number and batch number of the components as well as your name and address along with a detailed description of the fault and forward this information to us in written form.

FURTHER INFORMATION

If you have any complaints, suggestions or other points about the contents of this instruction for use or the use of the product, please contact the above address.

© HUMANTECH Spine GmbH. All rights reserved.

SYMBOLS

	Manufacturer
	Do not re-use
	Batch number
	Reference number
	Read the documentation
	Keep dry
	Do not re-sterilize
	Do not use when packaging is damaged
	Date of expiry
	Sterilization with ethylene oxide
	sterilized by irradiation
	Non sterile
	Attention

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© HUMANTECH Spine GmbH. Všetky práva vyhradené.