SAMSON®

VERTEBRAL BODY REPLACEMENT

Gebrauchsanweisung	DE	2 - 5
Instruction for use	EN	6 - 9
Manual de instrucciones	ES	10 - 13
Instruções de utilização	PT	14 - 17
Návod k použití	CZ	18 - 21
Inštrukcie na používanie	SK	22 - 25
Οδηγίες χρήσης	EL	26 - 30
Revision		31

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IMPORTANT INFORMATION

This instruction for use must be read carefully and its content must be adhered to.

SCOPE

The instruction for use is applicable for all sterile and non sterile delivered implants of the SAMSON®-System.

INTENDED USE

The SAMSON® Vertebral Body Replacement System (VBR) is an implant system for long-term use for anterior stabilization of the upper thoracic to lower lumbar spine as a replacement for one or more vertebral bodies in patients in whom general skeletal growth has ceased.

BASIC STRUCTURE

The SAMSON® Vertebral Body Replacement System consists of implant base bodies and implant endplates in various dimensions. The implant base bodies are available in different heights for different defect heights. The implant endplates are available in different widths and angles for adaptation to different vertebral body structures and for mapping the existing or achievable lordotic or kyphotic curvature of the spine. The unique anatomy of the individual patient can be taken into account by the possibility of connecting the different implant components with each other via plate holding screws.

MATERIAL

All components of the SAMSON® vertebral body replacement system are made of titanium alloy Ti6Al4V according to ISO 5832-3.

INDICATIONS FOR USE

The SAMSON® system is intended for the surgical treatment of tumorous, inflammatory and traumatic diseases which lead to instabilities around the anterior support, or compression of the neural structures, or for diseases which require treatment of an infection.

The SAMSON® system has been designed for use with an additional dorsal fixation system. Implants of SAMSON® system are not intended for stand-alone use.

SAMSON® is intended for patients whose general skeletal growth is finished.

When using a SAMSON® implant with base bodies Size 5 or Size 6, or when using 12 ° angled endplates mounted at the caudal and cranial position on the implant base body, ventral support (such as the VENUSnano®-Screw Rod System) is strongly advised.

SUITABLE PERSONNEL

The implants may only be implanted by surgeons who

have completed the necessary training in spinal surgery.

GENERAL CONDITIONS/ INFORMATIONS FOR USE

- The insertion of implants must be carried out in line with the surgical and medical indications, the potential dangers, and with the restrictions associated with this type of surgery. It should also be performed with an awareness of the contraindications, side effects, defined precautionary measures, and with regard for the structure and the physical, chemical and biological characteristics of the implant.
- The activities / movement behaviour of the patient has a significant influence on the service life of the implant. The patient must be informed that any activity increases the risk of loss, bending or fracture of the implant components. It is crucial to inform the patient about restrictions in activities in the post-operative phase and to monitor the patient post-operatively in order to assess the development of the fusion and the condition of the implant. Even when a solid bone fusion has occurred, implant components may still bend, break or loosen. Therefore, the patient must be informed that implant components can also bend, break or loosen if the restrictions in activities are properly followed.
- The information in this instruction for use is necessary but not sufficient for the use of this system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a clearly defined post-operative treatment regimen, and to conduct scheduled follow-up examinations.
- In case of complications, the surgeon has to decide whether a revision of the implant should be carried out, taking into account the patient's condition and the possible risks involved.
- Components of the SAMSON® system may not be replaced by components / products from other systems from another source or from a different manufacturer. Furthermore, no direct connection of components / products of the system to components of other systems may be established. If this is not complied with or if the products are otherwise used or used improperly, HumanTech Spine GmbH assumes no responsibility.
- Never reuse the implants. Even if the implant appears to be intact following revision, alterations within

the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break. Handle removed implants in such a way that their reuse is not possible.

- Implants that have already come into contact with a patient's body fluids or tissues or have been soiled must not be reused.
- Handle removed implants in such a way that their reuse is not possible.
- Care should be taken to maintain the integrity of the endplates of the vertebral bodies. Damage to the endplates of the vertebral bodies or excessive partial abrasion of the endplates of the vertebral bodies can lead to sintering of the implant and to a loss of segmental stability.
- The endplates of the implant should be selected in a sufficiently large size in order to make use of the maximum contact surface of the vertebral body. This will ensure the greatest possible stability and counteract implant sintering. However, in order to avoid injury to the adjacent structures, the endplates of the implant should not protrude beyond the vertebral body.
- The implant must be implanted under prestress in order to avoid dislocation. Expansion of the implant must be performed under continuous x-ray control in order to avoid over-distraction, damage to the implant endplates or other implant components, e.g. of the posterior fixation.
- The instructions in the operation instructions (Surgical Technique) are to be observed. You can download these from the internet from www.humantech-spine.de or directly from your HumanTech representative. It is advised to use the instruments described in the surgical technique, which are intended to be used with the system. Complications may occur if the implant is inserted with or without the intended instruments. No liability can be accepted for the use of external instruments.
- Proceed with extreme caution in the region of the spinal cord and the roots of the nerves, since damage to the nerves can lead to the impairment of neurological functions.
- Breakage, slippage or incorrect use of the instruments or implants can injure the patient or the operating staff or result in more time being required for surgery.
- Bone cement must not be used, as the removal of the components becomes difficult or impossible.
- Residues consisting of implant material and/ or not from implant material should be removed.
- Damaged implants must not be implanted.
- Safety and compatibility of the device in the setting of magnetic resonance (imaging) have not been evaluated. No thermal test or test of migration has been performed on the device in this setting.
- The endplates of the implant must be assembled correctly to the preassembled implant base body.
- Careful inspection of the assembled SAMSON®

implants is recommended before and during surgery to ensure that the implants are undamaged during assembly or previous procedures and have been correctly assembled according to Surgical Technique.

CONTRAINDICATIONS

Contraindications may be either relative or absolute. The selection of a particular implant must be weighed carefully against the overall assessment of the patient. The following conditions can have an adverse impact on the chances of successful surgery:

- Extreme malalignments, impairing the stability of the instrumentation
- Suspicion of a metal allergy or intolerance, and documented metal allergy or intolerance. Appropriate tests should be carried out.
- All cases in which the selected implant components are too large or too small to achieve a satisfactory result
- Every patient in which the use of the implant would interfere with anatomical structures
- Any neuromuscular disease that would place excess strain on the implant during the healing period.
- Any patient who is not willing or able to follow post-operative instructions
- Morbid obesity
- Pregnancy
- Use in combination with implants made of different metals or alloys
- All cases which are not described in the indications
- Any other medical or surgical condition which prevents possible improvement through the use of the implant

SIDE EFFECTS / COMPLICATIONS

The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

- Missing intact connection segments, bone resorption, osteopenia, osteomalacia and/or osteoporosis, as well as diffuse spinal tumours and infections, can, in addition to other causes, lead to the reduction of instrumentation stability, limit the degree of obtainable correction and stabilisation or increase the risk of sintering as well as implant loosening or migration.
- Failure to achieve the desired result of surgery, malpositioning of the implant
- Fracture, micro fracture, resorption, damage to or penetration of a vertebral body above or below the treated segment/s
- Reduction of bone density
- Change in the curvature and stiffness of the spine
- Loss or increase in spinal mobility or functions
- Physiological limitations, such as adjacent segment degeneration
- Bending, fracture or loosening of the implant components, as well as abrasion

- Delayed bone growth or no visible regeneration and pseudarthrosis, delayed healing or lack of healing
- Neurological deterioration with paralysis, functional failures, loss of sensation and pain, e.g. caused by scarring
- Soft-tissue lesions
- Pressure on surrounding tissue or organs
- The occurrence of laceration of the dura with risk of cerebrospinal fluid loss or a CSF fistula
- Radiculopathy
- Superficial or deep infection and inflammation as discitis, arachnoiditis, etc.
- Inability to perform one's daily tasks
- Dysfunctions of the bladder and rectum
- Temporary or permanent retrograde ejaculation in men
- Damage to the reproductive system, sterility and sexual dysfunctions
- Deep vein thrombosis, thrombophlebitis, pulmonary embolus
- Development of respiratory problems, including pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.
- Lung injuries
- Injuries of the large abdominal vessels
- Bleeding and/or haematomas
- Allergic reaction to the materials used
- Foreign body reaction to the implants, including possible tumour formation, autoimmune disease and or scarring
- Metallosis
- Bursitis
- Death
- as well as all common surgical risks.

PACKAGING, LABELING, TRANSPORT AND STORAGE

- The implant components must be handled, transported and stored with due care. Damages to the implant packaging or damages to the implant itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increase internal loads, which can result in a fracture of the implant.
- The implants and instruments should be stored at room temperature. Environmental factors, such as salt-laden air, humidity, chemicals, etc., must not be allowed to act on the implants.
- Before the surgery, a careful inspection of the components of the SAMSON® Vertebral Body Replacement system to be used must be carried out in order to avoid damages due to storage, transport or previous procedures. The sterilization containers, trays and the associated covers must not be damaged.
- The implants are delivered in both

NON STERILE and STERILE packaged condition:

Non sterile delivered implants:

NON STERILE delivered implants are labeled as NON STERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILIZATION).

The implants are supplied as an implant system in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

Sterile delivered implants:

STERILE delivered implants are sterilized by a validated gamma sterilization procedure and labeled as STERILE. Cleaning, disinfection and sterilization prior to use must not be performed.

The implants are supplied packaged in a blister-blister-combination and protected by a covering box. The implants may be used only when the label on the outer packaging and also the inner packaging are intact. If the packaging is damaged or already open, the sterility of the implant is not guaranteed and the implant may not be used.

The implants may not be used when the shelf life indicated has been exceeded.

Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or with damaged sterile packaging is not intended.

HumanTech Spine GmbH assumes no responsibility for the use of re-sterilized implants regardless of the person who performed the re-sterilization or the method used.

The rules of asepsis must be observed when removing the implant from the sterile packaging. The sterile packaging may only be opened immediately before the implant is inserted. It is recommended to always have a replacement implant available. The implant must be removed from the packaging using appropriate aseptic provisos.

CLEANING, DISINFECTION AND STERILIZATION

NON STERILE packaged delivered products of the SAMSON® Vertebral Body Replacement System must be cleaned, disinfected and sterilized prior to use. All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the instruction "Processing non sterile delivered Spine". The latest version can always be found on our homepage <https://www.humantech-spine.de/378-en-IFUs.html>, as well as on request from HumanTech Spine GmbH.

DISPOSAL

The products must be disposed of in accordance with the applicable local regulations, whereby the respective degree of contamination must be taken into account.

PRODUCT COMPLAINTS

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using

the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a SAMSON® implant ever exhibits “malfunctioning” (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, the representative should be informed immediately by phone, fax, or in writing.





Should you have a complaint, please provide the name, item number and batch number of the component along with your name, address and detailed description of the defect in writing.

FURTHER INFORMATION

In the event of complaints, suggestions or indications regarding the content of these instructions for use or the use of the product, please refer to the above-mentioned address.

© HUMANTECH Spine GmbH. All rights reserved.

SYMBOLS

	Manufacturer
	Do not re-use
	Batch number
	Reference number
	Read the documentation
	Keep dry
	Non sterile
	sterilized by irradiation
	Do not re-sterilize
	Date of expiry
	Do not use when packaging is damaged.
	Attention

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