VENUS®

POSTERIOR SPINAL IMPLANT SYSTEM

Gebrauchsanweisung	DE	2 - 5
Instruction for use	EN	6 - 9
Kullanım için talimatlar	TR	10 - 13
Manual de instrucciones	ES	14 - 17
Instructies voor gebruik	NL	18 - 21
Instruções de utilização	PT	22 - 25
Lietošanas instrukcija	LV	26 - 29
Naudojimo instrukcija	LT	30 - 32
Istruzioni per l'uso	IT	34 - 37
Návod k použití	CZ	38 - 41
Οδηγίες χρήσης	EL	42 - 45
Revision		46

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IMPORTANT NOTE

This instruction for use must be read carefully and its content must be adhered to.

SCOPE

The instruction for use is applicable for all sterile and unsterile delivered implants of the VENUS®-System.

BASIC STRUCTURE

The VENUS®-System, which is used as an implant system for spinal surgery and consists of rods, pedicle screws, transverse connectors, rod connectors and extension units as well as different types of hooks. Various forms and sizes of these implants are available so that adaptations can be utilized to take into account the unique anatomy of individual patients. Fenestrated screws contain a series of fenestrations which allows polymethylmethacrylate (PMMA) bone cement to be injected into the treated site.

MATERIAL

Components of the VENUS-System are made out of the titanium alloy Ti6Al4V, according to ISO 5832-3. There are also rods included, that are made of cobalt-chromium alloy CoCr28Mo6 according to ISO 5832-12. The Cementadapter CPS is made out of PEEK material.

INDICATIONS FOR USE

The VENUS®-System is intended for the surgical treatment of diseases and injuries of the spine such as instabilities, degenerative disc disease, degenerative spondylolisthesis, degenerative stenosis, scoliosis, kyphosis, spinal fractures, and spondylitis, tumor as well as all conditions that necessitate revision surgery. The VENUS®mini Subsystem, which is part of the VENUS®-System is suitable for percutaneous surgeries.

LEVEL OF FIXATION

Levels of fixation are for the thoracic, lumbar, and sacral spine.

GENERAL CONDITIONS / INFORMATIONS FOR USE

- The implants must be implanted only by surgeons having undergone the necessary training in spinal surgery. The insertion of implants must be carried out in accordance with the surgical and medical indications, the potential risks, and limitations related to this type of surgery; the contraindications, side effects, and precautions defined, and in the knowledge of the nature and physical, chemical, metallic, metallurgic and biological characteristics of the implants.
- This instruction for use provides essential information, but is not sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a clearly defined post-operative treatment regimen, and to conduct scheduled follow-up examinations.
- The activities and physical activity of the patient have a significant influence on the implant's useful life. Patients must be informed that every activity increases the risk of loss, loosening, dislocation, migration, bending or breakage of implant components. Informing patients about limitations to their activities in the postoperative phase and postoperative monitoring of the patients are crucial factors in assessing the development of the surgical result, the fusion and the condition of the implant. The above-mentioned effects may occur even if a bony fusion occurred and the implant is well integrated and the activity restrictions are complied with. The patient must be informed of this. If such an effect occurs, the surgeon must decide whether a revision of the implant or other measures should be taken, taking into account the well-being of the patient and the risks involved.
- Components of the VENUS® system may not be replaced by components / products from other systems from another source or from a different manufacturer. Furthermore, no direct connection may be made between components / products of the VENUS® system to components of other systems or to components/ products made of materials other than Ti6Al4V or CoCr28Mo6be established. If this is not complied with or if the products are otherwise used or used improperly, HumanTech Spine GmbH assumes no responsibility.
- Never reuse the implants. Even if the implant appears to be intact following revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break.
- Implants that have already come into contact with a patient's bodily fluids or tissues or have been soiled must not be reused.
- Handle removed implants in such a way that their reuse is not possible.
- The volumes of cement used for fenestrated / cannulated screws should ultimately be determined

by the surgeon based on the individual patient anatomy.

- Bending of the rods / revision connector affect the biomechanical properties of the implant. Bending in the area of the fixation of the rod in the Poly- or Monoaxial Screw can have negative influences on the fixation of the rod – Bending in this area has to be avoided.
- The rods of the VENUS® Mini-System have marks. It is not allowed to use the area between the ends of the rod and the marks for the fixation of the pedicle screw. In this area the rods may also not be bent. The system stability can be increased with a ventral support. In the lumbar area, especially when the use of pedicle screws with a diameter of 5,5mm or smaller can't be avoided, a ventral support is strongly recommended.
- The instructions in the operation instructions (Surgical Technique) are to be observed. You can download these from the internet from www.humantech-spine.de or directly from your HumanTech representative. It is advised to use the instruments described in the surgical technique, which are intended to be used with the system. Complications may occur if the implant is inserted with or without the intended instruments. No liability can be accepted for the use of external instruments
- Proceed with extreme caution in the region of the spinal cord and the roots of the nerves, since damage to the nerves can lead to the impairment of neurological functions.
- Breakage, slippage or incorrect use of the instruments or implants can injure the patient or the operating staff or result in more time being required for surgery.
- Residues consisting of implant material and/ or not of implant material should be removed.
- Damaged implants must not be implanted.
- Safety and compatibility of the device in the setting of magnetic resonance (imaging) has not been evaluated. No thermal test or test of migration has been performed on the device in this setting.

CONTRAINDICATIONS

Contraindications may be either relative or absolute. The selection of a particular implant must be weighed carefully against the overall assessment of the patient. The following conditions can have an adverse impact on the chances of successful surgery:

- acute infectious process or significant risk of infection (immunocompromise)
- Signs of local inflammation
- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Mental illness
- Grossly distorted anatomy caused by congenital abnormalities.
- Extreme malalignments, impairing the stability of the

instrumentation

- Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of congenital abnormalities, elevation of sedimentation rate unexplained by other diseases, elevation of white blood cell count (WBC) or a marked left shift in the WBC differential count.
- joint disease, bone absorption, osteopenia, osteomalacia, and/or osteoporosis. Osteoporosis or osteopenia is a relative contraindication since this condition may limit the degree of obtainable correction, stabilization, and/or the amount of mechanical fixation.
- Any neuromuscular disease that would place excess strain on the implant during the healing period.
- Suspected or documented allergy or intolerance to the material used. Appropriate tests should be carried out.
- Any case not needing a bone graft and fusion.
- Any case where the implant components selected for use would be too large or too small to achieve a successful result.
- All cases which require the use of components other than the metals or alloys used in this system.
- Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality.
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance.
- Any patient unwilling or not able to follow postoperative instructions.
- Any case not described in the indications.

SIDE EFFECTS / COMPLICATIONS

The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

- Delayed bone growth or no visible fusion and pseudarthrosis
- Neurological complications, paralysis, soft tissue lesions, and/or migration of the implant
- Breakage, loosening or deformation of the implant, as well as abrasion
- Implant breakage during implantation or when implanted
- Superficial or deep-set infection and inflammatory phenomena
- Allergic reaction to the implant material
- Reduction in bone density
- Neurological and spinal dura mater lesions from surgical trauma
- Genitourinary disorders, gastrointestinal disorders, vascular disorders including thrombus, bursitis, secondary bleeding, myocardial infarction, or death
- Disorder of anatomical structures.
- Fracture of a vertebra, the pedicle, and/or the sacrum.
- Presence of micro-particles around the implants, Metallosis

- Growth of the fused vertebra is altered
- Partial loss of the degree of correction achieved during surgery
- Modification of spinal curvature and stiffness of the vertebral column
- Pain, discomfort or abnormal sensations due to the presence of the implant.
- Pressure on the skin caused by components located in positions with insufficient tissue coverage over the implant, with potential penetration of the skin.
- Fracture, micro fracture, resorption, damage to or penetration of a vertebral body above or below the treated segment/s
- Physiological limitations, such as joint degeneration
- Bleeding and/or haematomas
- Revision surgery
- Development of respiratory problems, including pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.

PACKAGING, LABELING, TRANSPORT AND STORAGE

- The implant components must be handled, transported and stored with due care. Damages to the implant packaging or damages to the implant itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increase internal loads, which can result in a fracture of the implant.
- The implants and instruments should be stored at room temperature. Environmental influences such as salty air, humidity, chemicals etc. are not allowed to act on the implants.
- Before the surgery, a careful inspection of the components of the VENUS® system to be used must be carried out in order to avoid damages due to storage, transport or previous procedures. The sterilization containers, trays and the associated covers must not be damaged.

The implants are delivered in both NON STERILE and STERILE packaged condition:

Non sterile delivered implants:

NON STERILE delivered implants are labeled as NON STERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILIZATION). The implants are supplied as an implant system in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

Sterile delivered implants:

STERILE delivered implants are sterilized by a validated gamma sterilization procedure and labeled as STERILE. Cleaning, disinfection and sterilization prior to use must not be performed. The implants are supplied individually packaged in a blister and protected by a covering box.

The implants may be used only when the label on the outer packaging and also the inner packaging are intact. If the packaging is damaged or already open, the sterility of the implant is not guaranteed and the implant may not be used. The implants may not be used when the shelf life indicated has been exceeded. Processing, reprocessing, sterilization or re-sterilization of the products after opening the sterile packaging or with damaged sterile packaging is not intended. HumanTech Spine GmbH assumes no responsibility for the use of re-sterilized implants regardless of the person who performed the re-sterilization or the method used. The rules of asepsis must be observed when removing the implant from the sterile packaging. The sterile packaging may only be opened immediately before the implant is inserted. It is recommended to always have a replacement implant available. The implant must be removed from the packaging using appropriate aseptic provisos.

CLEANING, DISINFECTION AND STERILIZATION

NON STERILE packaged delivered products of the VENUS® System must be cleaned, disinfected and sterilized prior to use. All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the instruction "Processing unsterile delivered Spine". The latest version can always be found on our homepage <https://www.humantech-spine.de/378-en-IFUs.html>, as well as on request from HumanTech Spine GmbH.

DISPOSAL

The products must be disposed of in accordance with the applicable local regulations, whereby the respective degree of contamination must be taken into account.

PRODUCT COMPLAINTS

Any person involved in healthcare (e.g. customer or user of this product system) who has complaints of any kind or is dissatisfied using the product with respect to quality, identity, strength, durability, safety, effectiveness or function should notify the appropriate HumanTech representative. If a VENUS® implant should ever malfunction (i.e. fails to fulfill the performance specifications or does not function in the foreseen manner), or if the above is suspected, the HumanTech representative should be informed immediately.


If a HumanTech product should ever malfunction in a way that results in the death or serious injury of the patient or if it has contributed to the above, the representative is to be informed immediately by telephone, fax or in written form.

If you have any complaints then we ask you to state the name, article number and batch number of the components as well as your name and address along with a detailed description of the fault and forward this information to us in writing.

FURTHER INFORMATION

If you have any complaints, suggestions or other points about the contents of this instruction for use or the use of the product, please contact the above address.
© HumanTech Spine GmbH. All rights reserved.

SYMBOLS

	Manufacturer
	Do not re-use
	Batch number
	Reference number
	Read the documentation
	Keep dry
	Non sterile
	sterilized by irradiation
	Do not re-sterilize
	Date of expiry
	Do not use when packing is damaged
	Attention

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