# ADONIS®

## LUMBAR INTERBODY FUSION SYSTEM

Gebrauchsanweisung	DE	2 - 9
Instruction for use	EN	10 - 16
Návod k použití	CS	
Manual de instrucciones	ES	
οδηγίες χρήσης	EL	
Használati utasítás	HU	
Istruzioni per l'uso	IT	
Lietošanas instrukcija	LV	
Instructies voor gebruik	NL	
Instruções de utilização	PT	
Instrucțiuni de utilizare	RO	
Návod na použitie	SK	

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## IMPORTANT INFORMATION

This instruction for use and further product-specific information (e.g. surgical technique described in the brochure and separate IFU for instrumentation) must be read carefully and its content must be adhered to. For current product-specific information, refer to: [www.humantech-spine.de](http://www.humantech-spine.de)

This instruction for use provides essential information, but is not sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a clearly defined postoperative treatment regimen, and to conduct scheduled follow-up examinations. For improved readability, the tradename ADONIS® Lumbar Interbody Fusion System may be shortened to ADONIS® system throughout the IFU.

## SCOPE

The instruction for use is applicable for all sterile and unsterile delivered implants of the ADONIS® Lumbar Interbody Fusion System.

## INTENDED USE

The implants of the ADONIS® Lumbar Interbody Fusion System are intended to restore and maintain disc height and spinal stability in the lumbar and lumbosacral spine following partial or total discectomy due to diseases of the spine where intervertebral fusion is required.

The device is designed to facilitate bony fusion between two adjacent vertebral bodies by providing structural support.

It is intended for use via posterior, transforaminal, anterior or lateral surgical approaches, depending on the specific implant design.

Implants of the ADONIS® Lumbar Interbody Fusion System are not intended for stand-alone use and must be used with an additional dorsal fixation system (e.g. VENUS®).

## GENERAL PRODUCT DESCRIPTION

The implants of the ADONIS® Lumbar Interbody Fusion System are designed for surgical implantation into the intervertebral space of the lumbar or lumbosacral spine to replace a degenerated or damaged intervertebral disc. The primary functions of the implants are to restore the physiological disc height, provide structural stability, and help to reestablish the natural curvature of the spine following a discectomy.

The implants of the ADONIS® Lumbar Interbody Fusion System are indicated for use in skeletally mature patients and are designed to support fusion between two or more vertebral bodies. The implants serve as a structural support element to facilitate stabilization and promote bony fusion during the normal postoperative healing process after surgical correction of degenerative or other spinal pathologies.

The use of the ADONIS® Lumbar Interbody Fusion System should provide clinical benefit to patients by improving its quality of life through pain reduction, clinically relevant improvement in function, and/or stabilization to support the natural curvature.

Different types of implants are available for the various surgical approaches to the spine:

### posterior approaches:

- ADONIS® PLIF - posterior lumbar interbody fusion
- ADONIS® TLIF - transforaminal lumbar interbody fusion
- ADONIS® UniLIF - unilateral lumbar interbody fusion

### anterior approach:

- ADONIS® ALIF - anterior lumbar interbody fusion

The implants are available in different dimensions and angulations, so that the individual patient's anatomy could be taken into account.

The HumanTech instruments and accessories, intended for use with the implants of the ADONIS® Lumbar Interbody Fusion System, and the surgical technique are described in the brochure and are mandatory for implantation.

## INDICATIONS FOR USE

Implants of the ADONIS® Lumbar Interbody Fusion System are intended for use in patients in the following indications when structural stabilization and spinal fusion are indicated.

- Degenerative disc disease (DDD)
- Spondylolisthesis grade I or II
- Recurrent disc herniation
- Revision surgeries, including failed back surgery syndrome or pseudoarthrosis
- Spinal stenosis with segmental instability
- Spinal deformities (e.g., degenerative scoliosis)
- Structural defects or instability due to trauma
- Post-infectious or post-inflammatory conditions (e.g., spondylitis, spondylodiscitis)
- Structural defects due to tumor, trauma, or prior surgery, where the intervertebral disc space is intact or has been removed.

## CONTRAINDICATIONS

Contraindications may be either relative or absolute.

### Absolute contraindications:

- Suspected or documented allergy or intolerance to the materials used
- Medical or surgical condition which would preclude a potential benefit of spinal implant surgery, such as the presence of congenital abnormalities
- Severe destructive or inflammatory conditions (e.g. tumors, spondylitis, unstable fractures) requiring vertebral body replacement and additional stabilization due to significant structural damage to the vertebral body or motion segment.
- Anatomical structures expected to interfere with the usage of the implant (e.g. extreme anatomic abnormalities) or to impair the instrumentation
- If the implant is expected not being able to perform adequately after implantation
- If the implant components selected for use would be too large or too small or the selected angulation is insufficient to achieve a successful result
- All cases in which fusion is not required or fusion already exist in the segments treated
- Missing intact adjacent segments
- When the use of components of different metals or alloys is necessary
- Spondylolisthesis that cannot be reduced to the first degree.

### Relative contraindications:

- Relative contraindications are conditions that may decrease the surgical success and its outcomes.
- The responsible physician has to weigh the risks and benefits of an implantation in each individual case. If the implantation is performed despite the presence of a relative contraindication, special care or procedures may be required.
- Fever or leucocytosis
- Acute or chronic infection
- Significant risk of infection (e.g. in case of immune deficiencies)
- Signs of local inflammation
- Open wounds
- Wound healing disorders
- Pregnancy: Operations on pregnant women must be avoided if possible. If they are nevertheless performed, they require special care or procedures.
- Local bone tumours
- Inadequate tissue structure at the surgical site or inadequate bone or bone quality
- Known hereditary or acquired bone brittleness or calcification problems.

- Concomitant diseases or conditions that could negatively affect the function and success of the implant, e.g.
  - advanced age
  - smoking
  - diabetes
  - osteoporosis, osteopenia, osteomalacia, joint disease, bone absorption
  - obesity
  - Drug or alcohol abuse
  - chronic kidney disease
  - depression or conditions of reduced mental health
  - neuromuscular diseases that would place excess strain on the implant during the healing period
- If the patient is unwilling or not able to follow postoperative instructions
- Physical or medical conditions that would prevent beneficial surgical outcome.

## INTENDED PATIENT GROUP

All patients that are at least 18 years old and skeletally matured at the time of surgery.

## INTENDED USERS

The devices of the ADONIS® system may only be implanted by surgeons who have completed the necessary training in spinal surgery.

The

- professional judgement and/or clinical skills and experience of the surgeon with regard to careful patient selection
- preoperative planning and implant selection
- knowledge of the anatomy and biomechanics of the spine
- an understanding of the material and the mechanical properties of the used implants
- training and skills in spinal surgery and the use of the instruments required for inserting the implants
- the surgeon's ability to gain the patient's consent, adhere to a clearly defined post-operative treatment regimen, and to conduct scheduled follow-up examinations

is/are necessary for the use of the devices of the system.

The surgeons are supported in the operating theater environment by surgical assistants. General knowledge in surgical procedures, the surgical technique and special training and skills regarding the correct use of the devices of the system and the use of the instruments required for inserting the implants is/are necessary for the use of the devices.

Unsterile delivered products must be cleaned, disinfected and sterilized before use on the patient by persons responsible for processing procedures.

## INTENDED ENVIRONMENT

The devices have to be used in an operating theater.

A sterile environment and compliance with the general hygiene regulations that apply in operating theaters are basic requirements for using the devices.

It must be assured that the surgical field is fully illuminated with a view of the situs.

The implantation of the product in the patient takes place under general anesthesia. In order to be able to perform the operation under general anesthesia, anesthetic gas or other anesthetics are required.

It has to be assured that a C-Arm (X-ray) and the systemspecific instruments are available. Furthermore, general surgical instruments are required for the performance of the surgery.

## SIDE EFFECTS / ADVERSE EVENTS

### Possibly related to the implant:

- Fracture, bending or loosening of the implant or its components
- Implant migration
- Implant subsidence
- Instability of the instrumentation
- Allergic reaction to the implant material
- Foreign body reaction to the implants, including possible tumour formation, autoimmune disease and/ or scar formation
- Metallosis due to intra- or postoperative abrasion or corrosion
- Delayed bone growth or no visible fusion and pseudarthrosis, delayed bone healing or lack of bone healing
- Modification of spinal curvature and stiffness of the vertebral column
- Partial loss of the degree of correction achieved during surgery
- Reduction of bone density due to stress shielding
- Proximal Junctional Kyphosis (adjacent segment degeneration)
- Local pain, discomfort or unusual sensitivity
- Risk of infection or inflammation
- Fracture of a vertebra, micro fracture, damage to or penetration of a vertebral body above or below the treated segment/s
- Vertebral joint degeneration due to altered biomechanics
- Revision surgery.

### Possibly related to the surgical procedure:

- Malpositioning of the implant
- Loss or increase of spinal mobility or functions
- Herniated disc, disc destruction or degeneration in, above or below the treated segment(s)
- Fracture, resorption, necrosis or penetration adjacent vertebrae

- Pressure on or injury of surrounding tissue, nerves or organs
- Neurological complications or deterioration, e.g.:
  - Transient or permanent motor weakness, sensory loss, reflex changes, pain, or paralysis
  - Due to direct injury or irritation of the spinal cord, nerve roots (e.g., radiculopathy), or peripheral nerves
  - Secondary effects such as scarring, dural injury (durotomy), or postoperative inflammation
  - Heterotopic ossification (abnormal bone formation in soft tissue), potentially leading to pain, limited mobility, or nerve compression
  - Injury or lesion of the dura mater with potential cerebrospinal fluid (CSF) loss or formation of a CSF fistula
  - Complex regional pain syndrome (CRPS), potentially resulting in chronic pain, swelling, skin changes, and reduced limb function
- Neurological or neurogenic functional disorders, e.g.:
  - Urinary retention, loss of bladder control, or other urological impairments
  - Dysfunction of the rectum or gastrointestinal system
  - Sexual dysfunction, including retrograde ejaculation in men
  - Damage to the reproductive system or infertility
- Injury or lesion of surrounding soft tissue, blood vessels (e.g. aorta abdominalis, vena cava inferior or other vessels) or lymphatic vessels (including potential lymph leakage)
- Failure to achieve the desired result of surgery
- Reduced implant stability and correction potential due to reduction of bone density, diffuse spinal tumours or infections
- Complications at the bone transplant donor site.

### Possibly general surgical risk:

- Bleeding and/or haematomas
- Lung injuries
- Superficial or deep infections and inflammation as discitis, sepsis, arachnoiditis, etc.
- Change of mental state
- Back or leg cramps
- Headache
- Insomnia
- Vascular disorders including thrombus or myocardial infarction
- Development of respiratory problems, including pulmonary embolism, atelectasis, bronchitis, pneumonia, pneumothorax or transient pleural effusion
- Inability to perform daily tasks
- Death

## WARNINGS AND PRECAUTIONS

### General

- When using X-ray equipment for diagnosis or control, the national limit values for the radiation dose must be observed
- The ADONIS® System consist of implants and instruments. ADONIS® implants must be only used with the intended system-specific instruments and accessories described in the surgical technique.

### MR safety information

Implants of the ADONIS® Lumbar Interbody Fusion System have been examined in non-clinical-testings in accordance with ASTM norms and determined as “MR conditional”. Patients with ADONIS® implants can be safely examined via MRI on following conditions:

- Static magnetic field of 1.5 T to 3 T
- Spatial gradient magnetic field up to max. 1100 Gauss/cm or 11 T/m
- Whole body specific absorption rate (SAR) up to max. 2 W/kg in normal operating mode and at max. 15 min continuous scan duration

Implants of the ADONIS® Lumbar Interbody Fusion System can heat up by up to 3.5°C with continuous scanning time of 15 minutes and under the scanning conditions defined above. Under these conditions, a patient can be examined at low risk. To keep the risks of warming low, the scan duration should be kept as short as possible and the SAR as low as possible. The product has not been examined for artifacts in MR environment. MR image quality may be reduced if the examined area is in relative proximity to ADONIS® system implants.

A thorough risk-benefit assessment by the attending doctor must be carried out.

### Preoperative

- The use of implants for implantation must be decided upon in accordance with the surgical and medical indications, the potential risks and restrictions associated with this type of surgical procedure, as well as in knowledge of the contraindications, side effects and defined precautions, and in awareness of both the nature and the physical, chemical and biological characteristics of the implant.
- In older patients, the user must weigh the benefits against the risks, e.g. osteoporosis, slower healing and increased comorbidities such as heart disease or other age-related conditions that carry a higher risk of surgical complications such as infections or cardiovascular problems.
- The use of ADONIS implants during pregnancy should be avoided whenever possible, as it poses significant risks to both mother and fetus. These risks include infections, thromboembolic events,

implant-related complications, and potential harm to fetal development. Physicians are obligated to provide comprehensive counseling, carefully assess conservative treatment alternatives, and implement strategies to minimize blood loss, radiation exposure, and anesthesia-related risks. Careful planning and close monitoring throughout the preoperative, intraoperative, and postoperative phases are essential to avoid complications and to ensure the best possible outcomes for both mother and child.

- Polypharmacy, especially hyperpolypharmacy, can have a negative impact on the surgical outcome and the patient’s health in elderly patients. The number of prescribed medications in this patient group should be considered prior to surgery.
- The correct selection and placement of implants has to be ensured or to be checked based on suitable pre-, intra- and postoperative diagnostic procedures.
- The activities / movement behaviour of the patient has a significant influence on the service life of the implant. The patient must be informed that any activity increases the risk of loosening, dislocation, migration, bending or fracture of the implant or implant components. It is crucial to inform the patient about restrictions in activities in the post-operative phase and to monitor the patients postoperatively in order to assess the development of the fusion and the condition of the implant. Even when a solid bone fusion has occurred, implant components may still bend, break or loosen. Therefore, the patient must be informed that implant components can also bend, break or loosen if the restrictions in activities are properly followed.
- The patient has to be informed in detail about the risks of the procedure and the use of the implants, including any revisions that may be necessary.
- Errors in implant selection can lead to premature clinical implant failure. The number of segments to be served is to be determined carefully.

### Intraoperative

- The instructions in the surgical technique must be adhered to. You can download these online at [www.humantech-spine.de](http://www.humantech-spine.de) or obtain them directly from your HumanTech representative. It is mandatory to use the instruments described in the surgical technique, which are intended to be used with the system. No liability can be assumed for the use of third-party instruments.
- Damaged implants must not be implanted.
- Components of the ADONIS® system may not be replaced by components / products from other systems from another source or from a different manufacturer. Furthermore, no direct connection of components / products of the system to components of other systems may be established. If this is not complied with or if the products are otherwise used or used improperly, HumanTech Spine GmbH assumes no responsibility.

- Breakage, slippage or incorrect use of the instruments or implants can result in an injury of the patient or the surgical personnel or result in an extension of the duration of the surgery.
- Care should be taken to maintain the integrity of the endplates of the vertebral bodies. Damage to the endplates of the vertebral bodies or excessive partial abrasion of the endplates of the vertebral bodies can lead to sintering of the implant and to a loss of segmental stability.
- The sizes of the implant should be selected in a sufficiently large size in order to make use of the maximum contact surface of the vertebral body. This will ensure the greatest possible stability and counteract implant sintering.
- Careful inspection of the ADONIS® implants is recommended before and during surgery to ensure that the implants are undamaged during assembly or previous procedures and have been correctly assembled according to the instructions in the surgical technique.
- Proceed with extreme care in the area of the spinal cord and the nerve roots, as any damage to nerves can result in loss of neurological functions.
- Should an implant component or parts of an instrument loosen or break in situ, these residues must be completely removed.
- Implants that have already come into contact with a patient's body fluids or tissues or have been soiled must not be reused.
- Bone cement must not be used, as the removal of the components becomes difficult or impossible, particularly in revision surgeries. The heat produced by the hardening process can damage or deform the PEEK implants.
- When using PEEK implants particular care in preparing the intervertebral disc space is needed to keep the placement forces low.
- When inserting the PEEK-Ti coated implants, special care is required, as the coating for the improvement of the ingrowth behavior, consists of a rough surface.
- Implantation of the implant must be performed under continuous x-ray control in order to avoid misplacement.
- Implants are used to reconstruct substantial defects in the surgical area during a healing process lasting a maximum of two years. After spinal fusion has been achieved, the implants are part of the new formed body structure. Explantation of the implants is not intended unless complications, implant failure or a lack of fusion occur during the healing phase, which make implant removal necessary.
- A prolonged healing phase, lack of bone fusion or subsequent bone resorption can lead to overloading of the implant, which may result in implant failure
- In case of complications, the surgeon has to decide whether a revision of the implant should be carried out, taking into account the patient's condition and the possible risks involved.
- The implants are intended for single-use only, never reuse the implants. Even if the implant appears to be intact after the revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break.
- Handle removed implants in such a way that their reuse is not possible.
- An implant card must be given to the patient.

## MATERIAL

For products of the ADONIS® PEEK line, the base body is entirely made of PEEK Optima® according to ASTM F2026. To achieve a better radiographic contrast, radiographic markers made out of tantalum according to ASTM F560 have been incorporated into the implant.

For products of the ADONIS® R-PEEK-Ti line, the base body is also made of PEEK Optima® according to ASTM F2026. In addition, the surfaces of the implant have a coating out of pure titanium according to ASTM F1580. To achieve a better radiographic contrast, radiographic markers made out of tantalum according to ASTM F560 have been incorporated into the implant.

Products of the ADONIS® Ti line consist of the titanium alloy Ti6Al4V according to ISO 5832-3.

## PACKAGING, LABELING, TRANSPORT AND STORAGE

- The devices must be handled, transported and stored with care. Damages to the device packaging or damages to the device itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increased internal loads, which can result in a fracture of the device.
- The implants and corresponding instruments should be stored at room temperature. Environmental factors such as salty air, humidity, direct sunlight, chemicals, etc., must not be allowed to act on the devices.
- Before the surgery, a careful inspection of the components of the ADONIS® system to be used must be carried out in order to avoid damages caused by storage, transport or previous procedures. The steriliza-

## Postoperative

- The patient should be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area.
- The patient should be monitored if a change in the area where the implant(s) is/are located is noticed. The treating physician should assess the possibility of clinical implant failure and discuss with the patient the actions required to achieve further healing.
- The permanent stay of the implant in the body can lead to the previous implant-related complications and to possible unknown or unexpected long-term complications

tion containers, trays and the associated covers must not be damaged.

The devices are delivered in both UNSTERILE and STERILE packaged condition:

**Unsterile delivered devices:**

UNSTERILE delivered devices are labelled as UNSTERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILIZATION). The devices are supplied in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

**Sterile delivered devices:**

STERILE delivered devices are sterilized by a validated gamma or ethylene oxide sterilization procedure and labelled as STERILE. Cleaning, disinfection and sterilization prior to use need not be performed.

The implants are supplied individually packaged in a blister-pouch- or blister-blister-combination and protected by an outer carton. The devices may only be used when the label on the outer packaging and also the inner packaging is intact. If the packaging is damaged or already open, the sterility of the device is not guaranteed, and the device must not be used.

The devices must not be used if the indicated expiry date has been exceeded.

Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or if the sterile packaging is damaged is not intended.

HumanTech Spine GmbH assumes no responsibility for the use of resterilized devices regardless of the person who performed the resterilization or the method used.

The rules of asepsis must be observed when removing the devices from the sterile packaging. The sterile packaging may only be opened immediately prior to the insertion of the device. It is recommended that a replacement device is always kept available. The implant must be removed from the packaging using appropriate aseptic precautions.

**CLEANING, DISINFECTION AND STERILIZATION**

UNSTERILE packaged delivered products of the ADONIS® Lumbar Interbody Fusion System must be cleaned, disinfected and sterilized prior to use. For a correct execution of the processing procedure, general hygienic principles regarding the processing of medical devices and the knowledge of the defined steps of the processing instruction are necessary.

All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the instruction “Processing instructions for implants\_HuT Spine”. The

latest version can always be found on our homepage: “<https://humantech-spine.de/378-en-IFUs.html>”, as well as on request from HumanTech Spine GmbH.

**DISPOSAL**

Any HumanTech implant that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol.

Devices must be disposed of as a healthcare medical device in accordance with hospital procedures.

**PRODUCT COMPLAINTS**

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a ADONIS® implant ever exhibits “malfunctioning” (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, must be reported immediately to the manufacturer and the competent authority of the respective country by phone, fax, or in writing.

Should you have a complaint, please provide the name, item number and batch number of the component as well as your name, address, and detailed description of the fault as possible in written form.

**Information to be provided to the patient**

- The patient should be informed about contraindications, warnings, precautions, side effects and complications with HumanTech devices.
- The patient must be informed about MRI compatibility regarding the HumanTech products used.
- The patient should be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area
- An implant card must be given to the patient
- The patient has to be informed that the implant card must be shown to the treating physician prior to every planned examination, diagnostic intervention and surgical or non-surgical treatment.

## FURTHER INFORMATION



In the event of complaints, suggestions, or indications regarding the content of these instructions for use or the use of the product, please refer to address listed on the last page.

After the European Medical Devices Database is online, the summary of safety and clinical performance of ADONIS® implants will be available at the following address „<https://ec.europa.eu/tools/eudamed>“. Until then, you can find them on our homepage “<https://www.humantech-spine.de>”, as well as on request from HumanTech Spine GmbH.

If required, the instructions for use can be ordered in paper form from the address given on the last page. It will be sent by post within seven days.

© HUMANTECH Spine GmbH. All rights reserved.

## EXPLANATION OF SYMBOLS

	Manufacturer
	Date of manufacture
	Do not re-use
<b>LOT</b>	Lot number
<b>REF</b>	Catalogue number
<b>QTY</b>	Quantity
<b>MD</b>	Medical device
<b>UDI</b>	Unique Device Identifier
	Consult instructions for use
	Keep dry
	Non-sterile
<b>STERILE R</b>	Sterilized using irradiation
<b>STERILE EO</b>	Sterilized using ethylene oxide
	Single sterile barrier system with protective packaging inside
	Single sterile barrier system with protective packaging outside
	Do not re-sterilize
	Use-by date
	Do not use if package is damaged and consult instructions for use
	Caution
	MR conditional

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