BONE FILLER

HUMANTECH SPINE

Gebrauchsanweisung	DE	2 - 4
Instructions for use	EN	5 - 7

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IMPORTANT INFORMATION

This instruction for use must be read carefully and its content must be adhered to.

This instruction for use provides essential information but is not sufficient for the use of the Bone Filler. The use of the Bone Filler for applying cement must be decided according to the surgical and medical indications, the potential dangers, and the restrictions associated with this type of surgical intervention as well as in awareness of contraindications, side effects and defined precautionary measures.

SCOPE

The scope is limited to the application of bone cement to vertebral bodies using the Cement Adapter and Fenestrated Pedicle Screws of the VENUS Fixation System.

INTENDED USE

Bone Filler is a medical aid for applying bone cement to vertebral bodies using the Cement Adapter and the Fenestrated Pedicle Screws of the VENUS Fixation System.

GENERAL PRODUCT DESCRIPTION

Bone Filler is part of the VENUS System and intended for use only in conjunction with the Cement Adapter and the Fenestrated Pedicle Screws of the VENUS Fixation System. The Bone Filler is screwed onto the Cement Adapter of the VENUS Fixation System after being filled. The cement is applied by means of the push rod. The Bone Filler has a capacity of 1.5 ml. The Bone Filler is designed for one-time use on one patient.

INDICATIONS FOR USE

Bone Filler is indicated for applying bone cement to vertebral bodies using the Cement Adapter and the Fenestrated Pedicle Screws of the VENUS Fixation System.

CONTRAINDICATIONS

Contraindications are based on the implantable devices or procedure rather than the Bone Filler. Medical conditions that contraindicate the use of the Bone Filler include any medical condition that may contraindicate the surgical procedure itself.

Contraindications include but are not limited to:

- Infection in or around the operative site
- Signs of local inflammation
- Allergy or sensitivity to the materials used

- Use of incompatible materials from other systems
- Use of for anything other than their intended use
- The contraindications of the cement used

INTENDED PATIENT GROUP

At least 18 years old and skeletally mature at the time of surgery.

INTENDED USERS

The Bone Filler only be used by surgeons and competent surgical employees who have the necessary experience in spine surgery.

INTENDED ENVIRONMENT

The devices have to be used in an operating theater. A sterile environment and compliance with the general hygiene regulations that apply in operating theaters are basic requirements for using the devices. It must be ensured that the surgical field is fully illuminated with a view of the situs.

SIDE EFFECTS / ADVERSE EVENTS

Improper handling may render the Bone Filler unsuitable for its intended use. It may also result in component loosening, deformation and/or fracture, or injury to the patient or operating room personnel.

Possible complications include, but are not limited to, the side effects of the VENUS Fixation System, the side effects of the cement used and the following:

- Allergy or sensitivity to the materials used
- Nerve injury, paralysis, pain, injury to soft tissue, visceral organs or joints due to improper use or breakage of the Bone Filler
- Leakage of cement into surrounding tissues and burns associated with this while the cement is hardening

Due to the required mechanical properties, the Bone Filler is made of non-implantable materials. If the Bone Filler breaks, no fragments should remain in the patient. This can lead to postoperative complications such as allergy, infection, or biological complications related to the release of metal components and may require further surgery.

For additional adverse effects, please refer to related product specific IFUs (Implant IFUs).

WARNINGS AND PRECAUTIONS

General

The user has to ensure that the most current versions of the complete product materials, including Surgical Techniques, are available and that all aspects of the system instructions and instrument limitations are

thoroughly understood. These are available at: „www.humantech-spine.de“.

The Bone Filler must be used only for the intended purpose. Use of the Bone Filler for any purpose or in any manner other than those described here may cause the Bone Filler damage or failure, which could result in serious injury or death of the patient.

The Bone Filler should be used only by the intended users. The surgeon and all other persons involved in the use of the products are responsible for having adequate product knowledge based on the current state of the art in their field of activity. This will ensure proper use of the product and prevent risks to the health or safety of patients, users, or third parties.

It is not permitted to use the Bone Filler together with components of other spinal fusion systems from another supplier or another manufacturer. If this is the case, HumanTech Spine GmbH does not assume any liability.

Prior to use

The use of the Bone Filler for applying cement must be decided according to the surgical and medical indications, the potential dangers, and the restrictions associated with this type of surgical intervention as well as in awareness of contraindications, side effects and defined precautionary measures.

The amount of cement used must be chosen by the surgeon depending on the individual anatomy of the patient and the cement used.

Caution: The operating manual of the cement must also be observed in every situation.

The instructions in the operation instructions (Surgical Technique of the used VENUS system) are to be observed. These are also available at: „www.humantech-spine.de“.

An inspection and an assembly attempt are recommended before the operation in order to ensure that the products interact with each other correctly.

If the Bone Filler is to be used with imaging equipment, compatibility should be verified prior to use.

Handling

The Bone Filler must be handled with care and according to their intended use. It must not be overly stressed by twisting or using it as a lever as this can result in damage or breakage of the Bone Filler components, which could affect the safety of the patient, the user or others.

Failure to use the Bone Filler for its intended use, or improper handling may result in serious injury or death to the patient.

Observe general precautionary measures in order to avoid contact with pathogens in blood.

End of life

Bone Filler is distributed STERILE and must not be used if the expiry date given has been exceeded.

The Bone Filler is intended for single use in one patient. It may not be reused, resterilized nor reprocessed. Reusing, resterilizing or reprocessing may lead to the structural integrity of the product being compromised or to the product malfunctioning, which could lead to injury or illness of the patient. Furthermore, reusing, resterilizing or reprocessing may lead to an increased risk of contamination of the product and/or to the patient becoming infected, including but not limited to passing on infectious diseases from one patient to another.

MATERIAL

Different materials are used for the Bone Filler. For example, mechanically stressed components are made of various stainless steels, while plastics are used for gripping surfaces, or as a seal.

PACKAGING, LABELING, TRANSPORT AND STORAGE

- Bone Filler must be handled, transported and stored with care. Damage to the product can significantly reduce its stability. It can lead to cracks and/or increased internal tensions which can result in the rupture of the product.
- Bone Fillers should be stored between +5°C and +40°C.
- Environmental factors such as salt-laden air, moisture, chemicals etc. must not affect the products.

Bone Filler is distributed STERILE and must not be used if the expiry date given has been exceeded. If the packaging or the inner packaging is damaged or the packaging has already been opened, the sterility of the Bone Filler is not guaranteed and it must not be used.

DISPOSAL

The bone filler must never be used again. The bone filler and its packaging must be disposed of in accordance with the applicable local regulations and environmental legislation considering the respective degree of contamination.

PRODUCT COMPLAINTS

Any person operating within the healthcare system (e.g. customer or user of Bone Filler) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant

HumanTech representative.

If the Bone Filler ever exhibits “malfunctioning” (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If the Bone Filler ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, this must be reported immediately to the manufacturer and the competent authority of the respective country by phone, fax, or in writing. Should you have a complaint, please provide the batch number of the Bone Filler as well as your name, address, and as detailed as possible description of the fault in written form.



FURTHER INFORMATION

In the event of complaints, suggestions, or indications regarding the content of these instructions for use or the use of the product, please refer to address listed on the last page.

If required, the instructions for use can be ordered in paper form from the address given on the last page. It will be sent by post within seven days.

© HUMANTECH Spine GmbH. All rights reserved.

EXPLANATION OF SYMBOLS

	Manufacturer
	Date of manufacture
	Use-by date
	Do not re-use
	Lot number
	Catalogue number
	Quantity
	Medical device
	Unique Device Identifier
	Consult instructions for use
	Keep dry
	Sterilized using ethylene oxide
	Single sterile barrier system
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Recommended storage temperature: +5°C to +40°C
	Phthalate-free
	Non-pyrogenic

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