HERO[®]

ANTERIOR CERVICAL PLATE SYSTEM

Gebrauchsanweisung	DE	2 - 7
Instruction for use	EN	8 - 12
Revision		13

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IMPORTANT INFORMATION

This instruction for use must be read carefully and its content must be adhered to.

This instruction for use provides essential information but is not sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a clearly defined postoperative treatment regimen, and to conduct scheduled follow-up examinations.

SCOPE

The instruction for use is applicable for all sterile and unsterile delivered implants of the HERO®-System.

INTENDED USE

The HERO® Anterior Cervical Plate System is an implant system for long-term use for temporary anterior operative immobilization, stabilization, fixation and correction of the cervical spine (C2 to C7) while bone fusion occurs.

Devices of the HERO® Anterior Cervical Plate System are intended for stand-alone use.

Devices of the HERO® Anterior Cervical Plate System can also be used as additional anterior support for cervical interbody fusion systems (e.g. TRISTAN®) and cervical vertebral body replacement systems (e.g. LEANDER®).

GENERAL PRODUCT DESCRIPTION

The HERO® Anterior Cervical Plate System consists of bone plates and screws. The HERO® Cervical plates are single-use surgical implantable bone plates available in multiple lengths, suitable for 1-, 2-, 3-, and 4-level surgeries. The single-use screws serve as fixation devices for the plates and are available in several lengths and different configurations. The different implant sizes allow surgeons to accommodate individual anatomical features.

The HERO® Anterior Cervical Plate System is intended for an anterior surgical approach to the cervical spine.

INDICATIONS FOR USE

HERO® Anterior Cervical Plate System is intended for the surgical treatment of diseases and injuries of the cervical spine from C2 to C7 in skeletally mature patients, in particular for indications such as:

- Degenerative disc diseases, e.g. osteochondrosis

- Trauma (i.e. fractures or dislocations)
- Spondylolisthesis
- Tumors
- Pseudarthrosis
- Spinal stenosis
- Cervical myelopathy
- Deformities (i.e. kyphosis, lordosis and/or scoliosis)
- Failed previous fusion

CONTRAINDICATIONS

Contraindications may be either relative or absolute. The selection of a particular implant must be weighed carefully against the overall assessment of the patient.

The following conditions can have an adverse impact on the chances of successful surgery. Contraindications include, but are not limited to, the following:

- Acute infections or significant risks of infections (weakened immune system)
- All destructive and inflammatory diseases of the cervical vertebral body or movement segments such as tumors, spondylitis and fractures with a high degree of instability, with pronounced static and structural alterations of the vertebral body and segment structure, requiring the replacement of the vertebral body and other stabilizing measures for the reconstruction and stabilization of the segment
- Indications of local inflammation
- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Physical diseases
- Drug or alcohol abuse
- Severe anatomical deformities caused by congenital abnormalities
- Any other medical or surgical condition that prevents possible improvement through the use of the implant, such as the presence of congenital abnormalities, fracture near the surgical site, increase in deposition rate not indicated by other clinical pictures, increase in leucocyte count (white blood cells), or marked left shift in the white blood cell differential blood count
- Bone resorption, osteopenia, osteomalacia and/or osteoporosis are relative contraindications, as these can limit the degree of achievable correction/stabilization
- Any neuromuscular disease that would place excess strain on the implant during the healing period.
- Known hereditary or acquired bone brittleness or calcification problems
- Suspected allergy or intolerance, as well as documented allergy or intolerance, to the material used. Appropriate tests must be carried out
- All cases where the use of components of different metals or alloys is necessary
- All cases in which no fusion is required
- All cases in which the implant components selected for use are too large or too small to achieve a satisfactory result
- Any patient with inadequate tissue structure at the surgical site or with inadequate bone stock or bone quality
- All patients in whom the use of the implant would interfere with anatomical structures or restrict

- physiological performance
- All patients not willing to follow postoperative instructions
- All cases not described in the indications

INTENDED PATIENT GROUP

All patients that are at least 18 years old and skeletally matured at the time of surgery.

INTENDED USERS

The devices of the HERO® system may only be implanted by surgeons who have completed the necessary training in spinal surgery.

The

- professional judgement and/or clinical skills and experience of the surgeon with regard to careful patient selection
 - preoperative planning and implant selection
 - knowledge of the anatomy and biomechanics of the spine
 - understanding of the material and the mechanical properties of the used implants
 - training and skills in spinal surgery and the use of the instruments required for inserting the implants
 - the surgeon's ability to gain the patient's consent, adhere to a clearly defined post-operative treatment regimen, and to conduct scheduled follow-up examinations
- is/are necessary for the use of the devices of the system.

The surgeons are supported in the operating theater environment by surgical assistants. General knowledge in surgical procedures, the surgical technique and special training and skills regarding the correct use of the devices of the system and the use of the instruments required for inserting the implants is/are necessary for the use of the devices.

INTENDED ENVIRONMENT

The devices have to be used in an operating theater.

A sterile environment and compliance with the general hygiene regulations that apply in operating theaters are basic requirements for using the devices.

It must be ensured that the surgical field is fully illuminated with a view of the situs.

The implantation of the product in the patient takes place under general anesthesia. In order to be able to perform the surgery under general anesthesia, anesthetic gas or other anesthetics are required.

It has to be assured that C-Arm (X-ray) and the system specific instruments are available. Furthermore, general surgical instruments are required for the performance of the surgery.

SIDE EFFECTS / ADVERSE EVENTS

The side effects and adverse events listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

- Implant loosening
- Implant migration
- Implant breakage, respectively detached parts of the implant and wear
- Malpositioning of the implant

- Instability of the instrumentation
- Bending, fracture or loosening of implant components
- Signs of abrasion, for example occurrence of microparticles in the area of the implant (metallosis)
- Revision surgery or requirement of further surgeries
- Missing intact connection segments, bone resorption, osteopenia, osteomalacia and/or osteoporosis, as well as diffuse spinal tumours and infections, can, in addition to other causes, lead to the reduction of instrumentation stability, limit the degree of obtainable correction and stabilisation or increase the risk of implant loosening or migration.
- Changes in the curvature and stiffness of the spine
- Loss or increase of spinal mobility or functions
- Partial loss of the degree of correction achieved during the surgery
- Physiological restrictions, such as negative impairment or degeneration of adjacent segment(s)
- Delayed bone growth or no visible fusion and pseudarthrosis, delayed healing or lack of healing
- Herniated disc, disc destruction or degeneration in, above or below the treated segment(s)
- Changed growth of the fused spine
- Fracture, microfracture, resorption, damage to, necrosis or penetration of a vertebral body above or below the treated segment(s)
- Reduction of bone density
- Neurological deterioration with numbness, paralysis, functional failures, loss of sensitivity, reflex changes and pain, e.g. caused by scarring or durotomy
- Tactile hyperesthesia
- Neurological complications like injury or irritation of the spinal cord, the root nerves (radiculopathy) or peripheral nerves
- Soft tissue complications like lesions or heterotopic ossification
- Pressure on or injury of surrounding tissue, nerves or organs
- Neurological or spinal lesion of the dura mater due to surgical trauma with risk of loss of cerebrospinal fluid (CSF) or CSF fistula
- Injury to the lymphatic vessels / lymph leakage
- Bleeding and/or hematomas
- Injuries of the large abdominal and/or other blood vessels
- Development of respiratory problems, for example pulmonary embolism, atelectasis, bronchitis, or pneumonia
- Injury to the oesophagus and coughing and vomiting as well as difficulty swallowing (dysphagia) or hoarseness (dysphonia)
- Superficial or deep-seated infections and inflammation such as discitis, sepsis, arachnoiditis, etc.
- Inability to perform daily tasks
- Deep vein thrombosis, thrombophlebitis
- Urinary retention or loss of bladder control, or other types of impairment of the urological system; dysfunction of the rectum
- Temporary or permanent retrograde ejaculation in men
- Damage to the reproductive system, sterility and sexual dysfunction
- Change of mental state

- Allergic /hyper-sensitivity reactions to the materials used
- Foreign body reactions to the implant, including possible tumor formation, autoimmune disease and/or scar formation
- Failure to achieve the desired result of surgery, implantation failure
- Pain (e.g. neck pain, upper extremity pain, upper back pain)
- Motor dysfunction
- Pneumothorax
- Sensory dysfunction
- Complex regional pain syndrome (CRPS)
- Wound healing disorders
- Back or leg cramps
- Upper extremity tremor or weakness
- Muscle spasms
- Severe postoperative delirium
- Headache
- Insomnia
- Stroke
- Cardiac response
- Myocardial infarction
- Bursitis
- Death

As well as all general surgical risks, like high blood pressure, coagulation disorder (coagulopathy), blunt trauma, fever, pleural effusion, anesthesia and blood transfusion risks as well as positioning damage.

WARNINGS AND PRECAUTIONS

General

- The user has to ensure that the most current versions of the complete product materials are on hand and considered. These are also available at:
[„https://www.humantech-spine.de/“](https://www.humantech-spine.de/)
- When using X-ray equipment for diagnosis or control, the national limit values for the radiation dose must be observed.

MR safety information

HERO® implant system implants have been examined in non-clinical-testings in accordance with ASTM standards and determined as “MR conditional”. Patients with HERO® implants can be safely examined via MRI on following conditions:

- Static magnetic field of 1.5 T
- Spatial gradient magnetic field up to max. 1100 Gauss/cm or 11 T/m
- Whole body specific absorption rate (SAR) up to max. 2 W/kg in normal operating mode and at max. 8 min continuous scan duration

HERO® implant system implant can heat up to 6.0 °C with continuous scanning time of 8 min. and under the scanning conditions defined above.

Or:

- Static magnetic field of 3 T
- Spatial gradient magnetic field up to max. 1100 Gauss/cm or 11 T/m
- Whole body specific absorption rate (SAR) up to max. 2 W/kg in normal operating mode and at max. 15 min

continuous scan duration

HERO® implant system implant can heat up to 3.5 °C with continuous scanning time of 15 min. and under the scanning conditions defined above.

Under these conditions, a patient can be examined at low risk. To keep the risks of warming low, the scan duration should be kept as short as possible and the SAR as low as possible.

The product has not been examined for artifacts in MR environment. MR image quality may be reduced if the examined area is in relative proximity to HERO® system implants.

A thorough risk-benefit assessment by the attending doctor must be carried out.

Preoperative

- The use of implants for implantation must be decided upon in accordance with the surgical and medical indications, the potential risks and restrictions associated with this type of surgical procedure, as well as in knowledge of the contraindications, side effects and defined precautions, and in awareness of both the nature and the physical, chemical and biological characteristics of the implant.
- The correct selection and placement of implants has to be ensured or to be checked based on suitable pre-, intra- and postoperative diagnostic procedures.
- The activities / movement behavior of the patient has a significant influence on the service life of the implant. The patient must be informed that any activity increases the risk of loosening, dislocation, migration, bending or fracture of the implant or implant components. It is crucial to inform the patient about restrictions in activities in the postoperative phase and to monitor the patients postoperatively in order to assess the development of the fusion and the condition of the implant. Even when a solid bone fusion has occurred, implant components may still bend, break or loosen. Therefore, the patient must be informed that implant components can also bend, break or loosen if the restrictions in activities are properly followed.
- The patient has to be informed in detail about the risks of the procedure and the use of the implants, including any revisions that may be necessary.
- Errors in implant selection can lead to premature clinical implant failure. The number of segments to be served is to be determined carefully.

Intraoperative

- The instructions in the surgical technique must be adhered to. You can download these online at: [„https://www.humantech-spine.de“](https://www.humantech-spine.de/) or obtain them directly from your HumanTech representative. It is advised to use the instruments described in the surgical technique, which are intended to be used with the system. No liability can be assumed for the use of third-party instruments.
- Damaged implants must not be implanted.
- Components of the HERO® system may not be replaced by components / products from other systems from another source or from a different manufacturer. Furthermore, no direct connection of components /

products of the system to components of other systems may be established. If this is not complied with or if the products are otherwise used or used improperly, HumanTech Spine GmbH assumes no responsibility.

- Breakage, slippage or incorrect use of the instruments or implants can result in an injury of the patient or the surgical personnel or result in an extension of the duration of the surgery.
- Careful inspection of the HERO® implants is recommended before and during surgery to ensure that the implants are undamaged during assembly or previous procedures and have been correctly assembled according to the instructions in the surgical technique.
- Proceed with extreme care in the area of the spinal cord and the nerve roots, as any damage to nerves can result in loss of neurological functions.
- Should an implant component or parts of an instrument loosen or break in situ, these residues must be completely removed.
- Implants that have already come into contact with a patient's body fluids or tissues or have been soiled must not be reused.
- Bone cement must not be used, as the removal of the components becomes difficult or impossible.
- Bending of the plates can have in an adverse effect on the biomechanical properties of the implant. Deformations in the vicinity of locating holes can adversely affect the anchoring of the system. Any bending in this area is to be avoided.
- Choosing the right size, type, and design of implant contributes significantly to the potential for satisfactory fixation. However, the size, shape, and condition of the patient's bone place limitations on the size, type, and strength of implants.

Postoperative

- The patient should be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area.
- The patient should be monitored if a change in the area where the implant(s) is/are located is noticed. The treating physician should assess the possibility of clinical implant failure and discuss with the patient the actions required to achieve further healing.
- In case of complications, the surgeon has to decide whether a revision of the implant should be carried out, taking into account the patient's condition and the possible risks involved.
- Never reuse the implants. Even if the implant appears to be intact after the revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break.
- Handle removed implants in such a way that their reuse is not possible.
- An implant pass must be given to the patient.

MATERIAL

All implants of the HERO® Anterior Cervical Plate System including all components of preassembled assemblies are

made of titanium alloy Ti6Al4V according to ISO5832-3.

PACKAGING, LABELING, TRANSPORT AND STORAGE

- The devices must be handled, transported and stored with care. Damages to the device packaging or damages to the device itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increased internal loads, which can result in a fracture of the device.
- The implants and corresponding instruments should be stored at room temperature. Environmental factors such as salty air, humidity, direct sunlight, chemicals, etc., must not be allowed to act on the devices.
- The implants and corresponding instruments should be stored at room temperature. Environmental factors such as salty air, humidity, direct sunlight, chemicals, etc., must not be allowed to act on the devices.
- Before the surgery, a careful inspection of the components of the HERO® system to be used must be carried out in order to avoid damages caused by storage, transport or previous procedures. The sterilization containers, trays and the associated covers must not be damaged.

The devices are delivered in both NON-STERILE and STERILE packaged condition:

Unsterile delivered devices:

UNSTERILE delivered devices are labelled as UNSTERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILIZATION).

The devices are supplied in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

Sterile delivered devices:

STERILE delivered devices are sterilized by a validated gamma s terilization procedure and labelled as STERILE. Cleaning, disinfection and sterilization prior to use need not be performed.

The implants are supplied individually packaged in a blister-pouch- or blister-blister-combination and protected by an outer carton. The devices may only be used when the label on the outer packaging and also the inner packaging is intact. If the packaging is damaged or already open, the sterility of the device is not guaranteed, and the device must not be used.

The devices must not be used if the indicated expiry date has been exceeded. Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or if the sterile packaging is damaged is not intended.

HumanTech Spine GmbH assumes no responsibility for the use of resterilized devices regardless of the person who performed the resterilization or the method used. The rules of asepsis must be observed when removing the devices from the sterile packaging. The sterile packaging

may only be opened immediately prior to the insertion of the device. It is recommended that a replacement device is always kept available. The implant must be removed from the packaging using appropriate aseptic precautions.

CLEANING, DISINFECTION AND STERILISATION

UNSTERILE packaged delivered products of the HERO® System must be cleaned, disinfected and sterilized prior to use. All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the instruction "Processing instructions for implants_HuT Spine". The latest version can always be found on our homepage: "<https://www.humantech-spine.de/378-en-IFUs.html>", as well as on request from HumanTech Spine GmbH

DISPOSAL

Any HumanTech implant that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol.

Devices must be disposed of as a healthcare medical device in accordance with hospital procedures.

PRODUCT COMPLAINTS

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a product of the HERO® implant system ever exhibits "malfunctioning" (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, must be reported immediately to the manufacturer and the competent authority of the respective country by phone, fax, or in writing. Should you have a complaint, please provide the name, item number and batch number of the component as well as your name, address, and detailed description of the fault as possible in written form.

Information to be provided to the patient

- The patient should be informed about contraindications, warnings, precautions, side effects and complications with HumanTech devices.
- The patient must be informed about MRI compatibility regarding the HumanTech products used.
- The patient should be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area
- An implant pass must be given to the patient.

FURTHER INFORMATION



In the event of complaints, suggestions, or indications regarding the content of these instructions for use or the use of the product, please refer to address listed on the last page.

After the European Medical Devices Database is online, the summary of safety and clinical performance of HERO® implants will be available at the following address „<https://ec.europa.eu/tools/eudamed>“.

If required, the instructions for use can be ordered in paper form from the address given on the last page. It will be sent by post within seven days.

© HUMANTECH Spine GmbH. All rights reserved.

EXPLANATION OF SYMBOLS

	Manufacturer
	Date of manufacture
	Do not re-use
	Batch number
	Reference number
	Quantity
	Medical device
	Unique Device Identification
	Read the documentation
	Keep dry
	Unsterile
	Sterilized by irradiation
	Single sterile barrier system with protective packaging inside
	Do not re-sterilize
	Date of expiry
	Do not use when packaging is damaged
	Attention
	MR conditional

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