# HERO®

## ANTERIOR CERVICAL PLATE SYSTEM

|                           |    |         |
|---------------------------|----|---------|
| Gebrauchsanweisung        | DE | 2 - 9   |
| Instructions for use      | EN | 10 - 16 |
| Návod k použití           | CS |         |
| Manual de instrucciones   | ES |         |
| οδηγίες χρήσης            | EL |         |
| Használati utasítás       | HU |         |
| Istruzioni per l'uso      | IT |         |
| Lietošanas instrukcija    | LV |         |
| Instructies voor gebruik  | NL |         |
| Instruções de utilização  | PT |         |
| Instrucțiuni de utilizare | RO |         |
| Návod na použitie         | SK |         |

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## IMPORTANT INFORMATION

These instructions for use and further product-specific information (e.g. Surgical Technique described in the brochure and separate IFU for instrumentation) must be read carefully and its content must be adhered to. For current product-specific information, refer to:  
www.humantech-spine.de

These instructions for use provide essential information, but are not sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection, preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the implants used; the surgeon's training and skills in spinal surgery and in the use of the instruments required for implantation; his ability to obtain the patient's consent, to ensure the patient's cooperation, and to implement a clearly defined postoperative treatment programme and scheduled follow-up examination.

For improved readability, the tradename HERO™ Anterior Cervical Plate System may be shortened to HERO™ system throughout the IFU.

## SCOPE

The instructions for use are applicable for all sterile and unsterile delivered implants of the HERO™ Anterior Cervical Plate System.

## INTENDED USE

The implants of the HERO™ Anterior Cervical Plate System are intended for anterior stabilization of the cervical spine (C3 to C7) during the spinal fusion process following procedures of decompression and/or restoration of spinal segments. They are designed to provide mechanical support to maintain spinal corrections once achieved and to facilitate healing by fixation and immobilization of spinal segments.

They can be used as additional anterior support for cervical interbody fusion systems (e.g. TRISTAN™) and cervical vertebral body replacement systems (e.g. LEANDER™).

The implants are intended for long-term use, with mono- and multi-segmental treatments possible.

## GENERAL PRODUCT DESCRIPTION

The HERO™ Anterior Cervical Plate System is intended to provide anterior segmental stability to support bony fusion and the conservation of decompression and/or spinal correction in the cervical spine. The plates are fixed to the vertebral bodies (C3-C7) by screws in an anterior surgical

approach immobilizing the treated segments.

The plates are available in various sizes for one to four levels and curvatures to match the patient's individual anatomy. The HERO™ screws are available in different sizes and can be inserted at variable angles. Mono- or bicortical screw fixation is possible. HERO™ Single and Expansion Screws shall only be used for monocortical screw fixation. Expansion Screws are designed for bone of lower quality and revisions and shall not be used in hard bone. Locking Screws are intended to prevent HERO™ Standard and Expansion Screws from screw back-out, whereas Single Screws are locked in place via an integrated screw/plate interlocking mechanism.

The HERO™ Anterior Cervical Plate System is designed to contribute indirectly to pain reduction and functional improvement, thereby improving the patient's quality of life, due to its stabilizing supporting role of maintaining previous surgical achievements and facilitating healing.

The HumanTech instruments and accessories, intended for use with the implants of the HERO™ Anterior Cervical Plate System, and the Surgical Technique are described in detail in the brochure and are mandatory for implantation.

## INDICATIONS FOR USE

A treatment of the following diseases and injuries of the cervical spine (C3-C7) with implants of the HERO™ Anterior Cervical Plate System is indicated if (additional) anterior mechanical stabilization is radiographically or clinically required to support fusion, healing or neurological protection. Typically the devices are used in the context of ACDF (Anterior Cervical Discectomy and Fusion) or ACCF (Anterior Cervical Corpectomy and Fusion).

- Degenerative Disc Disease
- Stenosis (spinal canal stenosis, neuroforaminal stenosis)
- Myelopathy and/or radiculopathy
- Deformities
- Trauma
- Tumor
- Infection
- Spondylolisthesis
- Instabilities
- Failed previous fusion, pseudarthrosis, revisions

## CONTRAINDICATIONS

Contraindications may either be relative or absolute:

### Absolute contraindications

- Suspected or documented allergy or intolerance to the materials used.
- Untreated local infection in the operative field or in the cervical spine.
- Systemic sepsis or uncontrolled infection.
- Severe, decompensated medical illness.

- Medical or surgical conditions (e.g. severe anatomical abnormalities) that are expected to preclude any potential benefit from implantation and to interfere with the use of the implants or instruments.
- Inadequate bone stock or bone quality making implant anchoring impossible.
- Conditions which make an anterior approach impossible.
- If the implant components selected for use are insufficient to achieve a satisfactory result (e.g. selected plates or screws too long or too short, inappropriate bending of the plate, selection of wrong screw design).
- If fusion is not required in the affected segment.

### Relative contraindications

Relative contraindications are conditions that may decrease the surgical success and its outcomes.

The responsible physician has to weigh the risks and benefits of an implantation in each individual case. If the implantation is performed despite the presence of a relative contraindication, special care or procedures may be required.

- Fixation of HERO™ at the vertebral body C2 or T1.
- Fever or leucocytosis.
- Chronic infection.
- Acute infection which is not in the operative field.
- Significant risk of infection (e.g. in case immunocompromising chronic infections).
- Signs of local inflammation.
- Local bone tumors.
- Pre-operative dysphagia.
- Pre-operative airway compromise (e.g. COPD, chronic obstructive pulmonary disease).
- Open wounds.
- Inadequate cervical soft tissue structure or coverage (e.g. poor skin condition) at the surgical site.
- Pregnancy: Operations on pregnant women must be avoided if possible. If they are nevertheless performed, they require special care or procedures.
- If materials (e.g. metals, alloys) or components other than those available in this system are required.
- Concomitant diseases or conditions that could negatively affect the function and success of the implant, e.g.
  - advanced age
  - smoking
  - diabetes
  - osteoporosis, osteopenia, osteomalacia, joint disease, bone resorption
  - morbid obesity
  - chronic kidney disease
  - neuromuscular diseases that would place excess strain on the implant during the healing period
  - psychological disorders or other conditions of reduced mental health

- drug or alcohol abuse
- wound healing disorders
- If the patient is unwilling or not able to follow postoperative instructions.

### INTENDED PATIENT GROUP

All patients that are at least 18 years old and skeletally matured at the time of surgery.

### INTENDED USERS

The devices of the HERO™ Anterior Cervical Plate System may only be implanted by surgeons who have completed the necessary training in spinal surgery.

The

- professional judgement and/or clinical skills and experience of the surgeon with regard to careful patient selection.
- preoperative planning and implant selection.
- knowledge of the anatomy and biomechanics of the spine.
- understanding of the material and the mechanical properties of the used implants.
- training and skills in spinal surgery and the use of the instruments required for inserting the implants.
- surgeon's ability to gain the patient's consent to adhere to a clearly defined post-operative treatment regimen and to conduct scheduled follow-up examinations.

is/are necessary for the use of the devices of the system.

The surgeons are supported in the operating theatre environment by surgical assistants. General knowledge in surgical procedures, the Surgical Technique and special training and skills regarding the correct use of the devices of the system and the use of the instruments required for inserting the implants is/are necessary for the use of the devices.

Unsterile delivered products must be cleaned, disinfected and sterilized before use on the patient by persons responsible for processing procedures.

### INTENDED ENVIRONMENT

The devices have to be used in an operating theatre.

A sterile environment and compliance with the general hygiene regulations that apply in operating theatres are basic requirements for using the devices.

It must be ensured that the surgical field is fully illuminated with a view of the situs.

The implantation of the product in the patient takes place under general anesthesia. In order to be able to perform the surgery under general anesthesia, anesthetic gas or other anesthetics are required.

It has to be assured that C-Arm (X-ray) and the system specific instruments are available. Furthermore, general surgical instruments are required for the performance of the surgery.

## SIDE EFFECTS / ADVERSE EVENTS

The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

### Possibly related to the implant:

- Unintended modification of spinal curvature and stiffness of the vertebral column.
- Delayed bone growth / bone healing / osseointegration, non-union and pseudarthrosis.
- Partial loss of the degree of correction achieved during the surgery.
- Screw pull-out or back-out.
- Migration.
- Reduced bone density due to stress shielding.
- Hardware failure, e.g. breakage, deformation of the implant, bending, loosening of implant or implant components.
- Proximal Junctional Kyphosis (adjacent segment degeneration).
- Vertebral joint degeneration due to altered biomechanics in the spine.
- Local pain, discomfort or abnormal sensitivity.
- Allergic reaction to the implant material.
- Risk of infection or inflammation.
- Autoimmune reaction.
- Metallosis due to intra- or postoperative abrasion or corrosion with possible secondary chronic-inflammatory foreign body reactions as scar formation, local tissue necrosis, formation of neoplasia/tumours and/ or osteolysis.
- Long-term dysphagia due to implant pressure (e.g. on the esophagus).
- Revision surgery or requirement of further surgeries.
- Due to the permanent presence of an implant/implants in the body, there are postoperative risks such as fracture, microfracture, resorption, necrosis, damage to, penetration or subluxation of a vertebral body within, above or below the treated segment(s).

### Possibly related to the surgical procedure:

- Malpositioning of the implant.
- Superficial infection (wound), deep site infection and inflammatory phenomena.
- Dysphagia (difficulty in swallowing).
- Hoarseness (dysphonia).
- Airway compromise.
- Arm or neck cramps.
- Loss or increase of spinal mobility.
- Pain, persistent pain.
- Swelling and edema.
- Herniated disc, disc destruction or disc degeneration in the intervertebral space above or below the treated segment(s).
- Neurological complications or deterioration, e.g.:

- Transient or permanent motor weakness, loss or change of sensitivity or functions, reflex changes, numbness, pain, or paralysis
- Due to direct injury or irritation of the spinal cord, nerve roots (e.g. radiculopathy), or peripheral nerves
- Secondary effects such as scarring, dural injury (durotomy), or postoperative inflammation
- Injury or lesion of the dura mater with potential cerebrospinal fluid (CSF) loss or formation of a CSF fistula.
- Neurological or neurogenic functional disorders, e.g.:
  - Urinary retention, loss of bladder control, or other urological impairments
  - Dysfunction of the rectum or gastrointestinal system
  - Sexual dysfunction
  - Damage to the reproductive system or infertility.
- Heterotopic ossification of soft tissue.
- Injury of or pressure on surrounding bone structures, soft tissue, nerves or organs (e.g. lung, trachea, esophagus).
- Injury or lesion of blood vessels (including aorta, carotid artery, or jugular vein) or lymphatic vessels (including potential lymph leakage).

### Possible general surgical risk:

- Delayed wound healing.
- Superficial wound dehiscence.
- Infection.
- Discomfort.
- Inability to perform daily tasks.
- Bleeding and/or haematomas.
- Vascular disorders including e.g. thrombus, deep vein thrombosis, thrombophlebitis, arrhythmias, high blood pressure and blood pressure fluctuations, myocardial infarction, stroke or death.
- Development of respiratory problems, including e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, pneumothorax, transient pleural effusion or death.
- Temporary kidney dysfunction.
- All risks associated with anesthesia.
- All risks associated with blood transfusion.

## WARNINGS AND PRECAUTIONS

### General

- When using X-ray equipment for diagnosis or control, the national limit values for the radiation dose must be observed.
- The HERO™ System consists of implants and instruments. HERO™ implants must only be used with the intended system-specific instruments and accessories described in the Surgical Technique.

## MR safety information

Implants of the HERO™ Anterior Cervical Plate System have been examined in non-clinical-testings in accordance with ASTM norms and determined as “MR conditional”. Patients with HERO™ implants can be safely examined via MRI under following conditions:

- Static magnetic field of 1.5 T
- Spatial gradient magnetic field up to max. 1100 Gauss/cm or 11 T/m
- Whole body specific absorption rate (SAR) up to max. 2 W/kg in normal operating mode and at max. 8 min continuous scan duration.

Implants of the HERO™ Anterior Cervical Plate System can heat up by up to 6.0 °C with continuous scanning time of 8 min and under the scanning conditions defined above.

Or:

- Static magnetic field of 3 T
- Spatial gradient magnetic field up to max. 1100 Gauss/cm or 11 T/m
- Whole body specific absorption rate (SAR) up to max. 2 W/kg in normal operating mode and at max. 15 min continuous scan duration.

Implants of the HERO™ Anterior Cervical Plate System can heat up by up to 3.5 °C with continuous scanning time of 15 min and under the scanning conditions defined above. Under these conditions, a patient can be examined at low risk. To keep the risk of warming low, the scan duration shall be kept as short as possible and the SAR as low as possible.

The product has not been examined for artifacts in MR environment. MR image quality may be reduced if the examined area is in relative proximity to HERO™ system implants.

The treating doctor must carry out a thorough risk-benefit assessment.

## Preoperative

- The use of implants for implantation must be decided upon in accordance with the surgical and medical indications, the potential risks and restrictions associated with this type of surgical procedure, as well as in knowledge of the contraindications, side effects and defined precautions, and in awareness of both the nature and the physical, chemical and biological characteristics of the implant.
- A decision for an (elective) implantation must be made in agreement between the physician and the patient taking into account the general health condition including comorbidities and the level of activity of the patient. The physician has to inform the patient about the consequences of the implantation, especially that the treatment results in segmental stiffness, and about the possibility that the desired surgical result cannot be achieved.

- The patient has to be informed about any possible side effects of the implants and the risks of the procedure, including any revisions or further surgeries that may be necessary.
- The physician shall inform the patient that several causes, such as adjacent segment degeneration, bone resorption, diffuse spinal tumours or infections, can lead to the reduction of structural support of the cervical plate, increasing the risk of subsidence, implant loosening and migration and making correction and stabilisation difficult or impossible.
- To exclude a suspected allergy against the implant materials, appropriate tests shall be carried out.
- For elderly patients, the responsible physician must carefully weigh the expected benefits against the risks, taking into account the patient's general health condition, increased comorbidities such as osteoporosis, slower healing, and the age-related increased risk of general surgical complications, such as infections or cardiovascular problems.
- The use of HERO™ implants during pregnancy shall be avoided whenever possible, as the intervention poses significant risks to both mother and fetus. These risks include infections, thromboembolic events, implant-related complications, and potential harm to fetal development. Physicians are obliged to provide comprehensive counseling, carefully assess conservative treatment alternatives, and implement strategies to minimize blood loss, radiation exposure, and anesthesia-related risks. Careful planning and close monitoring throughout the preoperative, intraoperative, and postoperative phases are essential to avoid complications and to ensure the best possible outcomes for both mother and child.
- Polypharmacy, especially hyperpolypharmacy, can have a negative impact on surgical outcomes and the patient's overall health in elderly patients. The type and number of prescribed medicines and their potential interactions in this patient group shall be considered prior to surgery.
- The activities / movement behaviour of the patient has a significant influence on the service life of the implant. The patient must be informed that any activity increases the risk of loss, bending or fracture of the implant components. It is crucial to inform the patient about restrictions in activities in the post-operative phase and to monitor the patient post-operatively in order to assess the development of the fusion and the condition of the implant. Even when bone fusion is completed, implant components may still bend, break or loosen. The patient must be informed that implant components can also bend, break or loosen if the restrictions in activities are properly followed.
- The treating physician is responsible for the accurate diagnosis and decision-making for or against surgery as well as the appropriate planning and conductance of the surgery. This implies the patient-specific choice of size and design of the implants, taking into account

the individual pathology and anatomy.

- The correct selection and placement of implants has to be ensured and checked based on suitable pre-, intra- and postoperative imaging procedures.
- Errors in implant selection can lead to premature clinical implant failure. The number of segments to be treated has to be determined carefully. The surgeon has to ensure the availability of sufficient different implant sizes and the required instruments for the surgical procedure.
- The physician shall inform the patient that, if bone quality is low or cannot be assessed pre- or intraoperatively or in case of multisegmental treatments, additional stabilization e.g. with dorsal constructs may be necessary and can carry additional risks. Risks as e.g. dysphagia may increase with the number of treated segments. The treatment of more than 3 segments is generally associated with significantly higher risk of non-fusion and other complications and shall only be performed after careful risk-benefit consideration.
- HumanTech Spine accepts no liability for negative consequences resulting from incorrect diagnosis or indication, an inappropriate choice of implant, incorrect use of implant components and/or an improper Surgical Technique or asepsis.
- In treatments of more than two segments, additional posterior instrumentation may be necessary to decrease the complication rate.

### Intraoperative

- The instructions in the Surgical Technique must be adhered to. You can download these online at [www.humantech-spine.de](http://www.humantech-spine.de) or obtain them directly from your HumanTech representative. It is mandatory to use the instruments described in the Surgical Technique, which are intended to be used with the system. No liability can be assumed for the use of third-party instruments.
- Implants showing any signs of damage or contamination must not be implanted.
- Careful inspection of the implants is recommended before and during surgery to ensure that the implants have been assembled correctly according to the instructions in the Surgical Technique and that they remain undamaged during assembly, previous procedures and insertion.
- Components of the HERO™ system must not be replaced by components / products from other systems from another source or from a different manufacturer. Furthermore, a direct connection of components / products of the HERO™ system to components of other systems must not be established. If this is not observed or if the products are used improperly or otherwise than according to IFU and the mandatory Surgical Technique, HumanTech Spine GmbH assumes no responsibility.
- Breakage of implants or incorrect use of the instru-

ments or implants can result in injuries of the patient or the surgical personnel or in prolonged surgery duration.

- Care shall be taken to maintain the integrity of the vertebral bodies. Damage to the vertebral bodies can lead to insufficient grip and loosening of the screws and to a loss of segmental stability.
- Do not exert too much force to implants or instruments
- Do not bend the plate back and forth. Overbending or repeated bending must be avoided.
- The type and size of screws must be chosen adequately depending on the application (e.g. mono- or bicortical use) and the bone quality. Oversized plates may increase risks such as dysphagia, undersized plates may not cover the defect.
- Implantation of the implants must be performed under X-ray control, as indicated in the Surgical Technique, in order to avoid misplacement. Proceed with extreme caution in the area of the spinal cord and nerve roots, as any damage to nerves can result in the negative impairment or loss of neurological functions. Special caution is required during bicortical screw implantation which shall be continuously monitored by suitable imaging due to the increased risk of spinal cord injuries.
- Extreme caution is required to avoid injuries of the aorta, carotid artery, or jugular vein because this can be life-threatening if not treated immediately.
- Should an implant component or parts of an instrument loosen or break in situ, these residues must be removed completely.
- Implants that have already come into contact with a patient's body fluids or tissues or have been soiled must not be reused.
- The HERO™ implant system is not intended to be used with bone cement.

### Postoperative

- The patient shall be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area.
- The patient shall be monitored if a change in the area of the implant(s) is noticed. The treating physician should assess the possibility of clinical implant failure and discuss with the patient the actions required to achieve further healing.
- A prolonged healing phase, lack of bone fusion or subsequent bone resorption can lead to overloading of the implant, which may result in implant failure.
- The implants are used to provide support and stabilisation during a healing process lasting up to two years following the reconstruction of spinal defects. After spinal fusion has been achieved, the implants are part of the newly formed body structure. Explantation of the implants is not intended unless complications, implant failure or a lack of fusion occur, which can require implant removal. In case of complications, the surgeon has to decide whether a revision of the implant should

be carried out, taking into account the patient's condition and the possible risks involved, e.g. that the removal of an implant already fused with the adjacent vertebral bodies is associated with considerable trauma.

- The implants are intended for single-use only. Never reuse the implants. Even if the implant appears to be intact after the revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break.
- Handle removed implants in such a way that their reuse is impossible.
- An implant card must be given to the patient.

## MATERIAL

All implants of the HERO™ Anterior Cervical Plate System are made out of the titanium alloy Ti6Al4V according to ISO5832-3.

## PACKAGING, LABELING, TRANSPORT AND STORAGE

- The devices must be handled, transported and stored with care. Damages to the device packaging or damages to the device itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increased internal loads, which can result in a fracture of the device.
- The implants and corresponding instruments must be stored at room temperature. Environmental factors such as salty air, humidity, direct sunlight, chemicals, etc., must not be allowed to act on the devices.
- Before the surgery, a careful inspection of the components of the HERO™ System to be used must be carried out in order to avoid damages caused by storage, transport or previous procedures. The sterilization containers, trays and the associated covers must not be damaged.

The devices are delivered in both NON-STERILE and STERILE packaged condition:

### Unsterile delivered devices:

UNSTERILE delivered devices are labelled as UNSTERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILIZATION).

The devices are supplied in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

### Sterile delivered devices:

STERILE delivered devices are sterilized by a validated gamma sterilization procedure and labelled as STERILE. Cleaning, disinfection and sterilization prior to use need not be performed.

The implants are supplied individually packaged in a blister-blister-combination and protected by an outer carton. The implants may only be used when the label on the outer packaging and also the inner packaging is intact. If the packaging is damaged or already open, the sterility of the device is not guaranteed, and the device must not be used. The implants must not be used if the indicated expiry date has been exceeded. Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or if the sterile packaging is damaged is not intended.

HumanTech Spine GmbH assumes no responsibility for the use of resterilized devices regardless of the person who performed the resterilization or the method used. The rules of asepsis must be observed when removing the devices from the sterile packaging. The sterile packaging may only be opened immediately prior to the insertion of the device. It is recommended that a replacement device is always kept available. The implant must be removed from the packaging using appropriate aseptic precautions.

## CLEANING, DISINFECTION AND STERILISATION

UNSTERILE packaged delivered products of the HERO™ System must be cleaned, disinfected and sterilized prior to use. For a correct execution of the processing procedure, general hygienic principles regarding the processing of medical devices and the knowledge of the defined steps of the processing instruction are necessary.

All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the according processing instruction. The latest version can always be found on our homepage: "<https://humantech-spine.de/378-en-IF-Us.html>", as well as on request from HumanTech Spine GmbH.

## DISPOSAL

Any HumanTech implant that has been contaminated by blood, tissue, and/or bodily fluids/matter must never be used again and should be handled according to hospital protocol.

Devices must be disposed of as a healthcare medical device in accordance with hospital procedures.

## PRODUCT COMPLAINTS

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a product of the HERO™ implant system ever exhibits “malfunctioning” (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, must be reported immediately to the manufacturer and the competent authority of the respective country by phone, fax, or in writing.

Should you have a complaint, please provide the name, item number and batch number of the component as well as your name, address, and detailed description of the fault as possible in written form.

#### Information to be provided to the patient

- The patient must be informed about contraindications, warnings, precautions, side effects and complications with HumanTech devices.
- The patient must be informed about MRI compatibility regarding the HumanTech products used.
- The patient must be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area.
- An implant card must be given to the patient.
- The patient has to be informed that the implant card must be shown to the treating physician prior to every planned examination, diagnostic intervention and surgical or non-surgical treatment.

#### FURTHER INFORMATION

In the event of complaints, suggestions, or indications regarding the content of these instructions for use or the use of the product, please refer to address listed on the last page.

After EUDAMED (European Database on Medical Devices) is online, the summary of safety and clinical performance of HERO™ implants will be available at the following address „<https://ec.europa.eu/tools/eudamed>“. Until then, you can find them on our homepage “<https://www.humantech-spine.de>”, as well as on request from HumanTech Spine GmbH.

If required, the instructions for use can be ordered in paper form from the address given on the last page. It will be sent by post within seven days.

© HUMANTECH Spine GmbH. All rights reserved.

#### EXPLANATION OF SYMBOLS

|                                                                                     |                                                                   |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------|
|    | Manufacturer                                                      |
|    | Date of manufacture                                               |
|    | Do not re-use                                                     |
|    | Lot number                                                        |
|    | Catalogue number                                                  |
|    | Quantity                                                          |
|    | Medical device                                                    |
|    | Unique Device Identifier                                          |
|    | Consult instructions for use                                      |
|    | Keep dry                                                          |
|   | Non-sterile                                                       |
|  | Sterilized using irradiation                                      |
|  | Single sterile barrier system with protective packaging inside    |
|  | Do not resterilize                                                |
|  | Use-by date                                                       |
|  | Do not use if package is damaged and consult instructions for use |
|  | Caution                                                           |
|  | MR conditional                                                    |

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