HUMANTECH SPINE **INSTRUMENTATION**

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DE 2 - 6

Instructions for use

EN 7 - 11

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IMPORTANT INFORMATION

This instruction for use must be read carefully and its content must be adhered.

This instruction for use provides essential information but is not sufficient for the use of the instrumentation. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants.

SCOPE

The scope of surgical instruments is the spine. The handles and torque limiters are only to be used in combination with a HumanTech instrument and enable the handling of the instrument and, if necessary, limiting the torque to a specified value.

INTENDED USE

The spinal instruments and corresponding accessories from HumanTech Spine GmbH are always part of a system and intended for use only in conjunction with the associated implant systems from HumanTech Spine GmbH according to the corresponding surgical technique of each implant system.

GENERAL PRODUCT DESCRIPTION

The surgical instruments are used as instrument systems for various applications in spine surgery. They are available in various designs and sizes.

INDICATIONS FOR USE

The spinal instruments and corresponding accessories from HumanTech Spine are part of our implant systems and thus intended and necessary for the surgical treatment of diseases and injuries of the spine. Except for explicitly labeled products, the surgical instruments are reusable and intended for use with the following implant system:

- VENUS® Spinal Fixation System
- ADONIS® Lumbar Interbody Fusion System
- SAMSON® Vertebral Body Replacement System
- TRISTAN® Cervical Interbody Fusion System
- HERO® Anterior Cervical Plate System
- LEANDER® Cervical Vertebral Body Replacement System

These are only intended for manual operation. Use with active or motorized Class IIa propulsion system is not allowed.

CONTRAINDICATIONS

Contraindications are based on the implantable devices or procedure rather than the instruments.

Medical conditions that contraindicate the use of these instruments include any medical condition that may contraindicate the surgical procedure itself.

Contraindications include but are not limited to:

- Infection in or around the operative site
- Signs of local inflammation
- Allergy or sensitivity to instrument materials
- Use of incompatible materials from other systems
- Use of for anything other than their intended use

The instruments must not be used or processed again after contact with patients who suffer from Creutzfeldt-Jakob disease (CJD or vCJD), prion diseases or are infected by the HI-Virus.

INTENDED PATIENT GROUP

- At least 18 years old and skeletally mature at the time of surgery (except VENUSnano Instrumentation).
- As a Guideline: for patients aged from 5 to 11 years and for patients weighing less than 45 kg at the time of surgery. (for VENUSnano Instrumentation only).

INTENDED USERS

The products should only be used by surgeons and competent surgical employees who have the necessary experience in spine surgery. The use of instruments for the implantation must be decided upon in accordance with the surgical and medical indications, the potential risks, and limitations related to this type of surgery; the contraindications, side effects, and precautions defined, and in the knowledge of the nature and metallic, metallurgic and biological characteristics of the implants and instruments.

The surgeons are supported in the operating theater environment by surgical assistants. General knowledge in surgical procedures, the surgical technique and special training and skills regarding the correct use of the devices of the system and the use of the instruments required for inserting the implants is/are necessary for the use of the devices.

Unsterile delivered products must be cleaned, disinfected and sterilized before use on the patient by persons responsible for processing procedures.

INTENDED ENVIRONMENT

The devices have to be used in an operating theater. A sterile environment and compliance with the general hygiene regulations that apply in operating theaters are basic requirements for using the devices. It must be ensured that the surgical field is fully illuminated

with a view of the situs.

In order to be able to perform the surgery under general anesthesia, anesthetic gas or other anesthetics are required.

It has to be assured that C-Arm (X-ray) and the system specific instruments are available. Furthermore, general surgical instruments are required for the performance of the surgery.

SIDE EFFECTS / ADVERSE EVENTS

Improper maintenance, cleaning, or handling may render the instruments unsuitable for their intended use. It may also result in corrosion, loosening of components, deformation and/or fracture, or injury to the patient or operating room personnel.

Possible complications include, but are not limited to, the following:

- Allergy or sensitivity to instrument materials
- Nerve injury, paralysis, pain, injury to soft tissue, visceral organs or joints due to improper use or breakage of instruments
- Infection due to failure to properly clean and sterilize instruments
- Neurological or spinal lesion of the dura mater due to surgical trauma with risk of loss of cerebrospinal fluid (CSF) or CSF fistula, compression of vessels, damage to nerves or adjacent organs due to slippage or misplacement of damaged instruments
- Damage from inadvertent release of springs on certain instruments
- Damage from instruments used to bend or cut in situ due to excessive force applied during use
- Cuts to the gloves or skin of operating room personnel
- Tissue damage to the patient or OR personnel and/or increased OR time due to disassembly of instruments during surgery
- Bone cracks or fractures or inadvertent piercing of bone

Due to the required mechanical properties, most instruments are made of non-implantable materials. If an instrument breaks, no fragments must remain in the patient. This can lead to postoperative complications such as allergy, infection, or biological complications related to the release of metal components and may require further surgery.

For additional adverse effects, please refer to related product specific IFUs (Implant IFUs).

WARNINGS AND PRECAUTIONS

General

The user has to ensure that the most current versions of the complete product materials are on hand and considered. These are also available at: „www.humantech-spine.de”

The HumanTech spinal instruments and accessories that accompany those must be used only for the intended purpose. Use of this instrument for any purpose or in any manner other than those described here may cause instrument damage or failure, which could result in serious injury or death to the patient.

The Instruments should be used only by the intended users. The surgeon and all other persons involved in the use of the products are responsible for having adequate product knowledge based on the current state of the art in their field of activity. This will ensure proper use of the product and prevent risks to the health or safety of patients, users, or third parties.

HumanTech reusable instruments must be discarded after contact with patients with Creutzfeldt-Jakob disease (CJD or vCJD), prion diseases, or HIV infection. These instruments have not been tested or designed to withstand the chemical and thermal exposures recommended to eradicate prions.

Sterile or Single use instruments and accessories must not be reprocessed or reused. Those are clearly marked as sterile or single use.

Prior to use

HumanTech Reusable instruments and related storage are supplied NON-STERILE and must be removed from the transport packaging and must be cleaned, disinfected lubricated, and sterilized prior to initial use. Before each further use, the reusable instruments and accessories must again be cleaned, disinfected, inspected, lubricated, and sterilized.

Before and after each use, the instruments must be thoroughly inspected for visible soiling, moisture or damage such as fractures, corrosion (rust, pitting), discoloration, excessive scratches, flaking, cracks and wear, and for functional reliability. It is extremely important to inspect all working parts. Instruments that show signs of damage or corrosion should be repaired or replaced before further use. Scratches, dents, or other damage may result in instrument breakage or tissue injury. Long, narrow lumens, blind holes, moving and intricate parts require special attention during cleaning and inspection. Sharp instruments should be handled and cleaned with care. Special attention should be paid to cutting edges, tips, joints, ratchets, and any moving parts, as well as missing or removed part numbers.

The torque limiters should be checked to make sure that it disengages with an audible „click“.

If the instruments are to be used with imaging equipment, compatibility should be verified prior to use.

General handling of the instruments in conjunction with the corresponding storage

Contaminated or used reusable Instruments should not be directly placed in a tray and cleaned in a mechanical washer. Contaminated reusable Instruments are such that have blood, tissue and/or bodily fluid/matter in or on its surface. Those must be cleaned separately from HumanTech instrument trays and containers.

Cannulas, blind holes and intricate parts require particular attention during cleaning.

After cleaning and during mechanical washing and thermal disinfection standard hospital trays may be used for instrument positioning.

Do not stack trays of instruments in a mechanical washer.

Steam sterilization must be done within the storage trays provided by HumanTech. The trays have specially designed lips, adapted to each instrument, where the instruments are held.

The trays with their corresponding instruments must be placed in the sterilization containers provided by HumanTech, which have to be closed. Afterwards steam sterilization must be done.

The trays and containers and thus their handling comply with the usual hospital standards.

The accessories from HumanTech are intended for use in combination with the instruments from HumanTech and must not be used for products from other systems.

Handling

The instruments must be handled with care and according to their intended use. They must not be overly stressed by twisting or using them as a lever as this can result in damage or breakage of instrument parts that can affect patient safety. This applies specially to delicate parts of the instruments.

Surgical instruments are subject to corrosion and become functionally impaired following contact with aggressive substances, such as solutions containing iodine or high chlorine content. Do not use steel wool or abrasive cleaners.

Failure to use instruments for their intended use, or to perform functional and visual checks, or improper cleaning and handling, may result in serious patient injury or death.

End of life

End of life of HumanTech reusable Instruments is normally

determined by wear and tear and damage due to use and depends on the frequency of its use and may vary depending on the care with which the product is handled, cleaned and stored. All reusable instruments are exposed to repeated stresses associated with bone contact, impaction, repeated cleaning and sterilization processes. Repeated processing cycles that include ultrasonic, mechanical washing and sterilization have minimal effects on HumanTech reusable Instruments.

Bending, wear, cracks, corrosion (rust, pitting, etc.), excessive scratches, and missing parts may indicate damage and wear to a reusable instrument. Improperly functioning instruments, instruments with unrecognizable markings, missing or removed part numbers, or wear, corrosion, deformation, porosity, or other damage, and excessively worn instruments should not be used. If necessary, they must be sent to HumanTech customer service for exchange or repair.

The same applies to the accessories accompanying the Instruments.

MATERIAL

Different materials are used for instruments, depending on the application. For example, mechanically stressed instruments are made of various stainless steels, while plastics are used for gripping surfaces, sliding surfaces, or tissue protection.

PACKAGING, LABELING, TRANSPORT AND STORAGE

- The devices must be handled, transported and stored with care. Damages to the device packaging or damages to the device itself can significantly reduce the performance, stability and durability of the products. It can lead to cracks and/or increased internal loads, which can result in a fracture of the device.
- The products should be stored at room temperature. Environmental factors such as salty air, humidity, direct sunlight, chemicals, etc., must not be allowed to act on the devices.
- Prior to surgery, a careful inspection of the products to be used must be carried out in order to avoid damages caused by storage, transport or previous procedures. The sterilization containers, trays, and the associated covers must not be damaged.
- Defective instruments and damaged accessories must be separated immediately and returned to HumanTech Spine Customer Service for repair or replace. Instruments and accessories must be cleaned, disinfected and sterilized before being returned to HumanTech Spine. Proof of such treatment must be attached to the outer packaging. A decontamination certificate must always accompany the defective or returned product.

The devices are delivered in both UNSTERILE and STERILE packaged condition:

Unsterile delivered devices:

UNSTERILE delivered devices are labelled as UNSTERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILIZATION).

The devices are supplied in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

Sterile delivered devices:

STERILE delivered devices are sterilized by a validated sterilization procedure and labelled as STERILE. Cleaning, disinfection and sterilization prior to use need not be performed.

The devices are supplied individually packaged and protected by an outer carton. The devices may only be used when the label on the outer packaging and also the inner packaging is intact. If the packaging is damaged or already open, the sterility of the device is not guaranteed, and the device must not be used.

The devices must not be used if the indicated expiry date has been exceeded. Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or if the sterile packaging is damaged is not intended.

HumanTech Spine GmbH assumes no responsibility for the use of resterilized devices regardless of the person who performed the resterilization or the method used. The rules of asepsis must be observed when removing the devices from the sterile packaging. The sterile packaging may only be opened immediately prior to the insertion of the device. It is recommended that a replacement device is always kept available. The products must be removed from the packaging using appropriate aseptic precautions.

CLEANING, DISINFECTION AND STERILIZATION

UNSTERILE packaged delivered products must be cleaned, disinfected and sterilized prior to use. All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the instruction “Processing instructions HuT Spine Instr.”. The latest version can always be found on our homepage: “www.humantech-spine.de“ as well as on request from HumanTech Spine GmbH

DISPOSAL

All HumanTech products that have been in contact with patients with Creutzfeldt-Jakob disease (CJD or vCJD), prion diseases or HIV infection must never be reused and must be disposed of in accordance with the applicable local regulations and environmental legislation.

PRODUCT COMPLAINTS

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a HumanTech product ever exhibits “malfunctioning” (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, this must be reported immediately to the manufacturer and the competent authority of the respective country by phone, fax, or in writing. Should you have a complaint, please provide the name, item number and batch number of the component as well as your name, address, and as detailed as possible description of the fault in written form.

FURTHER INFORMATION

In the event of complaints, suggestions, or indications regarding the content of these instructions for use or the use of the product, please refer to address listed on the last page.

If required, the instructions for use can be ordered in paper form from the address given on the last page. It will be sent by post within seven days.

© HUMANTECH Spine GmbH. All rights reserved.

EXPLANATION OF SYMBOLS

	Manufacturer
	Date of manufacture
	Do not re-use
	Lot number
	Catalogue number
	Quantity
	Medical device
	Unique Device Identifier
	Consult instructions for use
	Keep dry
	Non-sterile
	Sterilized using irradiation
	Single sterile barrier system with protective packaging inside
	Do not re-sterilize
	Use-by date
	Do not use if package is damaged and consult instructions for use
	Caution

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