LEANDER[®]

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IMPORTANT INFORMATION

This instruction for use and further product-specific information (e.g. surgical technique described in the brochure and separate IFU for instrumentation) must be read carefully and its content must be adhered to. For current product-specific information, refer to:

www.humantech-spine.de

This instruction for use provides essential information, but is not sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a clearly defined postoperative treatment regimen, and to conduct scheduled follow-up examinations.

For improved readability, the tradename LEANDER® Cervical Vertebral Body Replacement System may be shortened to LEANDER® system throughout the IFU.

SCOPE

The instruction for use is applicable for all sterile delivered implants of the LEANDER® Cervical Vertebral Body Replacement System.

INTENDED USE

The implants of the LEANDER® Cervical Vertebral Body Replacement System are intended to surgically reconstruct one or more vertebral bodies at the anterior cervical spine (C3-C7) that have been substantially damaged or removed.

Implants of the LEANDER® system are intended for an anterior approach and for long-term use with an additional anterior cervical plate system (e.g. HERO®). Implants of the LEANDER® system are not intended for stand-alone use.

GENERAL PRODUCT DESCRIPTION

The LEANDER® Cervical Vertebral Body Replacement System (VBR) is an implant system to surgically replace one or more cervical vertebral bodies through an anterior approach. The essential function of LEANDER® is to restore and maintain structural integrity, height and natural curvature of the spine following cervical corpectomy. It provides immediate anterior column support and facilitates the re-establishment of sagittal alignment, thereby protecting nerve structures by decompression and maintaining physiological load transfer. The anterior column support offers mechanical stability for spinal fusion.

The implants of the LEANDER® Cervical Vertebral Body Replacement System consist of multiple components: a central base body and two endplates. The base bodies are available in various heights to accommodate different defect sizes. The endplates are offered in different widths and angulations to match individual vertebral anatomy and to support the restoration or preservation of the spine's natural lordotic curvature. The system's modular assembly using plate-holding screws enables intraoperative adaptation to the patient's individual anatomical conditions.

The use of the LEANDER® Cervical Vertebral Body Replacement System should provide clinical benefit to the patient by improving quality of life through pain reduction, clinically meaningful functional improvement, and/or spinal stabilization to support the physiological sagittal alignment. The HumanTech instruments and accessories, intended for use with the implants of the LEANDER® Cervical Vertebral Body Replacement System, and the surgical technique are described in the brochure and are mandatory for implantation.

INDICATIONS FOR USE

The implants of the LEANDER® Cervical Vertebral Body Replacement System are indicated for Anterior Cervical Corpectomy and Fusion (ACCF) procedures in the cervical spine (C3–C7) to reconstruct one or more vertebral bodies that have been resected or damaged due to trauma, fracture, inflammation, infection, tumor, or degenerative disease resulting in cervical instability or neural compression as well as revision procedures.

CONTRAINDICATIONS

Contraindications may be either relative or absolute. The selection of a specific implant should be based on a comprehensive assessment of the patient.

Absolute contraindications are:

- Suspected or documented allergy or intolerance to the materials used
- Untreated local infection in the operative field or in the cervical spine
- Systemic sepsis or uncontrolled infection
- Severe, decompensated medical illness
- Anatomic structures expected to interfere with the usage of the implant (e.g. extreme anatomic deformities) or to impair the instrumentation or to restrict physical performance of the patient
- Inadequate bone stock or bone quality making implant anchoring impossible
- Conditions which make an anterior approach impossible
- If the selected implant components are too large or too small or if the selected angulations of the endplates are insufficient to achieve a successful result
- If corpectomy and fusion is not required or if fusion already exists in the affected segment.

Relative contraindications are conditions that may decrease the surgical success and its outcomes.

The responsible physician has to weigh the risks and benefits of an implantation in each individual case. If the implantation is performed despite the presence of a relative contraindication, special care or procedures may be required.

- Fever or leucocytosis
- Chronic infection
- Acute infection which is not in the operative field
- Significant risk of infection (e.g. in case immunocompromising chronic infections)
- Signs of local inflammation or bone tumors in the adjacent segments
- Pre-operative dysphagia
- Pre-operative airway compromise (e.g. COPD, chronic obstructive pulmonary disease)
- Open wounds
- Inadequate cervical soft tissue structure or coverage (e.g. poor skin condition) at the surgical site
- Pregnancy: Operations on pregnant women must be avoided if possible. If they are nevertheless performed, they require special care or procedures
- If the use of other components, metals or alloys available in this system is required
- Concomitant diseases or conditions that could negatively affect the function and success of the implant, e.g.
 - advanced age
 - smoking
 - diabetes
 - osteoporosis, osteopenia, osteomalacia, joint disease, bone resorption
 - morbid obesity
 - chronic kidney disease
 - neuromuscular diseases that would place excess strain on the implant during the healing period
 - psychological disorders or other conditions of reduced mental health
 - drug or alcohol abuse
 - wound healing disorders
- If the patient is not willing or able to follow postoperative instructions

INTENDED PATIENT GROUP

All patients that are at least 18 years old and skeletally matured at the time of surgery.

INTENDED USERS

The devices of the LEANDER® system may only be implanted by surgeons who have completed the necessary training in spinal surgery.

The

- professional judgement and/or clinical skills and experience of the surgeon with regard to careful patient selection

- preoperative planning and implant selection
- knowledge of the anatomy and biomechanics of the spine
- understanding of the material and the mechanical properties of the used implants
- training and skills in spinal surgery and the use of the instruments required for inserting the implants
- surgeon's ability to gain the patient's consent to adhere to a clearly defined post-operative treatment regimen and to conduct scheduled follow-up examinations is/are necessary for the use of the devices of the system.

The surgeons are supported in the operating theater environment by surgical assistants. General knowledge in surgical procedures, the surgical technique and special training and skills regarding the correct use of the devices of the system and the use of the instruments required for inserting the implants is/are necessary for the use of the devices.

Unsterile delivered products must be cleaned, disinfected and sterilized before use on the patient by persons responsible for processing procedures

INTENDED ENVIRONMENT

The devices have to be used in an operating theater.

A sterile environment and compliance with the general hygiene regulations that apply in operating theaters are basic requirements for using the devices.

It must be ensured that the surgical field is fully illuminated with a view of the situs.

The implantation of the product in the patient takes place under general anesthesia. In order to be able to perform the surgery under general anesthesia, anesthetic gas or other anesthetics are required.

It has to be assured that C-Arm (X-ray) and the system specific instruments are available. Furthermore, general surgical instruments are required for the performance of the surgery.

SIDE EFFECTS / ADVERSE EVENTS

The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

Possibly related to the implant:

- Unintended modification of spinal curvature and stiffness of the vertebral column
- Delayed bone growth/bone healing/osseointegration, non-union and pseudarthrosis
- Partial loss of the degree of correction achieved during the surgery
- Subsidence
- Migration
- Reduced bone density due to stress shielding
- Hardware failure, e.g. breakage, deformation of the

- implant, loosening of implant or implant components
- Proximal Junctional Kyphosis (adjacent segment degeneration)
- Vertebral joint degeneration due to altered biomechanics in the spine
- Local pain, discomfort or abnormal sensitivity
- Allergic reaction to the implant material
- Risk of infection or inflammation
- Autoimmune reaction
- Metallosis due to intra- or postoperative abrasion or corrosion with possible secondary chronic-inflammatory foreign body reactions as scar formation, local tissue necrosis, formation of neoplasia/tumours and/or osteolysis
- Long-term dysphagia due to implant pressure
- Revision surgery or requirement of further surgeries
- Due to the permanent presence of an implant/implants in the body, there are postoperative risks such as fracture, microfracture, resorption, necrosis, damage to, penetration or subluxation of a vertebral body above or below the treated segment(s)

Possibly related to the surgical procedure:

- Malpositioning of the implant
- Superficial infection (wound), deep site infection and inflammatory phenomena
- Dysphagia (difficulty in swallowing)
- Hoarseness (dysphonia)
- Airway compromise (coughing, shortness of breath, risk of suffocation)
- Arm or neck cramps
- Loss or increase of spinal mobility
- Pain, persistent pain
- Swelling and edema
- Herniated disc, disc destruction or disc degeneration in the intervertebral space above or below the treated segment(s)
- Neurological complications or deterioration, e.g.:
 - Transient or permanent motor weakness, loss or change of sensitivity or functions, reflex changes, numbness, pain, or paralysis
 - Due to direct injury or irritation of the spinal cord, nerve roots (e.g. radiculopathy), or peripheral nerves
 - Secondary effects such as scarring, dural injury (durotomy), or postoperative inflammation
 - Injury or lesion of the dura mater with potential cerebrospinal fluid (CSF) loss or formation of a CSF fistula
- Neurological or neurogenic functional disorders, e.g.:
 - Urinary retention, loss of bladder control, or other urological impairments
 - Dysfunction of the rectum or gastrointestinal system
 - Sexual dysfunction
 - Damage to the reproductive system or infertility
- Heterotopic ossification of soft tissue
- Injury of or pressure on surrounding bone structures,

soft tissue, nerves or organs (e.g. lung, trachea, esophagus)

- Injury or lesion of blood vessels (including aorta, carotid artery, or jugular vein) or lymphatic vessels (including potential lymph leakage)

Possible general surgical risk:

- Delayed wound healing
- Superficial wound dehiscence
- Infection
- Discomfort
- Inability to perform daily tasks
- Bleeding and/or haematomas
- Vascular disorders including e.g. thrombus, deep vein thrombosis, thrombophlebitis, arrhythmias, high blood pressure and blood pressure fluctuations, myocardial infarction, stroke or death
- Development of respiratory problems, including e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, pneumothorax, transient pleural effusion or death
- Temporary kidney dysfunction
- All risks associated with anesthesia
- All risks associated with blood transfusion

WARNINGS AND PRECAUTIONS

General

- When using X-ray equipment for diagnosis or control, the national limit values for the radiation dose must be observed.
- The LEANDER® System consist of implants and instruments. LEANDER® implants must be only used with the intended system-specific instruments and accessories described in the surgical technique.

MR safety information

Implants of the LEANDER® Cervical Vertebral Body Replacement System have been examined in non-clinical-testings in accordance with ASTM norms and determined as “MR conditional”. Patients with LEANDER® implants can be safely examined via MRI on following conditions:

- Static magnetic field of 1.5 T to 3 T
- Spatial gradient magnetic field up to max. 1100 Gauss/cm or 11 T/m
- Whole body specific absorption rate (SAR) up to max. 2 W/kg in normal operating mode and at max. 15 min continuous scan duration

Implants of the LEANDER® Cervical Vertebral Body Replacement System can heat up by up to 5.4 °C with continuous scanning time of 15 minutes and under the scanning conditions defined above. Under these conditions, a patient can be examined at low risk. To keep the risks of warming low, the scan duration should be kept as short as possible and the SAR as low as possible.

The product has not been examined for artifacts in MR environment. MR image quality may be reduced if the ex-

aminated area is in relative proximity to LEANDER® system implants.

A thorough risk-benefit assessment by the attending doctor must be carried out.

Preoperative

- The use of implants for implantation must be decided upon in accordance with the surgical and medical indications, the potential risks and restrictions associated with this type of surgical procedure, as well as in knowledge of the contraindications, side effects and defined precautions, and in awareness of both the nature and the physical, chemical and biological characteristics of the implant.
- A decision for an (elective) implantation must be made in agreement between the physician and the patient taking into account the general health condition including comorbidities and the level of activity of the patient. The physician has to inform the patient about the consequences of the implantation, especially that the treatment results in segmental stiffness, and about the possibility that the desired surgical result cannot be achieved.
- The patient has to be informed about any possible side effects of the implants and the risks of the procedure, including any revisions or further surgeries that may be necessary.
- The physician shall inform the patient that several causes, such as e.g. adjacent segment degeneration, bone resorption, diffuse spinal tumours or infections, can lead to the reduction of structural support of the cage in the intervertebral space, increasing the risk of subsidence, implant loosening and migration and limit or make impossible correction and stabilisation
- To exclude a suspected allergy against the implant materials, appropriate tests shall be carried out
- Particularly in older patients, the user must weigh the benefits against the risks, e.g. osteoporosis, slower healing and increased comorbidities such as heart disease or other age-related conditions that carry a higher risk of surgical complications such as infections or cardiovascular problems
- The use of LEANDER® implants during pregnancy should be avoided whenever possible, as it poses significant risks to both mother and fetus. These risks include infections, thromboembolic events, implant-related complications, and potential harm to fetal development. Physicians are obligated to provide comprehensive counseling, carefully assess conservative treatment alternatives, and implement strategies to minimize blood loss, radiation exposure, and anesthesia-related risks. Careful planning and close monitoring throughout the preoperative, intraoperative, and postoperative phases are essential to avoid complications and to ensure the best possible outcomes for both mother and child
- Polypharmacy, especially hyperpolypharmacy, can

have a negative impact on the surgical outcome and the patient's health in elderly patients. The number of prescribed medications in this patient group should be considered prior to surgery

- The activities / movement behaviour of the patient has a significant influence on the service life of the implant. The patient must be informed that any activity increases the risk of loss, bending or fracture of the implant components. It is crucial to inform the patient about restrictions in activities in the post-operative phase and to monitor the patient post-operatively in order to assess the development of the fusion and the condition of the implant. Even when a solid bone fusion has occurred, implant components may still bend, break or loosen. Therefore, the patient must be informed that implant components can also bend, break or loosen if the restrictions in activities are properly followed
- The treating physician is responsible for the accurate diagnosis and decision-making for or against surgery as well as the appropriate planning and conductance of the surgery. This implies the patient-specific choice of size, endplate angulation and design of the implants, taking into account the individual pathology and anatomy
- The correct selection and placement of implants has to be ensured and checked based on suitable pre-, intra- and postoperative diagnostic procedures
- Errors in implant selection can lead to premature clinical implant failure. The number of segments to be treated has to be determined carefully. The surgeon has to ensure the availability of sufficient different implant sizes and the required instruments for the surgical procedure
- HumanTech Spine accepts no liability for negative consequences resulting from incorrect diagnosis or indication, an inappropriate choice of implant, incorrect use of implant components and/or an improper surgical technique or asepsis

Intraoperative

- The instructions in the surgical technique must be adhered to. You can download these at www.humantech-spine.de or obtain them directly from your HumanTech representative. It is mandatory to use the instruments described in the surgical technique, which are intended to be used with the system. No liability can be assumed for the use of external/third-party instruments
- Implants showing any signs of damage or contamination must not be implanted
- Components of the LEANDER® system may not be replaced by components / products from other systems of another source or different manufacturer. Furthermore, no direct connection of components / products of the system to components of other systems may be established. If this is not complied with, or if the products are improperly used or applied, HumanTech Spine GmbH assumes no responsibility

- Breakage, slippage or incorrect use of the instruments or implants can result in an injury of the patient or the surgical personnel or result in an extension of the duration of the surgery
- Care should be taken to maintain the integrity of the endplates of the vertebral bodies. Damage to the endplates of the vertebral bodies or excessive partial abrasion of the endplates of the vertebral bodies can lead to sintering of the implant and to a loss of segmental stability
- The endplates of the implant should be selected in a sufficiently large size in order to make use of the maximum contact surface of the vertebral body. This will ensure the greatest possible stability and counteract implant sintering. However, in order to avoid injury to the adjacent structures, the endplates of the implant should not protrude beyond the vertebral body
- Implantation of the implant must be performed under x-ray control, as indicated in the surgical technique, in order to avoid misplacement
- The endplates of the implant must be assembled correctly to the implant base body
- Careful inspection of the implants is recommended before and during surgery to ensure that the implants have been assembled correctly according to the instructions in the surgical technique and that they remain undamaged during assembly, previous procedures and insertion
- The size of the implant components must be selected so that the implant is under appropriate tension once the operation is completed
- Proceed with extreme caution in the area of the spinal cord and the nerve roots, as any damage to nerves can result in the negative impairment or loss of neurological functions
- Extreme caution is required to avoid injuries of the aorta, carotid artery, or jugular vein because this can be life-threatening if not treated immediately
- Should an implant component or parts of an instrument loosen or break in situ, these residues must be removed completely
- Implants that have come into contact with a patient's body fluids or tissues or have been soiled must not be reused
- Bone cement may not be used as it makes the removal of components difficult or even impossible, particularly in revision surgeries
- A prolonged healing phase, lack of bone fusion or subsequent bone resorption can lead to overloading of the implant, which may result in implant failure
- The implants are used to reconstruct substantial defects in the surgical area during a healing process lasting a maximum of two years. After spinal fusion has been achieved, the implants are part of the newly formed body structure. Explantation of the implants is not intended unless complications, implant failure or a lack of fusion occur, which can require implant removal. In case of complications, the surgeon has to decide whether a revision of the implant should be carried out, taking into account the patient's condition and the possible risks involved, e.g. that the removal of an implant already fused with the adjacent vertebral bodies is associated with considerable trauma
- The implants are intended for single-use only. Never reuse the implants. Even if the implant appears to be intact after the revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break
- Handle removed implants in such a way that their reuse is not possible
- An implant card must be given to the patient

MATERIAL

All implants of the LEANDER® Cervical Vertebral Body Replacement System are made out of the titanium alloy Ti6Al4V according to ISO 5832-3.

PACKAGING, LABELING, TRANSPORT AND STORAGE

- The devices must be handled, transported and stored with care. Damages to the device packaging or damages to the device itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increased internal loads, which can result in a fracture of the device.
- The implants and corresponding instruments should be stored at room temperature. Environmental factors such as salty air, humidity, direct sunlight, chemicals, etc., must not be allowed to act on the devices.
- Before the surgery, a careful inspection of the components of the LEANDER® system to be used must be carried out in order to avoid damages caused by storage, transport or previous procedures.

The implants are delivered in STERILE packaged condition:

All implants are sterilized by a validated gamma sterilization procedure and labelled as STERILE. Cleaning, disinfection and sterilization prior to use need not be performed. The implants are supplied packaged in a blister-blister-combination and protected by an outer carton. The devices

Postoperative

- The patient should be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area
- The patient should be monitored if a change in the area of the implant(s) is noticed. The treating physician should assess the possibility of clinical implant failure and discuss with the patient the actions required to achieve further healing

may only be used when the label on the outer packaging and also the inner packaging are intact. If the packaging is damaged or already open, the sterility of the device is not guaranteed, and the device must not be used.

The implants must not be used if the indicated expiry date has been exceeded. Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or if the sterile packaging is damaged is not intended.

HumanTech Spine GmbH assumes no responsibility for the use of re-sterilized devices regardless of the person who performed the resterilization or the method used.

The rules of asepsis must be observed when removing the devices from the sterile packaging. The sterile packaging may only be opened immediately prior to the insertion of the device. It is recommended that a replacement device is always kept available. The implant must be removed from the packaging using appropriate aseptic precautions.

CLEANING, DISINFECTION AND STERILISATION

UNSTERILE delivered instruments of the LEANDER® System must be cleaned, disinfected and sterilized prior to use. All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the according processing instruction.

The latest version can always be found on our homepage: "<https://humantech-spine.de/378-de-IFUs.html>", as well as on request from HumanTech Spine GmbH.

DISPOSAL

Any HumanTech implant that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol.

Devices must be disposed of as a healthcare medical device in accordance with hospital procedures.

PRODUCT COMPLAINTS

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a LEANDER® implant ever exhibits "malfunctioning" (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, must be reported immediately to the manufacturer and the competent authority of the respective country by phone, fax, or in writing.

Should you have a complaint, please provide the name,

item number and batch number of the component as well as your name, address, and detailed description of the fault as possible in written form.

Information to be provided to the patient

- The patient should be informed about contraindications, warnings, precautions, side effects and complications with HumanTech devices
- The patient must be informed about MRI compatibility regarding the HumanTech products used
- The patient should be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area
- The patient has to be informed that the implant card must be shown to the treating physician prior to every planned examination, diagnostic intervention and surgical or non-surgical treatment
- An implant card must be given to the patient

FURTHER INFORMATION

In the event of complaints, suggestions, or indications regarding the content of these instructions for use or the use of the product, please refer to address listed on the last page.

After the European Medical Devices Database is online, the summary of safety and clinical performance of LEANDER® implants will be available at the following address „<https://ec.europa.eu/tools/eudamed>“. Until then, you can find them on our homepage "<https://www.humantech-spine.de>", as well as on request from HumanTech Spine GmbH.

If required, the instructions for use can be ordered in paper form from the address given on the last page. It will be sent by post within seven days.

© HUMANTECH Spine GmbH. All rights reserved.

EXPLANATION OF SYMBOLS

	Manufacturer
	Date of manufacture
	Do not re-use
	Lot number
	Catalogue number
	Quantity
	Medical device
	Unique Device Identifier
	Consult instructions for use
	Keep dry
	Sterilized using irradiation
	Single sterile barrier system with protective packaging inside
	Do not re-sterilize
	Use-by date
	Do not use if package is damaged and consult instructions for use
	Caution
	MR conditional

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