TRISTAN[®]

CERVICAL INTERBODY FUSION SYSTEM

Gebrauchsanweisung	DE	2 - 10
Instruction for use	EN	11 - 18
Návod k použití	CS	
Manual de instrucciones	ES	
οδηγίες χρήσης	EL	
Használati utasítás	HU	
Istruzioni per l'uso	IT	
Lietošanas instrukcija	LV	
Instructies voor gebruik	NL	
Instruções de utilização	PT	
Instrucțiuni de utilizare	RO	
Návod na použitie	SK	

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IMPORTANT INFORMATION

This instruction for use and further product-specific information (e.g. surgical technique described in the brochure and separate IFU for instrumentation) must be read carefully and its content must be adhered to. For current product-specific information, refer to:

www.humantech-spine.de

This instruction for use provides essential information but is not sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a clearly defined postoperative treatment regimen, and to conduct scheduled follow-up examinations.

For improved readability, the tradename TRISTAN® Cervical Interbody Fusion System may be shortened to TRISTAN® system throughout the IFU.

SCOPE

The instruction for use is applicable for all sterile and unsterile delivered implants of the TRISTAN® Cervical Interbody Fusion System.

INTENDED USE

The implants of the TRISTAN® Cervical Interbody Fusion System are intended to be used in an ACDF surgery (Anterior Cervical Discectomy and Fusion). They replace an intervertebral disc in the cervical spine (between C3 and C7) following partial or total anterior cervical discectomy. The implants facilitate bony fusion of the adjacent vertebral bodies by providing structural support during the natural healing process, thereby contributing to the stabilization of the cervical spine. The implants are designed to restore and maintain the original disc height and to correct the spinal profile.

Implants of the TRISTAN® Cervical Interbody Fusion System are intended for stand-alone use, with mono- and multi-segmental treatments possible.

For multi-segmental fusions or in cases of low bone quality, additional stabilization with an additional anterior cervical plate screw fixation system (e.g. HERO®) is strongly recommended.

GENERAL PRODUCT DESCRIPTION

The TRISTAN® Cervical Interbody Fusion System is intended to be implanted after partial or total anterior

cervical discectomy in an open surgical approach. The wedge-formed cages are designed to provide structural support between C3 and C7 during bony fusion of the two adjacent vertebral bodies with the implant and to allow long-term stability of the affected spinal segment. Large filling apertures in the cage are typically filled with autologous bone or bone substitute materials to promote bone growth and healing (fusion).

The use of the TRISTAN® Cervical Interbody Fusion System should lead to clinical benefit for the patient by improving its quality of life through pain reduction, clinically relevant improvement in function and stabilization of the spine.

Different cage designs are available:

- TRISTAN® Ti
- TRISTAN® PEEK-S
- TRISTAN® R-PEEK-Ti Coated
- TRISTAN®flex

The devices are delivered in both unsterile and sterile packaged condition.

Different sizes and angulations are available to meet the individual natural curvature of the spine. The base bodies of the cages consist of titanium alloy or PEEK; the Ti-coating facilitates the osseointegration of the cages. Spikes are designed to allow additional anchoring of TRISTAN® PEEK-S and TRISTAN®flex cages in the caudal and cranial vertebral body cover plates. After discectomy, an implant of appropriate size and angulation is selected and inserted, followed by compression of the adjacent vertebral bodies. The compression ensures the proper penetration of spikes into the vertebral endplates.

The expansion mechanism of TRISTAN®flex provides adjustment to different angulations, 5° (not expanded), 7.5° (middle position) and 10° (fully expanded) facilitating adaptation to the natural spinal curvature without changing the cage.

The HumanTech instruments and accessories, intended for use with the implants of the TRISTAN® Cervical Interbody Fusion System, and the surgical techniques for the specific applications are described in detail in the corresponding brochures and are mandatory for implantation.

INDICATIONS FOR USE

A treatment of the following diseases and injuries of the cervical spine with implants of the TRISTAN® Cervical Interbody Fusion System is only intended if Anterior Cervical Discectomy and Fusion is indicated:

- Instabilities
- Cervical Degenerative Disc Disease
- Disc herniation
- Spinal stenosis (spinal canal stenosis or neuroforaminal stenosis)
- Myelopathy and/or radiculopathy
- Ossification of the posterior longitudinal ligament

- Trauma
- Revision surgery
- Deformities

CONTRAINDICATIONS

Contraindications may either be relative or absolute:

Absolute contraindications of TRISTAN®

- Suspected or documented allergy or intolerance to the materials used
- Untreated local infection in the operative field or in the cervical spine
- Systemic sepsis or uncontrolled infection
- Severe, decompensated medical illness
- Medical or surgical condition which would preclude a potential benefit of spinal implant surgery, such as fracture near the surgical site or extreme anatomic deformities expected to interfere with the usage of the implant, to impair the instrumentation or to restrict physical performance of the patient (beyond the consequences of fusion)
- Inadequate bone stock or bone quality making implant anchoring impossible
- Conditions which make an anterior approach impossible
- Pathologies which require corpectomy with vertebral body replacement and other stabilizing measures for the reconstruction and stabilization of the segment
- If the implant components selected for use are too large or too small, or the selected angulation is insufficient to achieve a satisfactory result
- If fusion is not required or already exists in the affected segment

Relative contraindications of TRISTAN®

Relative contraindications are conditions that may decrease the surgical success and its outcomes.

The responsible physician has to weigh the risks and benefits of an implantation in each individual case. If the implantation is performed despite the presence of a relative contraindication, special care or procedures may be required.

- Anterior cervical discectomy and fusion (ACDF) in the segment C7/T1
- Fever or leucocytosis
- Chronic infection
- Acute infection which is not in the operative field
- Significant risk of infection (e.g. in case immunocompromising chronic infections)
- Signs of local inflammation
- Local bone tumors
- Pre-operative dysphagia
- Pre-operative airway compromise (e.g. COPD, chronic obstructive pulmonary disease)
- Open wounds

- Inadequate cervical soft tissue structure or coverage (e.g. poor skin condition) at the surgical site
- Pregnancy: Operations on pregnant women must be avoided if possible. If they are nevertheless performed, they require special care or procedures
- If the use of other components, metals or alloys available in this system is required
- Concomitant diseases or conditions that could negatively affect the function and success of the implant, e.g.
 - advanced age
 - smoking
 - diabetes
 - osteoporosis, osteopenia, osteomalacia, joint disease, bone resorption
 - morbid obesity
 - chronic kidney disease
 - neuromuscular diseases that would place excess strain on the implant during the healing period
 - psychological disorders or other conditions of reduced mental health
 - drug or alcohol abuse
 - wound healing disorders
- If the patient is unwilling or able to follow postoperative instructions

INTENDED PATIENT GROUP

All patients that are at least 18 years old and skeletally matured at the time of surgery.

INTENDED USERS

The devices of the TRISTAN® Cervical Interbody Fusion System may only be implanted by surgeons who have completed the necessary training in spinal surgery.

The

- professional judgement and/or clinical skills and experience of the surgeon with regard to careful patient selection
- preoperative planning and implant selection
- knowledge of the anatomy and biomechanics of the spine
- understanding of the material and the mechanical properties of the used implants
- training and skills in spinal surgery and the use of the instruments required for inserting the implants
- surgeon's ability to gain the patient's consent to adhere to a clearly defined post-operative treatment regimen and to conduct scheduled follow-up examinations

is/are necessary for the use of the devices of the system. The surgeons are supported in the operating theater environment by surgical assistants. General knowledge in surgical procedures, the surgical technique and special training and skills regarding the correct use of the devices of the system and the use of the instruments required for inserting the implants is/are necessary for the use of the devices.

Unsterile delivered products must be cleaned, disinfected and sterilized before use on the patient by persons responsible for processing procedures.

INTENDED ENVIRONMENT

The devices have to be used in an operating theater.

A sterile environment and compliance with the general hygiene regulations that apply in operating theaters are basic requirements for using the devices.

It must be ensured that the surgical field is fully illuminated with a view of the situs.

The implantation of the product in the patient takes place under general anesthesia. In order to be able to perform the surgery under general anesthesia, anesthetic gas or other anesthetics are required.

It has to be assured that C-Arm (X-ray) and the system specific instruments are available. Furthermore, general surgical instruments are required for the performance of the surgery.

SIDE EFFECTS / ADVERSE EVENTS

The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

Possibly related to the implant:

- Unintended modification of spinal curvature and stiffness of the vertebral column
- Delayed bone growth/bone healing/osseointegration, non-union and pseudarthrosis
- Partial loss of the degree of correction achieved during the surgery
- Subsidence
- Migration
- Reduced bone density due to stress shielding
- Hardware failure, e.g. breakage, deformation of the implant, loosening of implant or implant components
- Proximal Junctional Kyphosis (adjacent segment degeneration)
- Vertebral joint degeneration due to altered biomechanics in the spine
- Local pain, discomfort or abnormal sensitivity
- Allergic reaction to the implant material
- Risk of infection or inflammation
- Autoimmune reaction
- Metallosis due to intra- or postoperative abrasion or corrosion with possible secondary chronic-inflammatory foreign body reactions as scar formation, local tissue necrosis, formation of neoplasia/tumours and/or osteolysis
- Long-term dysphagia due to implant pressure
- Revision surgery or requirement of further surgeries
- Due to the permanent presence of an implant/implants in the body, there are postoperative risks such as fracture, microfracture, resorption, necrosis, damage to, penetration or subluxation of a vertebral body above or below the treated segment(s)

Possibly related to the surgical procedure:

- Malpositioning of the implant
- Superficial infection (wound), deep site infection and inflammatory phenomena
- Dysphagia (difficulty in swallowing)
- Hoarseness (dysphonia)
- Airway compromise
- Arm or neck cramps
- Loss or increase of spinal mobility
- Pain, persistent pain
- Swelling and edema
- Herniated disc, disc destruction or disc degeneration in the intervertebral space above or below the treated segment(s)
- Complications at the bone transplant donor site
- Neurological complications or deterioration, e.g.:
 - Transient or permanent motor weakness, loss or change of sensitivity or functions, reflex changes, numbness, pain, or paralysis
 - Due to direct injury or irritation of the spinal cord, nerve roots (e.g. radiculopathy), or peripheral nerves
 - Secondary effects such as scarring, dural injury (durotomy), or postoperative inflammation
 - Injury or lesion of the dura mater with potential cerebrospinal fluid (CSF) loss or formation of a CSF fistula
- Neurological or neurogenic functional disorders, e.g.:
 - Urinary retention, loss of bladder control, or other urological impairments
 - Dysfunction of the rectum or gastrointestinal system
 - Sexual dysfunction
 - Damage to the reproductive system or infertility
- Heterotopic ossification of soft tissue
- Injury of or pressure on surrounding bone structures, soft tissue, nerves or organs (e.g. lung, trachea, esophagus)
- Injury or lesion of blood vessels (including aorta, carotid artery, or jugular vein) or lymphatic vessels (including potential lymph leakage)

Possible general surgical risk:

- Delayed wound healing
- Superficial wound dehiscence
- Infection
- Discomfort
- Inability to perform daily tasks
- Bleeding and/or haematomas
- Vascular disorders including e.g. thrombus, deep vein thrombosis, thrombophlebitis, arrhythmias, high blood pressure and blood pressure fluctuations, myocardial infarction, stroke or death
- Development of respiratory problems, including e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, pneumothorax, transient pleural effusion or death

- Temporary kidney dysfunction
- All risks associated with anesthesia
- All risks associated with blood transfusion

WARNINGS AND PRECAUTIONS

General

- When using X-ray equipment for diagnosis or control, the national limit values for the radiation dose must be observed
- The TRISTAN® System consist of implants and instruments. TRISTAN® implants must only be used with the intended system-specific instruments and accessories described in the surgical techniques in the respective brochures

MR safety information

Implants of the TRISTAN® Cervical Interbody Fusion System have been examined in non-clinical-testings in accordance with ASTM norms and determined as “MR conditional”. Patients with TRISTAN® implants can be safely examined via MRI on following conditions:

- Static magnetic field of 1.5 T to 3 T
- Spatial gradient magnetic field up to max. 1100 Gauss/cm or 11 T/m
- Whole body specific absorption rate (SAR) up to max. 2 W/kg in normal operating mode and at max. 15 min continuous scan duration

Implants of the TRISTAN® Cervical Interbody Fusion System can heat up by up to 3.5 °C with continuous scanning time of 15 minutes and under the scanning conditions defined above. Under these conditions, a patient can be examined at low risk. To keep the risks of warming low, the scan duration should be kept as short as possible and the SAR as low as possible.

The product has not been examined for artifacts in MR environment. MR image quality may be reduced if the examined area is in relative proximity to TRISTAN® system implants.

A thorough risk-benefit assessment by the attending doctor must be carried out.

Preoperative

- The use of implants for implantation must be decided upon in accordance with the surgical and medical indications, the potential risks and restrictions associated with this type of surgical procedure, as well as in knowledge of the contraindications, side effects and defined precautions, and in awareness of both the nature and the physical, chemical and biological characteristics of the implant
- A decision for an (elective) implantation must be made in agreement between the physician and the patient taking into account the general health condition including comorbidities and the level of activity of the pa-

tient. The physician has to inform the patient about the consequences of the implantation, especially that the treatment results in segmental stiffness, and about the possibility that the desired surgical result cannot be achieved

- The patient has to be informed about any possible side effects of the implants and the risks of the procedure, including any revisions or further surgeries that may be necessary
- The physician shall inform the patient that several causes, such as e.g. adjacent segment degeneration, bone resorption, diffuse spinal tumours or infections, can lead to the reduction of structural support of the cage in the intervertebral space, increasing the risk of subsidence, implant loosening and migration and limit or make impossible correction and stabilisation
- To exclude a suspected allergy against the implant materials, appropriate tests shall be carried out
- For elderly patients, the responsible physician must carefully weigh the potentially benefits against the risks, taking into account the patient’s general health condition, increased comorbidities such as e.g. osteoporosis, slower healing and an age-related higher risk of general surgical complications such as infections or cardiovascular problems
- The use of TRISTAN® implants during pregnancy should be avoided whenever possible, as it poses significant risks to both mother and fetus. These risks include infections, thromboembolic events, implant-related complications, and potential harm to fetal development. Physicians are obligated to provide comprehensive counseling, carefully assess conservative treatment alternatives, and implement strategies to minimize blood loss, radiation exposure, and anesthesia-related risks. Careful planning and close monitoring throughout the preoperative, intraoperative, and postoperative phases are essential to avoid complications and to ensure the best possible outcomes for both mother and child
- Polypharmacy, especially hyperpolypharmacy, can have a negative impact on the surgical outcome and the patient’s health in elderly patients. The number and interactions of prescribed medications in this patient group should be considered prior to surgery
- The activities / movement behavior of the patient has a significant influence on the service life of the implant. The patient must be informed that any activity increases the risk of loss, bending or fracture of the implant components. It is crucial to inform the patient about restrictions in activities in the post-operative phase and to monitor the patient post-operatively in order to assess the development of the fusion and the condition of the implant. Even when bone fusion is completed, implant components may still bend, break or loosen. The patient must be informed that implant components can also bend, break or loosen if the restrictions in activities are properly followed
- The treating physician is responsible for the accurate

diagnosis and decision-making for or against surgery as well as the appropriate planning and conductance of the surgery. This implies the patient-specific choice of size, angulation and design of the implants, taking into account the individual pathology and anatomy

- The correct selection and placement of implants has to be ensured and checked based on suitable pre-, intra- and postoperative diagnostic procedures
- Errors in implant selection can lead to premature clinical implant failure. The number of segments to be treated has to be determined carefully. The surgeon has to ensure the availability of sufficient different implant sizes and the required instruments for the surgical procedure
- The physician shall inform the patient that, if bone quality is low or cannot be assessed pre- or intraoperatively or in case of multisegmental treatments, additional stabilization with an anterior cervical plate-screw fixation system (e.g. HERO®) may be necessary and can carry additional risks such as e.g. dysphagia. Risks as e.g. dysphagia or pseudarthrosis may increase with the number of treated segments. The treatment of more than 3 segments is associated with significantly higher risk of non-fusion and other complications and should only be performed after careful risk-benefit consideration
- HumanTech Spine accepts no liability for negative consequences resulting from incorrect diagnosis or indication, an inappropriate choice of implant, incorrect use of implant components and/or an improper surgical technique or a sepsis
- In treatments of more than two levels, additional posterior instrumentation may be necessary to decrease the complication rate

Intraoperative

- The instructions in the surgical technique must be adhered to. You can download these online at www.humantech-spine.de or obtain them directly from your HumanTech representative. It is mandatory to use the instruments described in the surgical technique, which are intended to be used with the system. No liability can be assumed for the use of third-party instruments
- Implants showing any signs of damage or contamination must not be implanted
- Careful inspection of the implants is recommended before and during surgery to ensure that the implants have been assembled correctly according to the instructions in the surgical technique and that they remain undamaged during assembly, previous procedures and insertion
- Components of the TRISTAN® System must not be replaced by components / products from other systems from another source or from a different manufacturer. Furthermore, a direct connection of components / products of the TRISTAN® System to components of other systems must not be established. If this is not obser-

ved or if the products are used improperly or otherwise than according to IFU and the mandatory surgical technique, HumanTech Spine GmbH assumes no responsibility

- Breakage, slippage or incorrect use of the instruments or implants can result in injuries of the patient or the surgical personnel or in prolonged surgery duration
- Care should be taken to maintain the integrity of the endplates of the vertebral bodies. Damage to the endplates of the vertebral bodies or excessive partial abrasion of the endplates of the vertebral bodies can lead to subsidence of the implant and to a loss of segmental stability
- The sizes of the implant should be selected in a sufficiently large size in order to make use of the maximum contact surface of the vertebral body. This will ensure the greatest possible stability and counteract implant sintering
- Implantation of the implant must be performed under x-ray control, as indicated in the surgical technique, in order to avoid misplacement. Proceed with extreme caution in the area of the spinal cord and the nerve roots, as any damage to nerves can result in the negative impairment or loss of neurological functions
- Extreme caution is required to avoid injuries of the aorta, carotid artery, or jugular vein because this can be life-threatening if not treated immediately
- Should an implant component or parts of an instrument loosen or break in situ, these residues must be removed completely
- Implants that have already come into contact with a patient's body fluids or tissues or have been soiled must not be reused
- Bone cement must not be used, as the removal of the implant becomes difficult or impossible. Moreover, the heat produced by the hardening process can damage or deform the PEEK implants
- When using PEEK implants particular care in preparing the intervertebral disc space is needed to keep the applied forces to the implant during insertion as low as possible
- When inserting the PEEK-Ti coated implants, special care is required as the coating for the improvement of the ingrowth behavior consists of a structured, rough surface
- To achieve a good primary fixation, the TRISTAN®flex and TRISTAN® PEEK-S cages have spiky and sharp spikes, therefore they must be handled carefully to avoid injuries of patient, users and third parties. Moreover, sufficient distraction of the intervertebral space is necessary prior to the insertion of the implant

Postoperative

- The patient should be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area
- The patient should be monitored if a change in the

area of the implant(s) is noticed. The treating physician should assess the possibility of clinical implant failure and discuss with the patient the actions required to achieve further healing

- A prolonged healing phase, lack of bone fusion or subsequent bone resorption can lead to overloading of the implant, which may result in implant failure
- The implants are used to reconstruct substantial defects in the surgical area during a healing process lasting a maximum of two years. After spinal fusion has been achieved, the implants are part of the newly formed body structure. Explantation of the implants is not intended unless complications, implant failure or a lack of fusion occur, which can require implant removal. In case of complications, the surgeon has to decide whether a revision of the implant should be carried out, taking into account the patient's condition and the possible risks involved, e.g. is the removal of the implant already fused with the adjacent vertebral bodies associated with considerable trauma
- The implants are intended for single-use only. Never reuse the implants. Even if the implant appears to be intact after the revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break
- Handle removed implants in such a way that their reuse is not possible
- An implant card must be given to the patient

MATERIAL

There are different implant product lines.

Product line	Material
TRISTAN® Ti	Base body: Ti6Al4V according to ISO 5832-3
TRISTAN® PEEK-S	Base body: PEEK Optima® according to ASTM F2026 Spikes: Ti6Al4V according to ISO 5832-3 X-Ray markers: Tantalum according to ASTM F560
TRISTAN® R-PEEK-Ti Coated	Base body: PEEK Optima® according to ASTM F2026 Coating: cpTi (commercially pure titanium) X-Ray markers: Tantalum according to ASTM F560
TRISTAN®flex	Base body and additional comp.: PEEK Optima® according to ASTM F2026 Certain comp. (e.g., spikes): Ti6Al4V according to ISO 5832-3 X-Ray markers: Tantalum according to ASTM F560

PACKAGING, LABELING, TRANSPORT AND STORAGE

- The devices must be handled, transported and stored with care. Damages to the device packaging or damages to the device itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increased internal loads, which can result in a fracture of the device.
- The implants and corresponding instruments should be stored at room temperature. Environmental factors such as salty air, humidity, direct sunlight, chemicals, etc., must not be allowed to act on the devices.
- Before the surgery, a careful inspection of the components of the TRISTAN® System to be used must be carried out in order to avoid damages caused by storage, transport or previous procedures. The sterilization containers, trays and the associated covers must not be damaged.

The devices are delivered in both NON-STERILE and STERILE packaged condition:

Unsterile delivered devices:

UNSTERILE delivered devices are labelled as UNSTERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILIZATION). The devices are supplied in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

Sterile delivered devices:

STERILE delivered devices are sterilized by a validated gamma or ethylene oxide sterilization procedure and labelled as STERILE. Cleaning, disinfection and sterilization prior to use need not be performed.

The implants are supplied individually packaged in a blister-pouch- or blister-blister-combination and protected by an outer carton. The devices may only be used when the label on the outer packaging and also the inner packaging is intact. If the packaging is damaged or already open, the sterility of the device is not guaranteed, and the device must not be used.

The devices must not be used if the indicated expiry date has been exceeded. Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or if the sterile packaging is damaged is not intended.

HumanTech Spine GmbH assumes no responsibility for the use of resterilized devices regardless of the person who performed the resterilization or the method used. The rules of asepsis must be observed when removing the devices from the sterile packaging. The sterile packaging

may only be opened immediately prior to the insertion of the device. It is recommended that a replacement device is always kept available. The implant must be removed from the packaging using appropriate aseptic precautions.

CLEANING, DISINFECTION AND STERILISATION

UNSTERILE packaged delivered products of the TRISTAN® System must be cleaned, disinfected and sterilized prior to use. For a correct execution of the processing procedure, general hygienic principles regarding the processing of medical devices and the knowledge of the defined steps of the processing instruction are necessary.

All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the according processing instruction. The latest version can always be found on our homepage: "<https://humantech-spine.de/378-de-IFUs.html>", as well as on request from HumanTech Spine GmbH.

DISPOSAL

Any HumanTech implant that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol.

Devices must be disposed of as a healthcare medical device in accordance with hospital procedures.

PRODUCT COMPLAINTS

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a product of the TRISTAN® implant system ever exhibits "malfunctioning" (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, must be reported immediately to the manufacturer and the competent authority of the respective country by phone, fax, or in writing.

Should you have a complaint, please provide the name, item number and batch number of the component as well as your name, address, and detailed description of the fault as possible in written form.

Information to be provided to the patient

- The patient should be informed about contraindications, warnings, precautions, side effects and complications with HumanTech devices
- The patient must be informed about MRI compatibility

regarding the HumanTech products used

- The patient should be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area
- An implant card must be given to the patient
- The patient has to be informed that the implant card must be shown to the treating physician prior to every planned examination, diagnostic intervention and surgical or non-surgical treatment

FURTHER INFORMATION

In the event of complaints, suggestions, or indications regarding the content of these instructions for use or the use of the product, please refer to address listed on the last page.

After the European Medical Devices Database is online, the summary of safety and clinical performance of TRISTAN® implants will be available at the following address "<https://ec.europa.eu/tools/eudamed>". Until then, you can find them on our homepage "<https://www.humantech-spine.de>", as well as on request from HumanTech Spine GmbH. If required, the instructions for use can be ordered in paper form from the address given on the last page. It will be sent by post within seven days.

© HUMANTECH Spine GmbH. All rights reserved.

EXPLANATION OF SYMBOLS

	Manufacturer
	Date of manufacture
	Do not re-use
	Lot number
	Catalogue number
	Quantity
	Medical device
	Unique Device Identifier
	Consult instructions for use
	Keep dry
	Non-sterile
	Sterilized using irradiation
	Sterilized using ethylene oxide
	Single sterile barrier system with protective packaging inside
	Single sterile barrier system with protective packaging outside
	Do not re-sterilize
	Use-by date
	Do not use if package is damaged and consult instructions for use
	Attention
	MR conditional

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