VENUS[®]

SPINAL FIXATION SYSTEM

Gebrauchsanweisung	DE	2 - 10
Instruction for use	EN	11 - 17
Instrucciones de uso	ES	
Instruções de Uso	PT-BR	
Használati utasítás	HU	
Návod na použitie	SK	
Инструкция по применению	RU	
οδηγίες χρήσης	EL	
Istruzioni per l'uso	IT	
Instrucțiuni de utilizare	RO	
Instructies voor gebruik	NL	
Návod k použití	CS	
Lietošanas instrukcija	LV	

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IMPORTANT INFORMATION

This instruction for use and further product-specific information (e.g. surgical technique described in the brochure and separate IFU for instrumentation) must be read carefully and its content must be adhered to. For actual product-specific information, refer to:

www.humantech-spine.de

This instruction for use provides essential information, but is not sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a clearly defined post-operative treatment regimen, and to conduct scheduled follow-up examinations. For improved readability, the tradename VENUS™ Spinal Fixation System may be shortened to VENUS™ system throughout the IFU.

SCOPE

The instruction for use is applicable for all sterile and unsterile delivered implants of the VENUS™ Spinal Fixation System.

INTENDED USE

The implants of the VENUS™ Spinal Fixation System are intended for posterior mono- or multisegmental fixation, after prior correction if necessary, and immobilization during bony fusion to achieve stabilization of the thoracic, lumbar and iliosacral spine. The VENUS™ Spinal Fixation System can be used in an open or a minimally invasive surgical approach. The implants are intended for long-term stand-alone use.

GENERAL PRODUCT DESCRIPTION

The VENUS™ Spinal Fixation System includes pedicle screws, rods, set screws, transverse connectors, rod connectors, and extension units to form appropriate spinal screw-rod constructs depending on the application.

The implant system is used for surgical stiffening of spinal segments by fixation. The pedicle screws of the VENUS™ Spinal Fixation System are placed in the left and right pedicles of the vertebrae and then connected to a rod on each side. The essential function of the resulting screw-rod construct is to offer a fixed scaffold that is able to take off strain of the vertebral bodies in the treated segment so that the natural process of bone fusion is as undisturbed as possible. Once the fusion is complete, the implants have no longer any active function but can remain in the body

for the long-term. The use of the VENUS™ Spinal Fixation System should lead to clinical benefit for the patient by improving its quality of life through pain reduction, clinically relevant improvement in function, and/or stabilization to support the natural curvature.

The VENUS™ Spinal Fixation System is supplemented by several kits of additional implants and/or instruments for specific applications, e.g. the extension of an already implanted VENUS™ Spinal Fixation System (VENUS™ Extension), minimally invasive surgery (VENUS™ mini 2.0, VENUS™ mini 2.0 Fracture), the treatment of complicate deformities as scoliosis (VENUS™ Scoliosis), facilitated rod-to-screw approximation under difficult anatomic or surgical conditions by using long-head screws (VENUS™ Reduction) or augmentation in cases of low bone quality (VENUS™ Augmentation). For each application, refer to the corresponding brochure.

The pedicle screws of the VENUS™ Spinal Fixation System includes polyaxial, monoaxial and iliac screws. The screw types are either solid, cannulated, or fenestrated and are available in various forms, different lengths, and diameters to meet the patients' individual pathologies. Polyaxial screws are further available in short and long-head configurations. Long-head screws (Reduction Screws) allow for facilitated approximation of the rods to the screws in difficult anatomic or surgical situations. Cannulated and fenestrated screws of the VENUS™ system can be placed via a guide wire. All cannulated and fenestrated screws (except VENUS™ Reduction screws) are suitable for minimally invasive surgery with the integrated VENUS™ mini set. Fenestrated screws contain a series of fenestrations which allows bone cement (e.g. polymethylmethacrylate (PMMA)) to be injected into the treated site to increase screw fixation in cases of reduced bone quality (augmentation).

VENUS™ Reduction screws cannot be cement augmented with the VENUS™ cement adapter and bone filler.

VENUS™ rods are available in titanium alloy or cobalt chrome. Cobalt chrome offers higher strength than the titanium alloy. Titanium alloy rods are available in straight or pre-curved configurations and in various lengths. To extend an existing construction, Inline Rod Connectors, Domino Connectors, Parallel Connectors or Revisionsconnector Cap can be used. Lateral Connectors can be used for connecting an Ilium Screw or other offset constructions. Further, transverse connector hooks are available which are placed at a 90° angle to the rod construct in order to increase the rotational stability of the entire frame construct.

The VENUS™ Spinal Fixation System can also be used as additional posterior support for interbody fusion systems (e.g. ADONIS™ Lumbar Interbody Fusion System) or vertebral body replacement systems (e.g. SAMSON™ Vertebral Body Replacement System).

The HumanTech instruments and accessories, intended for use with the implants of the VENUS™ Spinal Fixation System, and the surgical techniques for the specific applications are described in the corresponding brochures and are mandatory for implantation.

INDICATIONS FOR USE

A treatment of the following diseases and injuries with VENUS™ Spinal Fixation System is only intended if surgical fixation is indicated.

- instabilities
- degenerative disc disease
- (degenerative) spondylolisthesis
- (degenerative) stenosis (spinal canal stenosis, neuroforaminal stenosis)
- deformities as (idiopathic) scoliosis, kyphosis or lordosis
- pseudarthrosis
- spinal trauma
- tumours
- inflammation (spondylitis, spondylodiscitis)
- conditions requiring revision surgery

CONTRAINDICATIONS

Contraindications may be either relative or absolute.

Absolute contraindications

- Suspected or documented allergy or intolerance to the materials used.
- Medical or surgical condition which would preclude a potential benefit of spinal implant surgery, such as the presence of congenital abnormalities.
- Inadequate tissue coverage over the operative site.
- Inadequate bone stock or bone quality making implant anchoring impossible
- Anatomical structures expected to interfere with the usage of the implant (e.g. extreme anatomic abnormalities) or to impair the instrumentation
- If the implant is expected not being able to perform adequately after implantation.
- If the implant components selected for use would be too large or too small to achieve a successful result.

For minimally invasive/percutaneous procedures only:

- Severe neural deficits, which require canal decompression

Relative contraindications

Relative contraindications are conditions that may decrease the surgical success and its outcomes. The responsible physician has to weigh the risks and benefits of an implantation in each individual case. If the implantation is performed despite the presence of a relative contraindication, special care or procedures may be required.

- Fever or leucocytosis
- Acute or chronic infection
- Significant risk of infection (e.g. in case of immune deficiencies)
- Signs of local inflammation
- Open wounds
- Wound healing disorders

- Pregnancy: Operations on pregnant women must be avoided if possible. If they are nevertheless performed, they require special care or procedures.
- If the use of components other than the metals or alloys used in this system is required.
- Local bone tumors
- Concomitant diseases or conditions that could negatively affect the function and success of the implant, e.g.
 - advanced age
 - smoking
 - diabetes
 - osteoporosis, osteopenia, osteomalacia, joint disease, bone absorption
 - obesity
 - chronic kidney disease
 - depression or conditions of reduced mental health
 - neuromuscular diseases that would place excess strain on the implant during the healing period.
- If the patient is unwilling or not able to follow postoperative instructions.

For minimally invasive/percutaneous procedures only:

- All the situations causing difficulty in imaging such as rotational deformities, obesity and osteopenia
- Multilevel vertebral fractures
- Pedicle fractures
- Adjacent body fracture
- Patient with significant kyphosis or sagittal malalignment

INTENDED PATIENT GROUP

There are no general limitations of the patient population in addition to the defined indications and contraindications. The patient's age may lead to restrictions, e.g. due to immature skeleton in pediatric patients or limited availability of implant sizes. This requires the consideration by an experienced surgeon to assess the suitability of the implant with respect to incomplete growth phase in pediatric patients or infants and other conditions as e.g. the expected activity of the individual patient. Age-related or other comorbidities may limit the surgical success, which requires careful risk-benefit weighting by the responsible surgeon.

INTENDED USERS

The devices of the VENUS™ system may only be implanted by surgeons who have completed the necessary training in spinal surgery.

The

- professional judgement and/or clinical skills and experience of the surgeon with regard to careful patient selection
- preoperative planning and implant selection

- knowledge of the anatomy and biomechanics of the spine
 - understanding of the material and mechanical properties of the used implants
 - training and skills in spinal surgery and the use of the instruments required for inserting the implants
 - surgeon's ability to gain the patient's consent to adhere to a clearly defined postoperative treatment regimen and to conduct scheduled follow-up examinations
- is/are necessary for the use of the devices of the system.

The surgeons are supported in the operating theatre environment by surgical assistants. General knowledge in surgical procedures, the surgical technique and special training and skills regarding the correct use of the devices of the system and the use of the instruments required for inserting the implants is/are necessary for the use of the devices.

Unsterile delivered products must be cleaned, disinfected and sterilized before use on the patient by persons responsible for processing procedures.

INTENDED ENVIRONMENT

The devices have to be used in an operating theater. A sterile environment and compliance with the general hygiene regulations that apply in operating theaters are basic requirements for using the devices.

It must be assured that the surgical field is fully illuminated with a view of the situs.

The implantation of the product in the patient takes place under general anesthesia. In order to be able to perform the operation under general anesthesia, anesthetic gas or other anesthetics are required.

It has to be assured that C-Arm (X-ray) and the systemspecific instruments are available. Furthermore, general surgical instruments are required for the performance of the surgery.

SIDE EFFECTS / ADVERSE EVENTS

The side effects and complications for VENUS™ Spinal Fixation System listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

Possibly related to implant:

- Delayed bone growth or no visible fusion and pseudarthrosis, delayed bone healing or lack of bone healing
- Modification of spinal curvature and stiffness of the vertebral column
- Partial loss of the degree of correction achieved during surgery.
- Reduction of bone density due to stress shielding
- Proximal Junctional Kyphosis (adjacent segment degeneration)
- Altered growth of the fused spine
- Fracture of a vertebra, the pedicle, and/or the sacrum,

- micro fracture, damage to or penetration of a vertebral body above or below the treated segment/s
- Vertebral joint degeneration due to altered biomechanics in the spine
- Allergic reaction to the implant material
- Persistent skin irritations
- Risk of infection or inflammation
- Autoimmune reaction
- Metallosis due to intra- or postoperative abrasion or corrosion with possible secondary chronic-inflammatory foreign body reactions as scar formation, local tissue necrosis, formation of neoplasia/tumours and/or osteolysis
- Breakage, deformation, loosening or migration of the implant
- Local pain, discomfort or abnormal sensitivity
- Pressure sores on the skin caused by components located in positions with insufficient tissue coverage over the implant, with potential penetration of the skin
- Revision surgery
- Stopped growth of the chest wall and lungs, with associated consequences.

Possible risks of augmentation:

- Allergic reaction to the bone cement used
- Cement leakage (possibly leading to neurological impairment, cement embolism, local pain and/or inflammation, revision surgery depending on the anatomical location of the cement flow)

Possibly related to surgical procedure:

- Malposition of the implant
- Superficial infection (wound), deep site infection and inflammatory phenomena
- Neurological complications, e.g.
 - due to injury or irritation of the spinal cord, the root nerves (radiculopathy) or peripheral nerves, e.g. transient or permanent motor weakness, loss of motory or sensory functions or paralysis
 - Injury or lesion of the dura mater with risk of loss of cerebrospinal fluid (CSF) or CSF fistula
 - Urinary retention or loss of bladder control, or other types of impairment of the urological system, dysfunction of the rectum, genitourinary disorders, gastrointestinal disorders
 - Temporary or permanent retrograde ejaculation in men, damage to the reproductive system, sterility and sexual dysfunctions
- Injury or lesion of surrounding soft tissue, blood vessels (e.g. aorta abdominalis, vena cava inferior or other vessels) or lymphatic vessels (including potential lymph leakage)

Possible general surgical risks:

- Delayed wound healing
- Superficial wound dehiscence
- Bleeding and/or haematomas

- Vascular disorders including thrombus, myocardial infarction or death
- Development of respiratory problems, including pulmonary embolism, atelectasis, bronchitis, pneumonia, pneumothorax, transient pleural effusion or death
- Inability to perform daily tasks

WARNINGS AND PRECAUTIONS

General

- When using X-ray equipment for diagnosis or control, the national limit values for the radiation dose must be observed
- VENUS™ system consist of implants and instruments. VENUS™ implants must only be used with the intended system-specific instruments and accessories described in the surgical techniques
- ⚠ VENUS™ Reduction screws are not compatible with the VENUS cement adapter and bone filler and cannot be cement augmented this way!

MR safety information

VENUS™ implants have been examined in non-clinical-testings in accordance with ASTM norms and determined as “MR conditional”. Patients with VENUS™ implants can be safely examined via MRI on following conditions:

- Static magnetic field of 1.5 T to 3 T
- Spatial gradient magnetic field up to max. 1100 Gauss/cm or 11 T/m
- Whole body specific absorption rate (SAR) up to max. 2 W/kg in normal operating mode and at max. 15 min continuous scan duration

VENUS™ implants can heat up to 4.7 °C with continuous scanning time of 15 min. and under the scanning conditions defined above. Under these conditions, a patient can be examined at low risk. To keep the risks of warming low, the scan duration should be kept as short as possible and the SAR as low as possible.

The product has not been examined for artifacts in MR environment. MR image quality may be reduced if the examined area is in relative proximity to VENUS™ implants. A thorough risk-benefit assessment by the attending doctor must be carried out.

Preoperative

- The use of implants for implantation must be decided upon in accordance with the surgical and medical indications, the potential risks and restrictions associated with this type of surgical procedure, as well as in knowledge of the contraindications, side effects and defined precautions, and in awareness of both the nature and the physical, chemical and biological characteristics of the implant.
- A decision for an implantation must be made in

agreement between the doctor and the patient taking into account the general health condition including comorbidities of the patient. The physician has to inform the patient about the consequences of the implantation, especially that the treatment results in segmental stiffness and that the use of the implant in the area of the thoracic spine in children or adolescents can stop the growth of the chest wall and lungs, with associated consequences. As a result, multiple operations must be performed during the growth period in order to adapt the structure to the spine as it grows.

- Moreover, the patient has to be informed, that adapting the structure to a spine that is fully grown can lead to a drastic reduction in the expected service life of the implant. It is therefore advisable to replace the rod and set screw.
- The patient has to be informed in detail about any possible side effects of the implants and the risks of the procedure, including any revisions that may be necessary.
- To exclude a suspected allergy against the implant materials, appropriate tests shall be carried out.
- For elderly patients, the responsible physician must carefully weigh the potentially benefits against the risks, taking into account the patient’s general health condition, increased comorbidities such as e.g. osteoporosis, slower healing and an age-related higher risk of general surgical complications such as infections or cardiovascular problems.
- The use of VENUS™ implants during pregnancy should be avoided whenever possible, as it poses significant risks to both mother and fetus. These risks include infections, thromboembolic events, implant-related complications, and potential harm to fetal development. Physicians are obligated to provide comprehensive counseling, carefully assess conservative treatment alternatives, and implement strategies to minimize blood loss, radiation exposure, and anesthesia-related risks. Careful planning and close monitoring throughout the preoperative, intraoperative, and postoperative phases are essential to avoid complications and to ensure the best possible outcomes for both mother and child.
- Polypharmacy, especially hyperpolypharmacy, can have a negative impact on the surgical outcome and the patient’s health in elderly patients. The number and interactions of prescribed medications in this patient group should be considered prior to surgery.
- The activities / movement behaviour of the patient has a significant influence on the service life of the implant. The patient must be informed that any activity increases the risk of loss, bending or fracture of the implant components. It is crucial to inform the patient about restrictions in activities in the post-operative phase and to monitor the patient post-operatively in order to assess the development of the fusion and the condition of the implant. Even when bone fusion is

completed, implant components may still bend, break or loosen. Therefore, the patient must be informed that implant components can also bend, break or loosen if the restrictions in activities are properly followed.

- Surgical therapy should only be considered after conservative treatment options have failed. The treating physician is responsible for the accurate diagnosis and decision-making for or against surgery as well as the appropriate planning and conductance of the surgery. This implies the patient-specific choice of size and design of the implants, taking into account the individual pathology and anatomy.
- The correct selection and placement of implants is crucial for the surgical success and has to be ensured or to be checked based on suitable pre-, intra- and postoperative diagnostic procedures.
- Errors in implant selection can lead to premature clinical implant failure. The number of segments to be treated is to be determined carefully. The surgeon has to ensure the availability of sufficient different implant sizes and required instruments for the surgical procedure.
- Improper preoperative planning can have a negative impact on the surgical result and the patient's health.
- HumanTech Spine accepts no liability for negative consequences resulting from incorrect diagnosis or indication, an inappropriate choice of implant, incorrect use of implant components and/or an improper surgical technique or asepsis.

Intraoperative

- The instructions in the surgical technique must be adhered to. You can download these online at www.humantech-spine.de or obtain them directly from your HumanTech representative. It is mandatory to use the instruments described in the surgical technique, which are intended to be used with the implant system. No liability can be assumed for the use of third-party instruments.
- Implants showing any signs of damage or contamination must not be implanted.
- Components of the VENUS™ system must not be replaced by components / products from other systems from another source or from a different manufacturer. Furthermore, no direct connection of components / products of the system to components of other systems shall be established. If this is not observed or if the products are used improperly or otherwise than according to IFU and the mandatory surgical technique, HumanTech Spine GmbH assumes no responsibility.
- Bending of the rods / revision connector affect the biomechanical properties of the implant. Bending in the area of the fixation of the rod in the Poly- or Monoaxial Screw can have negative influences on the fixation of the rod – bending in this area has to be avoided.
- Breakage, slippage or incorrect use of the instruments

or implants can result in an injury of the patient or the surgical personnel or result in an extension of the duration of the surgery.

- Care should be taken to maintain the integrity of the pedicle and vertebral body structures. Damage to the pedicle can lead to a loosening of the implants and to a loss of the stability of the fixation.
- The diameter of the screws should be selected in a sufficiently wide size in order to make use of the maximum compression in the pedicle area. This will ensure the greatest possible stability and counteract implant loosening. However, in order to avoid injury to the pedicle structures, the screw width should not fit too tight into the pedicle.
- The rod must be tightened correctly to pedicle screw with the corresponding setscrew.
- Careful inspection of the implants is recommended before and during surgery to ensure that the implants remain undamaged during insertion.
- Implantation of the implant must be performed under continuous x-ray control in order to avoid misplacement
- Proceed with extreme caution in the area of the spinal cord and the nerve roots, as any damage to nerves can result in the negative impairment or loss of neurological functions.
- Extreme caution is required to avoid injuries of the aorta abdominalis or the vena cava inferior because this can be life-threatening.
- If an implant component or part of an instrument loosens or breaks in situ, all resulting residues – whether made of implant material or other materials – must be completely removed.
- Implants that have come into contact with a patient's body fluids or tissues or have been soiled must not be reused.
- The rods of the VENUS™ mini system have marks. It is not allowed to use the area between the ends of the rod and the marks for the fixation of the pedicle screw. In this area the rods must also not be bent. The system stability can be increased with a ventral support. In the lumbar area, especially when the use of pedicle screws with a diameter of 5,5mm or smaller can't be avoided, a ventral support is strongly recommended
- Augmentation: The IFU provided by the manufacturer of the bone cement used as well as the IFU "Bone Filler" and the surgical technique of HumanTech Spine must be adhered to. The volumes of cement used for fenestrated / cannulated screws should ultimately be determined by the surgeon based on the individual patient anatomy. Especially in patients with severe heart or lung impairment, the quantity of applied cement must be chosen as low as possible. Bone cement application has to be monitored by suitable imaging to avoid cement leakage as far as possible.
- For the usage of guide wires in minimally invasive procedures, the instructions in the surgical techniques must be adhered to and all the cautions have to be observed strictly. The correct positioning of the guide

wire must be monitored by suitable imaging during the whole surgical procedure to allow accurate screw insertion and to avoid vessel perforation or other tissue damage.

Postoperative

- The patient should be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area.
- The patient should be monitored if a change in the area where the implant(s) is/are located is noticed. The treating physician should assess the possibility of clinical implant failure and discuss with the patient the actions required to achieve further healing.
- A prolonged healing phase, lack of bone fusion or subsequent bone resorption can lead to overloading of the implant, which may result in implant failure
- The implants are used for internal fixation and immediate immobilization of the treated segments during the bone healing process lasting a maximum of two years. Once the fusion is complete, the implants have no longer any active function. Removal of the implant is possible, but can cause damage to surrounding tissue and loss of stability and should only be performed after a careful risk-benefit consideration by a qualified surgeon and the patient.
- In case of complications, the surgeon has to decide whether a revision of the implant should be carried out, taking into account the patient's health condition and the possible risks involved
- Risk of further injuries due to postoperative trauma
- Never reuse the implants. Even if the implant appears to be intact after the revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been exposed can cause the implant to break.
- Handle removed implants in such a way that their reuse is not possible.

MATERIAL

Components of the VENUS™ system are made out of the titanium alloy Ti6Al4V, according to ISO 5832-3. There are also rods included, that are made of cobaltchromium alloy CoCr28Mo6 according to ISO 5832-12. CoCr28Mo6 contains cobalt, which is classified as a CMR substance. A risk-benefit analysis must be performed prior to use.

PACKAGING, LABELING, TRANSPORT AND STORAGE

- The devices must be handled, transported and stored with care. Damages to the device packaging or damages to the device itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increased internal loads, which can result in a fracture of the device.
- The implants and corresponding instruments should

be stored at room temperature. Environmental factors such as salty air, humidity, direct sunlight, chemicals, etc., must not be allowed to act on the devices.

- Before the surgery, a careful inspection of the components of the VENUS™ system to be used must be carried out in order to avoid damages caused by storage, transport or previous procedures. The sterilization containers, trays and the associated covers must not be damaged.

The devices are delivered in both UNSTERILE and STERILE packaged condition:

Unsterile delivered devices:

UNSTERILE delivered devices are labelled as UNSTERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILISATION). The devices are supplied in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

Sterile delivered devices:

STERILE delivered devices are sterilized by a validated gamma sterilization procedure and labelled as STERILE. Cleaning, disinfection and sterilization prior to use need not be performed.

The devices are supplied packaged in a blister-blister combination and protected by an outer carton. The devices may only be used when the label on the outer packaging and also the inner packaging is intact. If the packaging is damaged or already open, the sterility of the device is not guaranteed, and the device must not be used.

The devices must not be used if the indicated expiry date has been exceeded. Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or if the sterile packaging is damaged is not intended.

HumanTech Spine GmbH assumes no responsibility for the use of re-sterilized devices regardless of the person who performed the resterilization or the method used.

The rules of asepsis must be observed when removing the devices from the sterile packaging. The sterile packaging may only be opened immediately prior to the insertion of the device. It is recommended that a replacement device is always kept available. The implant must be removed from the packaging using appropriate aseptic precautions.

CLEANING, DISINFECTION AND STERILISATION

UNSTERILE packaged delivered products of the VENUS™ system must be cleaned, disinfected and sterilized prior to use. All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the

instruction “Processing instructions for implants_HuT Spine”. The latest version can always be found on our homepage:

“<https://humantech-spine.de/378-en-IFUs.html>”, as well as on request from HumanTech Spine GmbH.

DISPOSAL

Any HumanTech implant that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol.

Devices must be disposed of as a healthcare medical device in accordance with hospital procedures.

PRODUCT COMPLAINTS

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a product of the VENUS™ system ever exhibits “malfunctioning” (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, must be reported immediately to the manufacturer and the competent authority of the respective country by phone, fax, or in writing.

Should you have a complaint, please provide the name, item number and batch number of the component as well as your name, address, and detailed description of the fault as possible in written form.

Information to be provided to the patient

- The patient should be informed about contraindications, warnings, precautions, side effects and complications with HumanTech devices.
- The patient must be informed about MRI compatibility regarding the HumanTech products used.
- The patient should be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area
- An implant card must be given to the patient
- The patient has to be informed that the implant card must be shown to the treating physician prior to every planned examination, diagnostic intervention and surgical or non-surgical treatment.

FURTHER INFORMATION

In the event of complaints, suggestions, or indications regarding the content of these instructions for use or the use of the product, please refer to address listed on the last page.







After the European Medical Devices Database is online, the summary of safety and clinical performance of VENUS™ implants will be available at the following address

„<https://ec.europa.eu/tools/eudamed>“. Until then, you can find them on our homepage www.humantech-spine.de, as well as on request from HumanTech Spine GmbH.

If required, the instructions for use can be ordered in paper form from the address given on the last page. It will be sent by post within seven days.

© HUMANTECH Spine GmbH. All rights reserved.

EXPLANATION OF SYMBOLS

	Manufacturer
	Date of manufacture
	Do not re-use
	Lot number
	Catalogue number
	Quantity
	Medical device
	Unique Device Identifier
	Consult instructions for use
	Keep dry
	Non-sterile
	Sterilized using irradiation
	Single sterile barrier system with protective packaging inside
	Do not resterilize
	Use-by date
	Do not use if package is damaged and consult instructions for use
	Caution
	MR conditional
	Contains hazardous substances

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