# VENUS<sup>®</sup>

## Bone Filler

Gebrauchsanweisung	DE	2 - 3
Instruction for use	EN	4 - 5
Manual de instrucciones	ES	6 - 7
Kullanım için talimatlar	TR	8 - 9
Инструкция за употреба	BG	10 - 11
Návod k použití	CZ	12 - 13
Οδηγίες χρήσης	EL	14 - 15
Istruzioni per l'uso	IT	16 - 17
Gebruiksaanwijzing	NL	18 - 19
Instrução de utilização	PT	20 - 21
Instrucțiuni de utilizare	RO	22 - 23
Инструкция по применению	RU	24 - 25
Inštrukcie na používanie	SK	26 - 27
Használati útmutató	HU	28 - 29
Lietošanas instrukcija	LV	30 - 31

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## IMPORTANT INFORMATION

This instruction for use must be read carefully and its content must be adhered to.

## DESCRIPTION

Bone Filler is a medical aid for applying bone cement to vertebral bodies using a cement adapter and the cannulated pedicle screws of the VENUS fixation system. The Bone Filler is screwed onto the cement adapter of the VENUS fixation system after being filled. The cement is applied by means of the push rod. The Bone Filler has a capacity of 1.5ml.

## INDICATIONS FOR USE

The Bone Filler is intended for applying bone cement to vertebral bodies using a cement adapter and the cannulated pedicle screws of the VENUS fixation system.

## AREA OF APPLICATION

The area of application is restricted to applying bone cement to vertebral bodies using a cement adapter and the cannulated pedicle screws of the VENUS fixation system.

- The Bone Filler may only be used by surgeons who have completed the necessary training in the area of spinal surgery. The use of the Bone Filler for applying cement must be decided according to the surgical and medical indications, the potential dangers, and the restrictions associated with this type of surgical intervention as well as in awareness of contraindications, side effects and defined precautionary measures.
- It is recommended that the Bone Filler is not used together with components of other spinal fusions systems from another supplier or another manufacturer. If this is the case, HumanTech Spine GmbH does not assume any liability.
- The Bone Filler is designed for one-time use on one patient. It may not be reused, resterilized nor reprocessed. Reusing, resterilizing or reprocessing may lead to the structural integrity of the product being compromised or to the product malfunctioning, which could lead to injury or illness of the patient. Furthermore, reusing, resterilizing or reprocessing may lead to an increased risk of contamination of the product and/or to the patient becoming infected, including but not limited to passing on infectious diseases from one patient to another.
- Observe general precautionary measures in order to avoid contact with pathogens in blood.
- Caution: The operating manual of the cement must also be observed in every situation.
- The amount of cement used must be chosen by the surgeon depending on the individual anatomy of the patient and the cement used.

Contraindications include, but are not restricted to the contraindications of the VENUS fixation system and the contraindications of the cement used. Side effects include, but are not restricted to the side effects of the VENUS fixation system and the side effects of the cement used:

- Leakage of cement into surrounding tissues and burns associated with this while the cement is hardening.
- Injuries resulting from the Bone Filler breaking.
- Bone Fillers are distributed STERILE. They may not be used if the inner packaging has been compromised. If the packaging is damaged or has already been opened, the sterility of the implant is not guaranteed and it must not be used.
- Bone Fillers must not be used if the expiry date given has been exceeded.
- Bone Fillers must not be resterilized.
- Bone Fillers must be handled and stored with care. Damage to the product can significantly reduce its stability. It can lead to cracks and/or increased internal tensions which can result in the rupture of the product.
- Bone Fillers should be stored between +5°C and +40°C.
- Environmental factors such as salt-laden air, moisture, chemicals etc. must not affect the products.
- An inspection and an assembly attempt are recommended before the operation in order to ensure that the products interact with each other correctly.

The instructions in the operation instructions (Surgical Technique of the used VENUS system) are to be observed. You can download these from the internet from [www.humantech-spine.de](http://www.humantech-spine.de) or directly from your HumanTech representative. It is advised to use the instruments described in the surgical technique, which are intended to be used with the system. No liability can be accepted for the use of external instruments.

## PRODUCT COMPLAINTS


Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who was dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/ or function should inform the relevant HumanTech representative. If a Bone Filler ever exhibits "malfunctioning" (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately. If a HumanTech product ever exhibits malfunctioning which has caused or contributed to the death or serious injury of the patient, the representative should be informed immediately by phone, fax or in writing. Should you have a complaint, we would ask you to send us the names, item numbers and the batch number of the components as well as your name and address together with a fault description with as much detail as possible in writing.

## FURTHER INFORMATION

In the event of complaints, suggestions or indications regarding the content of these instructions for use or the use of the product, please refer to the above-mentioned address.

© HumanTech Spine GmbH. All rights reserved.

## SYMBOLS

	Manufacturer
	Date of manufacture
	Expiry date
	Do not re-use
	Batch number
	Order number
	Read the documentation
	Sterilization with ethylene oxide
	Do not use if the packaging is damaged
	Phthalate-free
	Recommended storage temperature: +5°C to +40°C
	Pyrogen-free

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© HumanTech Spine GmbH. Visas tiesības aizsargātas.