VENUS®

POSTERIOR SPINAL IMPLANT SYSTEM

HumanTech Spine GmbH
Gewerbestraße 5
D-71144 Steinenbronn
Tel: +49 (0) 7157/5246-66
Fax: +49 (0) 7157/5246-71
e-mail: info@humantech-spine.de
www.humantech-spine.de

Stand: 25.10.2019

DE

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EN

IMPORTANT NOTE

This instruction for use must be read carefully and its content must be adhered.

SCOPE

The instruction for use is applicable for all sterile and unsterile delivered implants of the VENUS®-System.

BASIC STRUCTURE

The VENUS®-System, which is used as an implant system for spinal surgery and consists of rods, pedicle screws, transverse connectors, rod connectors and extension units as well as different types of hooks. Various forms and sizes of these implants are available so that adaptations can be utilized to take into account the unique pathology of individual patients.

Fenestrated Screws contain a series of fenestrations which allows polymethylmethacrylate (PMMA) bone cement to be injected into the treated site.

MATERIAL

Components of the VENUS-System are made out of the titaniumalloy Ti6Al4V, according to ISO 5832-3. There are also rods included, that are made of cobalt-chromium alloy CoCr28Mo6 according to ISO 5832-12. The Cementadapter CPS is made out of PEEK material.

INDICATIONS FOR USE

The VENUS®-System is intended for the surgical treatment of diseases and injuries of the spine such as instabilities, degenerative disc disease, degenerative spondylolisthesis, degenerative stenosis, scoliosis, kyphosis, spinal fractures, and spondylitis, tumor as well as all conditions that necessitate revision surgery. The VENUS® Mini Subsystem, which is part of the VENUS®-System is suitable for percutaneous surgeries.

LEVEL OF FIXATION

Levels of fixation are for the thoracic, lumbar, and sacral spine.

GENERAL CONDITIONS / INFORMATIONS FOR USE

- The implants must be implanted only by surgeons having undergone the necessary training in spinal surgery. The insertion of implants must be carried out in accordance with the surgical and medical indications, the potential risks, and limitations related to this type of surgery; the contraindications, side effects, and precautions defined, and in the knowledge of the nature and physical, chemical, metallic, metallurgical and biological characteristics of the implants.
- This instruction for use provides essential information, but is not sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a clearly defined post-operative treatment regimen, and to conduct scheduled follow-up examinations.
- The activities and kinesic behavior of the patient have a significant influence on the implant's useful life. Patients must be informed that every activity increases the risk of loss, loosening, dislocation, migration, bending or breakage of implant components. Informing patients about limitations to their activities in the postoperative phase and postoperatively monitoring of the patients are crucial factors in assessing the development of the surgical result, the fusion and the condition of the implant. The above-mentioned effects may occur even if a bony fusion occurred and the implant is well integrated and the activity restrictions are complied with. The patient must be informed of this. If such an

effect occurs, the surgeon must decide whether a revision of the implant or other measures should be taken, taking into account the well-being of the patient and the risks involved.

- It is recommended not to replace components of VENUS® system with components or products of other systems, or make a direct connection to components of other systems or to use components or products made out of material differing from Ti6Al4V or CoCr28Mo6 or to use with components of other spinal fusion systems from a different source, manufacturer, or material other than Ti6Al4V and CoCr28Mo6. HumanTech Spine GmbH shall accept no liability if this recommendation is not observed.

- Never reuse the implants. Even if the implant appears to be intact following revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break.
- Implants may also not be reused when they have been in contact with body fluids or tissues of third persons.
- Handle removed implants in such a way that their reuse is not possible.
- The volumes of cement used for fenestrated / cannulated screws should ultimately be determined by the surgeon based on the individual patient anatomy.
- Bending of the rods / revision connector affect the biomechanical properties of the implant. Bending in the area of the fixation of the rod in the Poly- or Monoaxial Screw can have negative influences to the fixation of the rod – Bending in this area has to be avoided.

- The rods of the VENUS® Mini-System do have marks. It is not allowed to use the area between the ends of the rod and the marks for the fixation of the pedicle screw. In this area the rods may also not be bended.
- The system stability can be increased with a ventral support. In the lumbar area, especially when the use of pedicle screws with a diameter of 5,5mm or smaller can't be avoided, a ventral support is strongly recommended.
- Fluted Rods shall not be used in combination with cannulated or fenestrated screws. With higher stress, metallosis can be caused due to the connection geometry.
- The instructions in the operation instructions (Surgical Technique) are to be observed. You can download these from the internet from www.humantech-spine.de or directly from your HumanTech representative. It is advised to use the instruments described in the surgical technique, which are intended to be used with the system. Complications may occur if the implant is inserted with or without the intended instruments. No liability can be accepted for the use of external instruments

- Proceed with extreme caution in the region of the spinal cord and the roots of the nerves, since damage to the nerves can lead to the impairment of neurological functions.
- Breakage, slippage or incorrect use of the instruments or implants can injure the patient or the operating staff or result in more time being required for surgery.
- Residues consisting of implant material and/ or not from implant material should be removed.
- Damaged implants must not be implanted.
- Safety and compatibility of the device in the setting of magnetic resonance (imaging) have not been evaluated. No thermal test or test of migration has been performed on the device in this setting.

CONTRAINDICATIONS

Contraindications may be either relative or absolute. The selection of a particular implant must be weighed carefully against the overall assessment of the patient. The following conditions can have an adverse impact on the chances of successful surgery:

- acute infectious process or significant risk of infection (immunocompromise)
- Signs of local inflammation
- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Mental illness
- Grossly distorted anatomy caused by congenital abnormalities.
- Extreme malalignments, impairing the stability of the instrumentation
- Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of congenital abnormalities, elevation of sedimentation rate unexplained by other diseases, elevation of white blood cell count (WBC) or a marked left shift in the WBC differential count.
- Joint disease, bone absorption, osteopenia, osteomalacia, and/or osteoporosis. Osteoporosis or osteopenia is a relative contraindication since this condition may limit the degree of obtainable correction, stabilization, and/or the amount of mechanical fixation.
- Any neuromuscular disease that would place excess strain on the implant during the healing period.
- Suspected or documented allergy or intolerance to the material used. Appropriate tests should be carried out.
- Any case not needing a bone graft and fusion.
- Any case where the implant components selected for use would be too large or too small to achieve a successful result.

- All cases which require the use of components other than the metals or alloys used in this system.
- Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality.
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance.
- Any patient unwilling or not able to follow postoperative instructions.
- Any case not described in the indications.

SIDE EFFECTS / COMPLICATIONS

The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

- Delayed bone growth or no visible fusion and pseudarthrosis
- Neurological complications, paralysis, soft tissue lesions, and/or migration of the implant
- Breakage, loosening or deformation of the implant, as well as abrasion
- Implant breakage during implantation or when implanted
- Superficial or deep-set infection and inflammatory phenomena
- Allergic reaction to the implant material
- Reduction in bone density
- Neurological and spinal dura mater lesions from surgical trauma
- Genitourinary disorders, gastrointestinal disorders, vascular disorders including thrombus, bursitis, secondary bleeding, myocardial infarction, or death
- Disorder of anatomical structures.
- Fracture of a vertebra, the pedicle, and/or the sacrum.
- Presence of micro-particles around the implants, Metallosis
- Growth of the fused vertebra is altered
- Partial loss of the degree of correction achieved during surgery
- Modification of spinal curvature and stiffness of the vertebral column
- Pain, discomfort or abnormal sensations due to the presence of the implant.
- Pressure on the skin caused by components located in positions with insufficient tissue coverage over the implant, with potential penetration of the skin.
- Fracture, micro fracture, resorption, damage to or penetration of a vertebral body above or below the treated segment/s
- Physiological limitations, such as joint degeneration
- Bleeding and/or haematomas
- Revision surgery
- Development of respiratory problems, including pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.

PACKAGING, LABELING AND STORAGE

- Handle and store the implant components carefully. Damage to the implant or cutting, bending or scratching the surface of the implant, can significantly reduce the strength and long-term stability of the implant system. Cracking and/or higher internal stresses can result, possibly resulting in the breakage of the implant.
- The implants and instruments should be stored at room temperature. Environmental factors, such as salt-laden air, humidity, chemicals, etc., must not be allowed to act on the implants.
- Thorough inspection is recommended before operating in order to ensure that the instruments or implants have not been damaged during storage or previous procedures.

The implants are delivered in both NON STERILE and STERILE packaged condition:

UNSTERILE delivered implants:

UNSTERILE delivered implants are labeled as UNSTERILE and must be cleaned and sterilized before use (see below REPROCESSING and STERILISATION).

The implants are supplied as a complete set or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilisation in the original pouch packaging is not allowed.

The holes of cannulated screws have to be rinsed and blown out after each cleaning process.

STERILE delivered implants:

STERILE delivered implants are sterilized by a validated gamma sterilization procedure and labeled as STERILE. Cleaning, processing and sterilization prior to use need not be performed.

The implants are supplied individually packaged in a blister-blister-combination or peel-pouches and protected by a covering box. The implants may be used only when the label on the outer packaging and also the inner packaging are intact. If the packaging is damaged or already open, the sterility of the implant is not guaranteed and the implant may not be used.

The implants may not be used when the shelf life indicated has been exceeded. Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or with damaged sterile packaging is not intended.

REPROCESSING

UNSTERILE delivered products of the VENUS®-System must be processed prior to sterilisation. The necessary steps are described in the reprocessing instructions for Ti implants.

STERILISATION

UNSTERILE delivered products of the VENUS®-System must be sterilized or resterilized prior to use by means of validated steam sterilisation in an autoclave. The minimum recommended sterilisation parameters are as follows:

Steam	Pre-vacuum	132-134°C	5 minutes
Steam	Pre-vacuum	121°C	30 minutes

PRODUCT COMPLAINTS

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a VENUS® implant ever exhibits "malfunctioning" (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, the representative should be informed immediately by phone, fax, or in writing.

Should you have a complaint, please provide the name, item number and batch number of the component along with your name, address and detailed description of the defect in writing.

FURTHER INFORMATION

In the event of complaints, suggestions or indications regarding the content of these instructions for use or the use of the product, please refer to the above-mentioned address.

© HumanTech Spine GmbH. All rights reserved.

	Manufacturer
	Do not re-use
	Batch number
	Reference number
	Read the documentation
	Keep dry
	Non-sterile
	Sterilized by irradiation
	Date of expiry
	Do not use when packing is damaged
	Attention

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