VENUSnano

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EN

IMPORTANT INFORMATION

This instruction for use must be read carefully and its content must be adhered to.

BASIC STRUCTURE

The VENUSnano Spinal System is intended for use as a posterior and anterior implant system for fixation in spinal surgery for children and adults of small stature. It consists of rods, pedicle screws, fasteners, flat washers and a variety of hooks.

The VENUSnano VDS instrument kit contains mono-axial screws, washers, rods and set screws for a ventral insertion using a double-rod treatment. A bicortical screw with a minimum diameter of 5 mm should be used for a ventral treatment.

These implants are available in different shapes and sizes, which means that adjustments can be made according to the individual pathology of each patient.

MATERIAL

The VENUSnano Spinal System is made from Ti6Al4V, a titanium alloy certified in accordance with ISO 5832-3. It also contains rods made from cobalt-chrome (CoCr) certified in accordance with ISO 5832-12.

INDICATIONS FOR USE

The VENUSnano Spinal System is designed for the surgical treatment of diseases and injuries of the thoracic, lumbar, and sacral spine, in particular for indications such as instability, degenerative disc disease, deformities such as scoliosis and kyphosis, spondylolisthesis, trauma, tumours, and revision surgery.

As a guideline, the implants are designed for patients aged from 5 to 11 years and for patients weighing less than 45 kg.

AREA OF APPLICATION

The area of application is the thoracic, lumbar and sacral spine.

GENERAL CONDITIONS FOR USE

- The implants must only be inserted by surgeons who have completed the requisite training in spinal surgery. The insertion of implants must be carried out in line with the surgical and medical indications, the potential dangers, and with the restrictions associated with this type of surgery. It should also be performed with an awareness of the contraindications, side effects, defined precautionary measures, and with regard for the structure and the metallic, metallurgical and biological characteristics of the implant.
- This leaflet provides essential information, but this is not sufficient to enable you to use the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a

clearly defined post-operative treatment regimen, and to conduct scheduled follow-up examinations.

- Use of the VENUSnano Spinal System in conjunction with components of other fusion systems from other suppliers or manufacturers, or which are made from a material other than Ti6Al4V and cobalt-chrome (CoCr), is not recommended. HumanTech Spine GmbH shall accept no liability in such cases.
- The implants must not be used more than once under any circumstances. Although the implant may appear intact after revision surgery, changes to the interior of the implant, or small defects caused by loads and stresses, can cause the implant to break.
- Implants may also not be reused when they have been in contact with body fluids or tissues of third persons.
- Use of the implant in the area of the thoracic spine in children can stop the growth of the chest wall and lungs, with associated consequences. As a result, multiple operations must be performed during the growth period in order to adapt the structure to the spine as it grows.
- Adapting the structure to a spine that is fully grown can lead to a drastic reduction in the expected service life of the implant. It is therefore advisable to replace the rod and set screw.
- Bending of the Rods affect the biomechanical properties of the implant. Bending in the area of the fixation of the rod in the Poly- or Monoaxial Screw can have negative influences to the Fixation of the Rod – Bending in this area has to be avoided.
- The instructions in the operation instructions (Surgical Technique) are to be observed. You can download these from the internet from www.humantech-solutions.de or directly from your HumanTech representative. It is advised to use the instruments described in the surgical technique, which are intended to be used with the system. No liability can be accepted for the use of external instruments

CONTRAINDICATIONS

Contraindications may be either relative or absolute. The selection of a particular implant must be weighed carefully against the overall assessment of the patient. The following conditions can have an adverse impact on the chances of successful surgery.

Contraindications include, but are not limited to, the following:

- Acute infections or significant risk of infection (weakened immune system).
- Signs of local inflammation.
- Open wounds.
- Fever or leukocytosis.
- Morbid obesity.
- Pregnancy.
- Psychological illness.
- If the anatomy is too deformed due to congenital abnormalities.
- Any other medical or surgical condition which prevents possible improvement through the use of implant systems, such as the presence of congenital abnormalities, increased white blood cell count (WBC) or a marked shift to the left in the WBC differential blood count.
- Joint diseases, bone shrinkage, osteopenia, osteomalacia and/or osteoporosis.
- Any neuromuscular disease that would place excess strain on the implant during the healing period.
- Suspicion of a metal allergy or intolerance, and documented metal allergy or intolerance. Appropriate tests should be carried out.
- All cases in which fusion is not required.
- All cases in which the selected implant components are too large or too small to achieve a satisfactory result.
- All cases which require the use of components other than the metals or alloys used in this system.
- Any patient with inadequate tissue structure on the operational side, or an inadequate osseous bed or bone quality.
- Any patient who is not willing to follow post-operative instructions.
- All cases which are not described in the indications.

COMPLICATIONS/ SIDE EFFECTS

Complications and side effects include, but are not limited to, the following:

- Delayed bone growth or no visible fusion and pseudarthrosis.
- Neurological complications, paralysis, soft-tissue lesions and/or migration of the implant.
- Breakage, loosening or deformation of the implant, as well as abrasion
- Implant breakage during implantation or when implanted
- Partial loss of the level of correction achieved during the operation.
- Superficial or deep infection and inflammation.
- Allergic reaction to the implant material.

- Reduction of bone density.
- Neurological or spinal lesion in the dura mater caused by surgical trauma.
- Genitourinary disorders, gastrointestinal disorders, vascular disorders including thrombus, respiratory tract disorders including embolism, bursitis, secondary bleeding, myocardial infarction, or death.
- Disorder of anatomical structures.
- Restricted physiological function.
- Fracture of a vertebra, the pedicle, and/or the sacrum.
- Occurrence of microparticles in the area of the implant (metallosis).
- Modified growth of the fused spine.
- Stopped growth of the chest wall and lungs, with associated consequences.
- Change in the curvature and stiffness of the spine.
- Pain, discomfort or abnormal sensations due to the presence of the implant.
- Pressure on the skin caused by components located in positions with insufficient tissue coverage over the implant, with potential penetration of the skin.
- Fracture, micro fracture, resorption, damage to or penetration of a vertebral body above or below the treated segment/s
- Physiological limitations, such as joint degeneration
- Bleeding and/or haematomas
- Revision surgery
- Development of respiratory problems, including pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.

PACKAGING, LABELLING AND STORAGE

- The implants are NON-STERILE. They must be cleaned and sterilised before use.
- The implants can be supplied as a complete set or individually. The packaging must be unopened at the time of delivery. The original bag packaging is not intended to be used for sterilization.
- These implants and instruments are kept in specially designed trays/containers which can be sterilised directly using approved sterilisation processes.
- The implant components must be handled and stored with care. Cutting, bending or scratching the surface of the implant can significantly reduce the strength and durability of the implant system. This may cause cracks and/or increased internal tension, which can cause the implant to break.
- Stored implants and instruments must be protected from corrosive environments such as salty air, moisture, etc. It is recommended that the equipment is inspected and a trial assembly carried out before surgery to ensure that the instruments or implants have not been damaged during storage or previous procedures.

STERILISATION

The components of the VENUSnano Spinal System must be sterilised or re-sterilised prior to use by means of validated steam sterilisation in an autoclave. Consult the autoclave manufacturer's instruction manual for details of the correct sterilisation temperature and cycle time. Steam sterilisation should be performed in accordance with EN ISO 17665-1. The minimum recommended sterilisation parameters are as follows:

Steam	Pre-vacuum	132-134°C	5 minutes
Steam	Pre-vacuum	121°C	30 minutes

IMPORTANT INFORMATION FOR THE SURGEON

- The service life of an implant is limited and is significantly influenced by the activities and movement patterns of the patient.
- The patient must be informed that any activity increases the risk of the loss, bending or breakage of implant components. It is imperative that the patient is informed about restrictions to his/her activities during the postoperative period, and that the patient is monitored after surgery in order to assess the development of the fusion and the condition of the implant components. Even if a solid bone fusion has occurred, implant components can still bend, break or loosen. The patient must therefore be informed that implant components may also bend, break or loosen even if he/she complies with the restriction of his/her activities. If an implant does break, the doctor must decide whether a revision of the implant should be performed, taking into account the wellbeing of the patient and the risks involved.
- The patient should be informed about the side effects.

PRODUCT RECLAMATIONS

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a component of the implanted VENUSnano Spinal System ever shows signs of a defect (i.e. does not meet the performance specifications or does not function as intended), or if such a defect is suspected, the HumanTech representative should be notified immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, the representative should be informed immediately by phone, fax, or in writing.

Should you have a complaint, please provide the name, item number and batch number of the component along with your name, address and detailed description of the defect in writing.

FURTHER INFORMATION

Please contact us at the above address if you have any complaints, suggestions or comments about the content of this manual or the use of the product.

© HUMANTECH Spine GmbH. All rights reserved.

	Manufacturer
	Only for one-time use
	Batch number
	Order number
	Read this documentation
	Keep dry
	Unsterile
	Do not use if the packaging is damaged
	Attention

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© HUMANTECH Spine GmbH. Reservados todos os direitos.