# HERO®

## CERVICAL PLATE SYSTEM

|                          |    |         |
|--------------------------|----|---------|
| Gebrauchsanweisung       | DE | 2 - 5   |
| Instruction for use      | EN | 6 - 9   |
| Mode d'emploi            | FR | 10 - 13 |
| Manual de instrucciones  | ES | 14 - 17 |
| Instruções de utilização | PT | 18 - 21 |
| Inštrukcie na používanie | SK | 22 - 25 |
| Οδηγίες χρήσης           | EL | 26 - 29 |
| Lietošanas instrukcija   | LV | 30 - 33 |
| Kullanım için talimatlar | TR | 34 - 37 |
| Használati utasítás      | HU | 38 - 41 |

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## IMPORTANT INFORMATION

This instruction for use must be read carefully and its content must be adhered to.

## SCOPE

The instruction for use is applicable for all sterile and unsterile delivered implants of the HERO® Cervical Plate System.

## BASIC DESIGN

The HERO® Cervical Plate System, which is used as an implant system for fixing in spinal surgery, consists of plates and screws of various sizes making adjustments possible, taking into account the unique anatomy of the individual patient.

## MATERIAL

All components are made of titanium alloy Ti6Al4V according to DIN EN ISO 5832-3.

## INDICATIONS FOR USE

The HERO® Cervical Plate System is designed for the surgical treatment of diseases and injuries of the cervical spine from C2 to C7 in patients in whom general skeletal growth has ceased, in particular for indications such as:

- degenerative diseases of the discs
- fractures
- tumors
- pseudarthrosis
- spinal canal stenosis
- cervical myelopathy
- deformities (i.e. cyphosis, lordosis and/or scoliosis)
- surgical revisions

## AREAS OF APPLICATION

The area of application is the cervical spine from C2 to C7.

## GENERAL CONDITIONS / INFORMATION FOR USE

- The implants may only be emplaced by surgeons who have successfully completed the necessary training in the field of spinal surgery. The decision to use implants must be made in accordance with the surgical and medical indications, taking into account the potential dangers and the limitations associated with this type of surgical intervention, as well as in knowledge of the contraindications, side effects and the defined precautionary measures and in awareness of the properties and the physical, chemical and biological characteristics of the implant.
- The activities and physical activity of the patient have a significant influence on the implant's useful life. Patients must be informed that every activity increases the

risk of loss, loosening, dislocation, migration, bending or breakage of implant components. Informing patients about limitations to their activities in the postoperative phase and postoperatively monitoring patients are crucial factors in assessing the development of the surgical result and the condition of the implant. The above-mentioned effects may occur even if the implant is well integrated and the activity restrictions are complied with. The patient must be informed of this.

- This instruction for use provides essential information, but is not sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a clearly defined post-operative treatment regimen, and to conduct scheduled follow-up examinations.
- For complications, the doctor must decide whether a revision of the implant or other measures should be taken, taking into account the well-being of the patient and the risks involved.
- Components of the HERO® system may not be replaced by components / products from other systems from another source or from a different manufacturer. Furthermore, no direct connection of components / products of the system to components of other systems may be established. If this is not complied with or if the products are otherwise used or used improperly, HumanTech Spine GmbH assumes no responsibility.
- The instructions in the operation instructions (Surgical Technique) are to be observed. You can download these from the internet from [www.humantech-spine.de](http://www.humantech-spine.de) or directly from your HumanTech representative. It is advised to use the instruments described in the surgical technique, which are intended to be used with the system. No liability can be accepted for the use of external instruments.
- Proceed with extreme care in the vicinity of the spinal cord and the nerve roots, as any damage to nerves can result in the loss of neurological functions.
- Fracture, slipping or incorrect use of instruments or implants can result in injury to the patient or surgical personnel.
- Bone cement may not be used as it makes the removal of components difficult or even impossible.
- Residues consisting of implant material and/ or not from implant material should be removed.
- Damaged implants must not be implanted.
- Under no circumstance may the implants be re-used. Even if the implant appears intact after the revision, there may be alterations within the implant or small defects, which are attributable to the action of loads and strain and which may result in fracture of the implant.

- Implants that have already come into contact with a patient's bodily fluids or tissues or have been soiled must not be reused.
- Once removed, implants are to be treated so as to render any re-use impossible.
- Safety and compatibility of the device in the setting of magnetic resonance (imaging) have not been evaluated. No thermal test or test of migration has been performed on the device in this setting.
- Bending of the plates can have an adverse effect on the biomechanical properties of the implant. Deformations in the vicinity of locating holes can adversely affect the anchoring of the system. Any bending in this area is to be avoided.
- all cases in which no fusion is required
- all cases in which the implant components selected for use are too big or too small to achieve a satisfactory result
- any patient with inadequate tissue structure at the surgical site or with inadequate bone stock or bone quality
- Any patient in whom the use of implant is hindered by anatomical structures or where these structures would restrict the physiological performance.
- any patient who is unwilling to follow postoperative instructions
- all cases not described under the indications

## CONTRAINDICATIONS

Contraindications may be either relative or absolute. The selection of a particular implant must be weighed carefully against the overall assessment of the patient. The following conditions can have an adverse impact on the chances of successful surgery:

- All destructive and inflammatory diseases of the cervical vertebrae and of the motor segments such as tumors, spondylitis and fractures with a high degree of instability, with pronounced static and structural alteration of the spine and the segment structure, which require a vertebral body replacement and other measures to reconstruct and stabilize the segment.
- acute infections or significant risks of infection (weakened immune system)
- signs of local inflammations
- fever or leukocytosis
- morbid obesity
- pregnancy
- physical diseases
- severe anatomical deformities caused by congenital abnormalities
- Any other medical or surgical condition, which prevents a possible improvement by the use of the implant, such as the presence of congenital abnormalities, a fracture close to the operation site, increase of the deposition rate which is not explained by other disease patterns, elevated white blood cell count (WBC) or a pronounced leftward shift in the differential WBC.
- Articular diseases, bone resorption, osteopenia, osteomalacia and/or osteoporosis are relative contraindications as these may limit the degree of achievable correction or stabilization.
- Any neuromuscular disease that would place excess strain on the implant during the healing period.
- known hereditary or acquired bone brittleness or calcification problems
- Suspicion of a metal allergy or intolerance, and documented metal allergy or intolerance. Appropriate tests should be carried out
- all cases in which the use of components made of various metals or alloys is required

## SIDE EFFECTS / COMPLICATIONS

The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

- loosening of the implant
- fracture of the implant
- foreign body reactions to the implants including possible formation of tumors, autoimmune disease or scar formation
- neurological complications
- paralysis
- lesions of soft tissue
- neurological or spinal lesion of the dura mater due to surgical trauma
- swallowing difficulties (dysphagia), hoarseness
- damage to the lymph vessels / lymph leakage
- changes to the curvature and rigidity of the spinal column
- failure to achieve the desired result of surgery, implantation failure
- partial loss of degree of correction achieved in operation
- superficial or internal infections and inflammations such as discitis, arachnoiditis
- delayed engraftment of bone or no visible fusion and pseudarthrosis
- fracture, microfracture, resorption, damage or penetration of a vertebra and/or bone transplant and/ or bone transplant extraction site above or below the treated segment
- bleeding and/or hematomas
- deep vein thrombosis, thrombophlebitis, lundenembolus
- complications at bone transplant donor site
- herniated disc, spinal disc destruction or degeneration, above or below the treated segment
- loss or increase of spinal mobility or functions
- damage to the reproductive system, sterility and sexual dysfunction
- development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia etc.
- changes in the mental condition
- changed growth of the fused spinal column
- allergic reaction to the Ti6Al4V alloy

- appearance of microparticles in the vicinity of the implant
- reduction of bone density
- bursitis
- death

## PACKAGING, LABELING, TRANSPORT AND STORAGE

- The implant components must be handled, transported and stored with due care. Damages to the implant packaging or damages to the implant itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increase internal loads, which can result in a fracture of the implant.
- The implants and instruments should be stored at room temperature. Environmental influences such as salty air, humidity, chemicals etc. may not be allowed to act on the implants.
- Before the surgery, a careful inspection of the components of the HERO® Cervical Plate system to be used must be carried out in order to avoid damages due to storage, transport or previous procedures. The sterilization containers, trays and the associated covers must not be damaged.

The implants are delivered in both NON STERILE and STERILE packaged condition:

### Non sterile delivered implants:

NON STERILE delivered implants are labeled as NON STERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILIZATION).

The implants are supplied as an implant system in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

### Sterile delivered implants:

STERILE delivered implants are sterilized by a validated gamma sterilization procedure and labeled as STERILE. Cleaning, disinfection and sterilization prior to use must not be performed.

The implants are supplied individually packaged in a blister and protected by a covering box. The implants may be used only when the label on the outer packaging and also the inner packaging are intact. If the packaging is damaged or already open, the sterility of the implant is not guaranteed and the implant may not be used.

The implants may not be used when the shelf life indicated has been exceeded.

Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or with damaged sterile packaging is not intended.

HumanTech Spine GmbH assumes no responsibility for the use of re-sterilized implants regardless of the person who performed the re-sterilization or the method used. The rules of asepsis must be observed when removing the implant from the sterile packaging. The sterile packaging may only be opened immediately before the implant is inserted. It is recommended to always have a replacement implant available. The implant must be removed from the packaging using appropriate aseptic provisos.

## CLEANING, DISINFECTION AND STERILIZATION

NON STERILE packaged delivered products of the HERO® Cervical Plate System must be cleaned, disinfected and sterilized prior to use. All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the instruction "Processing unsterile delivered Spine". The latest version can always be found on our homepage <https://humantech-spine.de/378-en-IFUs.html>, as well as on request from HumanTech Spine GmbH.

## DISPOSAL

The products must be disposed of in accordance with the applicable local regulations, whereby the respective degree of contamination must be taken into account.

## PRODUCT COMPLAINTS

Any person involved in healthcare (e.g. customer or user of this product system) who has complaints of any kind or is dissatisfied using the product with respect to quality, identity, strength, durability, safety, effectiveness or function should notify the appropriate HumanTech representative. If a HERO® implant should ever malfunction (i.e. fails to fulfill the performance specifications or does not function in the foreseen manner), or if the above is suspected, the HumanTech representative should be informed immediately.

If a HumanTech product should ever malfunction in a way that results in the death or serious injury of the patient or if it has contributed to the above, the representative is to be informed immediately by telephone, fax or in written form. If you have any complaints then we ask you to state the name, article number and batch number of the components as well as your name and address along with a detailed description of the fault and forward this information to us in writing.



## FURTHER INFORMATION

If you have any complaints, suggestions or other points about the contents of this instruction for use or the use of the product, please contact the above address.

© HUMANTECH Spine GmbH. All rights reserved.

## SYMBOLS

|   |                                      |
|---|--------------------------------------|
|    | Manufacturer                         |
|    | Do not re-use                        |
|    | Batch number                         |
|    | Reference number                     |
|    | Read the documentation               |
|    | Keep dry                             |
|    | Unsterile                            |
|    | sterilized by irradiation            |
|    | Do not re-sterilize                  |
|   | Date of expiry                       |
|  | Do not use when packaging is damaged |
|  | Attention                            |

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© HUMANTECH Spine GmbH. Minden jog fenntartva