# Distraction Pin<sup>®</sup>

## FOR CERVICAL APPLICATIONS

|                      |    |       |
|----------------------|----|-------|
| Gebrauchsanweisung   | DE | 2 - 3 |
| Instruction for use  | EN | 4 - 5 |
| Instrucciones de uso | ES | 6 - 7 |
| Használati utasítás  | HU | 8 - 9 |

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## IMPORTANT NOTE

This instruction for use must be read carefully and its content must be adhered to.

## SCOPE OF APPLICATION

The instruction for use is applicable for all sterile and non sterile delivered Distraction Pins.

## BASIC STRUCTURE

The Distraction Pin is used in the distraction of cervical vertebral bodies. It is screwed into the vertebral bodies, which are distracted using another instrument (e.g. Retrieval Body Retractor).

## MATERIAL

The Distraction Pin consists of stainless steel 1.4441 that corresponds to the requirements on implanted steels (ISO 5832-1).

## INDICATIONS FOR USE

The Distraction Pin is intended for anterior surgeries in the cervical region of the spine, for which distraction is required in order to insert a disc replacement implant.

## AREA OF APPLICATION

The area of application is the cervical spine.

## GENERAL CONDITIONS OF USE

- Distraction Pins must only be used by surgeons who have completed the requisite training in spinal surgery.
- It is recommended that the Distraction Pins are not used together with instruments from another supplier or another manufacturer. HumanTech Spine GmbH shall accept no liability if this recommendation is not observed.
- Following the notes in the operating instructions (surgical technique) is essential. You will find these in the Internet at [www.humantech-spine.de](http://www.humantech-spine.de) or you can obtain them directly from your HumanTech representative.
- Proceed with extreme caution in the region of the spinal cord and the roots of the nerves, since damage to the nerves can lead to the impairment of neurological functions.
- Breakage, slippage or incorrect use of the instruments or implants can injure the patient or the operating staff.
- Do not use bone cement, as this material makes the removal of the components difficult or impossible.
- Residues consisting of instrument material should be

removed.

- Damaged products must not be used.
- The Distraction Pins may not be used multiple times. Even if the instrument appears intact after revision, alterations within the instrument or small defects due to acting stresses and strains during use can cause the instrument to break when it is reused
- Handle removed instruments in such a way that their reuse is not possible.

## CONTRAINDICATIONS

All cases that are not described in the indications and all cases in which the bone quality is not suitable for their application.

## SIDE EFFECTS / COMPLICATIONS

Complications can arise if the Distraction Pins are placed with or without the intended instruments.

Apart from the side effects that arise due to the disc replacement implant and the operation required for this, the following side effects can also occur but are not limited to:

- Breaking of the pins and thus remainder of the tip left in the vertebral body
- Foreign body reaction to the Distraction Pin and its material
- neurological complications, paralysis, soft-tissue lesions
- neurological or spinal lesion in the dura mater caused by surgical trauma
- superficial or deep infection and inflammation
- fracture, micro fracture, resorption, damage to or penetration of vertebral bone
- bleeding and/or haematomas
- reduction of bone density

## PACKAGING, LABELING AND STORAGE

- The product must be handled and stored with care. Damage can significantly reduce rigidity and fatigue strength, which can result in breakage.
- The product must be stored at room temperature. Environmental factors such as salt-laden air, moisture, chemicals etc. must not affect the products.
- It is recommended that a careful inspection is carried out before surgery to ensure that the instruments have not been damaged during storage or previous procedures.

Distraction Pins are delivered in both NON STERILE and STERILE packaged condition:

### Non sterile delivered products:

NON STERILE delivered products are labeled as NON STERILE and must be cleaned and sterilized before use (see REPROCESSING and STERILIZATION).

The products are supplied combined with the associated

instruments or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original pouch packaging is not allowed.

**Sterile delivered products:**

STERILE delivered products are sterilized by a validated gamma sterilization procedure and labeled as STERILE. Cleaning, processing and sterilization prior to use need not be performed.

The products are supplied individually packaged in a blister and protected by a covering box. The products may be used only when the label on the outer packaging and also the inner packaging are intact. If the packaging is damaged or already open, the sterility of the product is not guaranteed and the product may not be used.

The product may not be used when the shelf life indicated has been exceeded.

Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or with damaged sterile packaging is not intended.

**REPROCESSING**

Non sterile delivered Distraction Pins must be processed prior to sterilization. The necessary steps are described in the reprocessing instructions for surgical instruments.

**STERILIZATION**

Non sterile delivered Distraction Pins must be sterilized or re-sterilized prior to use by means of validated steam sterilization in an autoclave. The minimum recommended sterilization parameters are as follows:

|       |            |       |             |
|-------|------------|-------|-------------|
| Steam | Pre-vacuum | 132°C | ≥ 5 minutes |
| Steam | Pre-vacuum | 121°C | 30 minutes  |

**PRODUCT COMPLAINTS**

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a Distraction Pin ever exhibits “malfunctioning” (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, the representative should be informed immediately by phone, fax, or in writing.


Should you have a complaint, please provide the name, item number and batch number of the component along with your name, address and detailed description of the defect in writing.

**FURTHER INFORMATION**

In the event of complaints, suggestions or indications regarding the content of these instructions for use or the use of the product, please refer to the above-mentioned address.

© HumanTech Spine GmbH. All rights reserved.

**SYMBOLS**

|   |                                      |
|---|--------------------------------------|
|    | Manufacturer                         |
|    | Do not re-use                        |
|    | Batch number                         |
|    | Reference number                     |
|    | Read the documentation               |
|    | Keep dry                             |
|   | Do not re-sterilize                  |
|  | sterilized by irradiation            |
|  | Date of expiry                       |
|  | Do not use when packaging is damaged |
|  | Attention                            |
|  | Non sterile                          |

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© HumanTech Spine GmbH. Minden jog fenntartva.