Distraction Pins

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1501010022 Distraction Pin14mm

1501010023 Distraction Pin16mm

1501010024 Distraction Pin18mm

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EN

IMPORTANT INFORMATION

This instruction for use must be read carefully and its content must be adhered to.

BASIC STRUCTURE

Distraction Pins are used in the distraction of cervical vertebral bodies. They are screwed into the vertebral bodies, which are distracted using another instrument (e.g. Retrival Body Retractor)

MATERIAL

The Distraction Pins consist of stainless steel that corresponds to the requirements on implanted steels (ISO 5832-1).

INDICATIONS FOR USE

The Distraction Pins are suitable for anterior operations in the cervical region, for which distraction is required in order to insert a replacement intervertebral body.

AREA OF APPLICATION

The area of application is the cervical spine.

GENERAL CONDITIONS FOR USE

- The Distraction Pins must only be used by surgeons who have completed the requisite training in spinal surgery.
- It is recommended that the Distraction Pins are not used together with instruments from another supplier or another manufacturer. HumanTech Spine GmbH shall accept no liability if this recommendation is not observed.
- The Distraction Pins may not be used multiple times. Even if the instrument appears intact, re-use can cause the Distraction Pin to break.

CONTRAINDICATIONS

- All cases that are not described in the indications and all cases in which the bone quality is not suitable for their application.

SIDE EFFECTS

Complications can arise if the Distraction Pins are placed with or without the intended instruments.

Apart from the side effects that arise due to the intervertebral body replacement and the operation required for this, the following side effects can also occur but are not limited to:

- Breaking of the pins and thus remainder of the tip left in the vertebral body
- Foreign body reaction to the Distraction Pins and their material

- Neurological complications, paralysis, soft-tissue lesions
- Neurological or spinal lesion in the dura mater caused by surgical trauma
- Superficial or deep infection and inflammation
- Fracture, micro fracture, resorption, damage to or penetration of vertebral bone
- Bleeding and/or haematomas
- Reduction of bone density

PACKAGING, LABELLING AND STORAGE

- The Distraction Pins are UNSTERILE. They must be cleaned and sterilised before use (see below)
- The Distraction Pins are supplied packaged. The packaging must be unopened at the time of delivery.
- The Distraction Pins are kept with the associated instruments in specially designed trays/containers which can be sterilised directly using approved sterilisation processes. The original bag packaging is not intended to be used for sterilization.
- Product must be handled and stored with care. Damage can significantly reduce rigidity and fatigue strength, which can result in breakage.
- Product should be stored at room temperature. Environmental factors such as salt-laden air, moisture, chemicals etc. should not effect the products.
- It is recommended that a careful inspection is carried out before surgery to ensure that the instruments have not been damaged during storage or previous procedures.

STERILISATION

The Distraction Pins must be sterilised or re-sterilised prior to use by means of validated steam sterilisation in an autoclave. The correct settings with respect to sterilization temperature and cycle time shall be adopted from the instructions of the autoclave manufacturer. Sterilization may be accomplished by steam according to EN ISO 17665-1. The minimum recommended sterilization parameters are as follows:

Steam	Pre-vacuum	132°C	≥ 5 minutes
Steam	Pre-vacuum	121°C	30 minutes

IMPORTANT INFORMATION FOR THE SURGEON

- Proceed with extreme caution in the region of the spinal cord and the roots of the nerves, since damage to the nerves can lead to the impairment of neurological functions.
- Breakage, slippage or incorrect use of the instruments or implants can injure the patient or the operating staff.
- Handle removed instruments in such a way that their reuse is not possible.

PRODUCT COMPLAINTS

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a Distraction Pin ever exhibits "malfunctioning" (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, the representative should be informed immediately by phone, fax, or in writing.

Should you have a complaint, please provide the name, item number and batch number of the component along with your name, address and detailed description of the defect in writing.

FURTHER INFORMATION

In the event of complaints, suggestions or indications regarding the content of these instructions for use or the use of the product, please refer to the above-mentioned address.

© HumanTech Spine GmbH. All rights reserved.

	Manufacturer
	Only for one-time use
	Batch number
	Order number
	Read this documentation
	Keep dry
	Unsterile
	Do not use if the packaging is damaged.

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