TRISTAN®

CERVICAL INTERBODY FUSION

Gebrauchsanweisung	DE	2 - 6
Instruction for use	EN	7 -10
Kullanım Talimatı	TR	11 - 14
Manual de instrucciones	ES	15 - 19
Istruzioni per l'uso	IT	20 - 24
Instructions d'utilisation	FR	25 - 29
Gebruiksaanwijzing	NL	30 - 34
Instrukcja użytkowania	PL	35 - 39
Инструкция по применению	RU	40 - 44
Instruções de utilização	PT	45 - 49
Οδηγίες χρήσης	EL	50 - 54
Inštrukcie na používanie	SK	55 - 58
Használati utasítás	HU	59 - 62
Lietošanas instrukcija	LV	63 - 66
Návod k použití	CS	67 - 70

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IMPORTANT NOTICE

This instruction for use must be read carefully and its content must be adhered to.

SCOPE

The instructions for use is valid for sterile delivered implants:

- TRISTAN® PEEK...
- TRISTAN® PEEK...-S
- TRISTAN® R-PEEK-Ti...
- TRISTANflex...
- TRISTAN® Ti... sterile,

and unsterile delivered implants:
TRISTAN® Ti...

BASIC DESIGN

The TRISTAN® Cervical Interbody Fusion, used as an implant system for the fusion of vertebral bodies in spinal surgery, is comprised of intervertebral implants in different sizes, angles and heights, so that adaptations taking into account the unique anatomy of individual patients are possible.

MATERIAL

There are different implant product lines.

Basic material PEEK Optima®:

TRISTAN® PEEK...

The base body of products in the TRISTAN® PEEK range is exclusively made from PEEK Optima®. To achieve a better radiographic contrast, radiographic markers made out of titanium alloy Ti6Al4V respectively tantalum have been incorporated into the implant.

TRISTAN® R-PEEK-Ti...

The base body of products in the TRISTAN® R-PEEK- Ti range also consists of PEEK Optima®. In addition, the surfaces of the implant have a titanium coating. To achieve a better radiographic contrast, radiographic markers made out of titanium alloy Ti6Al4V respectively tantalum have been incorporated into the implant.

TRISTANflex and TRISTAN® PEEK...-S

The base body of products in the TRISTANflex and TRISTAN® PEEK...-S ranges also consists of PEEK Optima®. To achieve a better radiographic contrast, radiographic markers made out of titanium alloy Ti6Al4V respectively tantalum have been incorporated into the implant. In addition, there are other components made out of titanium alloy Ti6Al4V included (e.g. spikes).

Basic material Titanium alloy:

TRISTAN® Ti and TRISTAN® Ti ... sterile

Products of the TRISTAN® Ti and TRISTAN® Ti...sterile range consist of the titanium alloy Ti6Al4V.

INDICATIONS FOR USE

TRISTAN® is intended for the operative treatment of diseases and injuries of the cervical spine in patients whose overall skeletal growth is fully developed, in particular for such indications as:

- disk herniation
- hard disk herniation
- mechanical instabilities
- calcification of the posterior longitudinal ligament
- osteochondrosis (for TRISTANflex reduced to Modic I and II)
- spinal canal stenosis

AREA OF APPLICATION

The area of application is the cervical spine from C3 to C7.

GENERAL CONDITIONS FOR USE

- Only surgeons with the required professional training in the field of spinal surgery may insert the implants. The use of the implants must be decided in keeping with the surgical and medical indications, the potential risks and the limitations relating to this type of surgical intervention, as well as with the knowledge of contraindications, side effects and known precautionary measures and the knowledge of the condition and the physical, chemical and biological properties of the implant.
- The activities and physical activity of the patient have a significant influence on the implant's useful life. Patients must be informed, that every activity increases the risk of loss, loosening, dislocation, migration, bending or breakage of implant or implant components. Informing patients about limitations to their activities in the postoperative phase and postoperatively monitoring patients are crucial factors in assessing the development of the surgical result and the condition of the implant. The above- mentioned effects may occur even if the implant is well integrated and the activity restrictions are complied with. The patient must be informed of this. If such an effect occurs, the surgeon must decide whether a revision of the implant or other measures should be taken, taking into account the well-being of the patient and the risks involved.
- This instruction for use provides essential information, but is not solely sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a

- clearly defined post-operative treatment regimen, and to conduct scheduled follow-up examinations.
- Components of the TRISTAN® system may not be replaced by components/ products from other systems, from another source or from a different manufacturer. Furthermore, no direct connection of components/ products of the systems to components of other systems may be established. If this recommendation is not complied with or if the products are otherwise used or used improperly, HumanTech Spine GmbH assumes no liability responsibility.
- Following the notes in the operating instructions (surgical technique) is essential. You will find these in the Internet at www.humantech-spine.de or you can obtain them directly from your HumanTech representative. It is advisable to use the instruments described there and provided by the manufacturer. There may be complications if the implant is used with or without the dedicated instrument. For the use of foreign instruments, no liability can be accepted.
- Proceed with extreme caution in the region of the spinal cord and the roots of the nerves, since damage to the nerves can lead to the impairment of neurological functions.
- Breakage, slippage or incorrect use of the instruments or implants can injure the patient or the operating staff.
- Do not use bone cement, as this material makes the removal of the components difficult or impossible. The heat produced by the hardening process can damage or deform the PEEK implants.
- Residues consisting of implant material and/ or not from implant material should be removed.
- Damaged implants must not be implanted.
- Never reuse the implants. Even if the implant appears to be intact after revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break.
- Implants that have already come into contact with a patient's bodily fluids or tissues or have been soiled must not be reused.
- Handle removed implants in such a way that their reuse is not possible.
- Safety and compatibility of the device in the setting of magnetic resonance (imaging) have not been evaluated. No thermal test or test of migration has been performed on the device in this setting.
- To achieve a good primary fixation, the TRISTANflex and TRISTAN® PEEK...-S Cages have spiky and sharp spikes, therefore they must be carefully handled and a sufficient distraction is necessary.
- It is recommended to use an additional anterior plating if the bone quality or the implant site can't be evaluated pre- or intraoperatively in question of primary fixation or if the implant site or the bone quality is expected not to ensure a proper primary fixation.
- When using PEEK implants particular care in preparing the intervertebral disc space is needed to keep

the placement forces low

- For a better ongrowing, the R-PEEK-Ti implants are coated with a titanium surface. Take care during insertion of these cages as the surface is structured and rough.

CONTRAINDICATIONS

Contraindications may be either relative or absolute. The selection of a particular implant must be weighed carefully against the overall assessment of the patient. The following conditions can have an adverse impact on the chances of successful surgery:

- all destructive and inflammatory diseases of the cervical vertebral body and vertebral motor segments such as tumors, spondylitis and fractures with a high degree of instability, with pronounced static and structural alterations of the vertebral body and segment structure, requiring the replacement of the vertebral body and other stabilizing measures for the reconstruction and stabilization of the segment
- acute infections or significant risks of infections (weakened immune system)
- indications of local inflammation
- fever or leucocytosis
- pathological obesity
- pregnancy
- psychological disorders
- pronounced anatomical deformities resulting from congenital abnormalities
- any other medical or surgical condition hindering possible improvement through the use of the implant, such as the existence of congenital abnormalities, fracture in the vicinity of the operating point, elevated deposition rate not indicated by other clinical pictures, elevated leucocyte count (white blood cells), or distinct leftward shift in the white blood cell differential blood count.
- Diseases of the joints, bone resorption, osteopenia, osteomalacia and/or osteoporosis are relative contraindications, since these can limit the degree of correction and stabilization that can be achieved.
- Any neuromuscular disease that would place excess strain on the implant during the healing period.
- known hereditary or acquired brittleness of the bones or calcification problems
- for TRISTANflex: Osteochondrosis Modic III
- spondylolisthesis that cannot be reduced to the first degree
- Suspicion of a metal allergy or intolerance, and documented metal allergy or intolerance. Appropriate tests should be carried out.
- all cases in which the use of components of different metals or alloys is necessary
- already existing fusion in the segments treated
- all cases in which fusion is not necessary
- all cases in which the implant component selected for use is too large or too small to achieve a satisfactory

- result
- all patients with inadequate tissue structure on the operative side or inadequate bone or bone quality
- all patients in which the use of the implant would interfere with anatomical structures or restrict physiological performance
- all patients not willing to follow postoperative instructions
- all cases not described in the indications

SIDE EFFECTS / COMPLICATIONS

The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

- loosening of the implant
- shifting of the implant
- breakage of the implant, respectively detached parts of the implant and wear
- foreign body reactions to the implant, including possible tumor formation, autoimmune disease and/ or scar formation
- scar formation with possible neurological impairments, nerve compression or pain
- neurological complications
- paralysis
- soft tissue lesions
- neurological or spinal lesions of the dura mater as a result of surgical trauma
- The migration of the implant in the posterior direction can cause a loss of neurological function, the occurrence of radiculopathy, dural tears and/or pain.
- Neurovascular impairment, including paralysis
- temporary or permanent retrograde ejaculation in men, or other severe injuries
- Discharge of cerebrospinal fluid
- inability to perform daily tasks
- urinary retention or loss of control over the bladder, or other types of urinary system impairment
- difficulty swallowing (dysphagia)
- alteration of the curvature and stiffness of the spine
- Failure to achieve the desired result of surgery, malpositioning of the implant
- partial loss of the degree of correction achieved during the operation
- superficial or deep-seated infections and inflammations, such as discitis or arachnoiditis
- bone fracture above or below the segments treated
- late osseointegration or no visible fusion and pseudarthrosis
- fracture, micro fracture, resorption, damage to or penetration of vertebral bone and/or the bone transplant and/or the bone transplant donor site above and below the segment treated
- adjacent segment disease
- bleeding and/or hematomas
- deep-seated venous thrombosis, thrombophlebitis, pulmonary embolism

- complications at the bone transplant donor site
- herniated disk, disk destruction or disk degeneration at, above or below the segment treated
- loss or enhancement of spinal mobility or functions
- damage to the reproductive system, sterility and sexual functional disturbances
- development of respiratory problems, for example pulmonary embolism, atelectasis, bronchitis, or pneumonia
- alteration of the mental state
- alteration in the growth of the fused spine
- allergic reaction to the materials used
- decrease in bone density
- bursitis
- death

PACKAGING, LABELING, TRANSPORT AND STORAGE

- The implant components must be handled, transported and stored with due care. Damages to the implant packaging or damages to the implant itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increase internal loads, which can result in a fracture of the implant.
- The implants and instruments should be stored at room temperature. Environmental influences such as salty air, humidity, chemicals etc. are not allowed to act on the implants.

The implants are delivered in both NON STERILE and STERILE packaged condition:

Non sterile delivered implants:

NON STERILE delivered implants are labeled as NON STERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILIZATION).

The implants are supplied as an implant system in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

Sterile delivered implants:

STERILE delivered implants are sterilized by a validated gamma sterilization or ethylene oxide procedure and labeled as STERILE. Cleaning, disinfection and sterilization prior to use must not be performed.

The implants are supplied individually packaged in a blister and protected by a covering box. The implants may be used only when the label on the outer packaging and also the inner packaging are intact. If the packaging is damaged or already open, the sterility of the implant is not guaranteed and the implant may not be used.

The implants may not be used when the shelf life indicated

has been exceeded.

Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or with damaged sterile packaging is not intended.

HumanTech Spine GmbH assumes no responsibility for the use of re-sterilized implants regardless of the person who performed the re-sterilization or the method used. The rules of asepsis must be observed when removing the implant from the sterile packaging. The sterile packaging may only be opened immediately before the implant is inserted. It is recommended to always have a replacement implant available. The implant must be removed from the packaging using appropriate aseptic provisos.

CLEANING, DISINFECTION AND STERILIZATION

NON STERILE packaged delivered products of the TRISTAN® Cervical Interbody Fusion System must be cleaned, disinfected and sterilized prior to use. All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the instruction "Processing unsterile delivered Spine". The latest version can always be found on our homepage <https://humantech-spine.de/378-en-IFUs.html>, as well as on request from HumanTech Spine GmbH.

DISPOSAL

The products must be disposed of in accordance with the applicable local regulations, whereby the respective degree of contamination must be taken into account.

PRODUCT COMPLAINTS

Any person involved in healthcare (e.g. customer or user of this product system) who has complaints of any kind or is dissatisfied using the product with respect to quality, identity, strength, durability, safety, effectiveness or function should notify the appropriate HumanTech representative. If a TRISTAN® implant should ever malfunction (i.e. fails to fulfill the performance specifications or does not function in the foreseen manner), or if the above is suspected, the HumanTech representative should be informed immediately.

If a HumanTech product should ever malfunction in a way that results in the death or serious injury of the patient or if it has contributed to the above, the representative is to be informed immediately by telephone, fax or in written form. If you have any complaints then we ask you to state the name, article number and batch number of the components as well as your name and address along with a detailed description of the fault and forward this information to us in writing.

FURTHER INFORMATION

If you have any complaints, suggestions or other points about the contents of this user manual or the use of the product, please contact the above address.

© HUMANTECH Spine GmbH. All rights reserved

SYMBOLS

	Manufacturer
	Do not re-use
	Batch number
	Reference number
	Read the documentation
	Keep dry
	Non sterile
	sterilized by irradiation
	Sterilization with ethylene oxide
	Do not resterilize
	Date of expiry
	Do not use when packaging is damaged
	Attention

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© HUMANTECH Spine GmbH. Všechna práva vyhrazena.