ADONIS® PEEK, PEEK-Ti, ADONIS®Ti sterile LUMBAR INTERBODY FUSION

HumanTech Germany GmbH
Gewerbestraße 5
D-71144 Steinheimbronn
Tel: +49 (0) 7157/5246-71
Fax: +49 (0) 7157/5246-33
e-mail:info@humantech-solutions.de
www.humantech-solutions.de

Stand: 12.04.2019

DE

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EN

IMPORTANT NOTICE

This instruction for use must be read carefully and its content must be adhered.

SCOPE

The instructions for use is valid for sterile delivered implants:

- ADONIS® PEEK ...
- ADONIS® R-PEEK-TI...
- ADONIS® Ti... sterile

The sterile delivered implants are packed in a blister and protected by a covering box. The label indicates the sterilization by the symbol "sterilized with ethylene oxide".

BASIC STRUCTURE

The ADONIS® Lumbar Interbody Fusion System is used as an implant system for the fusion of vertebrae bodies in spine column surgery (cages). The cages are available in different dimensions, so that adaptations taking into account the unique pathology of individual patients are possible.

MATERIAL

There are different implant product lines.

Basic material PEEK Optima®:

ADONIS® PEEK

For products of the ADONIS® PEEK line, the base body is entirely made of PEEK Optima®. To achieve a better radiographic contrast, radiographic markers have been incorporated into the implant.

ADONIS® R-PEEK-TI

For products of the ADONIS® R-PEEK-Ti line, the base body is also made of PEEK Optima®. In addition, the surfaces of the implant have a titanium coating. To achieve a better radiographic contrast, radiographic markers have been incorporated into the implant.

Basic material Titanium alloy:

ADONIS® Ti... sterile

Products of the ADONIS® Ti line consist of the titanium alloy Ti6Al4V. The titanium alloy products are available both sterile and non-sterile. This instruction manual is only valid for the sterile version.

INDICATIONS

ADONIS® is intended for the surgical treatment of diseases and injuries of the lumbar and lumbosacral spine in patients whose general skeletal growth is completed, especially for indications such as:

- Prolapsed intervertebral disc
- Hard prolapsed intervertebral disc
- Mechanical instabilities
- Calcification of the posterior longitudinal ligament
- Osteochondrosis
- Lumbar spinal stenosis

APPLICATION

The area of application is the lumbar and lumbosacral spine.

GENERAL TERMS AND CONDITIONS FOR USE

- The implants may only be implanted by surgeons who have completed the necessary training in spinal column surgery. The decision to use the implants for implantation must be made in accordance with the surgical and medical indications, the potential dangers, the constraints related to

this type of surgery, as well as in full knowledge of contraindications, side effects and defined precautions, and in full knowledge of both the type of implant and the physical, chemical and biological properties of the implant.

- It is not recommended to use ADONIS® in combination with implants of a different source or from another manufacturer. If this recommendation is adhered to, HumanTech Germany GmbH accepts no responsibility.
- Under no circumstances may the implants be used multiple times. Although the implant may appear intact after revision, changes in the interior of the implant or small defects, resulting from acting loads and stresses, can result in the rupture of the implant.
- Safety and compatibility of the device in the setting of magnetic resonance (imaging) have not been evaluated. No thermal test or test of migration has been performed on the device in this setting.

CONTRAINDICATIONS

Contraindications include but are not limited to:

- All destructive and inflammatory diseases of the vertebral body or motion segment such as tumours, spondylitis and fractures with a high degree of instability, with a strong static and structural change of the vertebral body and segment structure, which require vertebral body replacement and other stabilising measures for reconstruction and stabilisation of the segment.
- Acute infections or significant risks of infections (weakened immune system)
- Signs of local inflammations
- Fever or leukocytosis
- Adiposis
- Pregnancy
- Mental illness
- Any other medical or surgical condition which prevents the possible improvement through the application of the implant, such as the presence of congenital abnormalities, a fracture close to the site of operation, increase in the deposition rate which is not described by other disease patterns, increase in white blood cell counts (WBC) or a marked left shift in the WBC differential blood count.
- Diseases of the joints, bone resorption, osteopenia, osteomalacia and/or osteoporosis are relative contraindications, since these can limit the degree of correction and stabilization that can be achieved
- Known hereditary or acquired bone brittleness or calcification problems
- Spondylolisthesis, which cannot be reduced to the first degree.
- Suspicion of allergy or intolerance as well as documented allergy to or intolerance of the material used.
- All cases where the use of components of different metals or alloys is necessary
- Previous fusions in the segments to be treated
- All cases in which no fusion is required
- All cases in which the implant component selected for use is too large or too small to achieve a satisfactory result.
- Any patient with inadequate tissue structure on the operative side, or inadequate bone bed or bone quality
- Any patient in whom the application of the implant would interfere with anatomical structures or restrict physiological performance.
- Any patient who is not willing to adhere to postoperative instructions.
- All cases which are not described in the indications.

SIDE EFFECTS

Complications can occur if the implant is inserted with or without the armamentarium designated for this purpose.

Side effects include, but are not limited to:

- Implant loosening
- Implant migration
- Implant rupture
- Foreign body reaction to the implants including possible tumour formation, autoimmune disease and / or scarring
- Scarring with possible neurological impairment or nerve compression or pain
- Neurological complications, paralysis, soft-tissue lesions
- Neurological or spinal lesion of the dura mater through surgical trauma
- Pressure on surrounding tissue, nerves or organs; injury to surrounding tissue, nerves or organs
- Loss of neurological function, the appearance of radiculopathy, dura tears and / or pain. Neuro-vascular impairment including paralysis, temporary or permanent retrograde ejaculation in men or other severe injuries. Leakage of cerebrospinal fluid
- Inability to perform activities of daily life.
- Urinary retention or loss of bladder control or other types of impairment of the urologic system.
- Change in the curvature and the rigidity of the spinal column
- Failure to achieve the desired operation result, implant incorrectly positioned
- Partial loss of the degree of correction, which was achieved during surgery
- Superficial or deep infection and inflammation such as discitis, arachnoiditis, etc.

- Bone fracture above or below the treated segments
- Late bone growth or no visible fusion and pseudarthrosis
- Fracture, micro fracture, resorption, damage or penetration of a dorsal vertebra and / or bone graft and / or bone graft extraction location above or below the treated segment.
- Bleeding and / or haematomas
- Deep vein thrombosis, thromboembolism, pulmonary embolus
- Complications at the bone graft donor site
- Prolapsed intervertebral disc, intervertebral disc destruction or - degeneration on, above or below the treated segment
- Loss or increase in spinal mobility or function
- Damage to the reproductive system, sterility and sexual dysfunction.
- Development of respiratory problems such as pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.
- Change in mental state
- changed growth of the fused spinal column
- Allergic reaction to materials used
- Reduction of bone density, osteolysis
- Bursitis
- Death

PACKAGING, LABELLING AND STORAGE

- The implants are delivered in a STERILE condition. They may only be used if the outer packaging label as well as the inner packaging is intact. If the packaging is damaged or has already been opened, the sterility of the implant is not guaranteed and it may not be used.
- The implants may not be used if the specified expiration date is exceeded.
- The implant may not be re-sterilised.
- The implant components must be handled and stored with care. Damage to the implant can significantly reduce the stability and the fatigue strength of the implant system. Cracks and / or increased internal stresses can occur, which could result in rupture of the implant.
- The implants and instruments should be stored at room temperature. Environmental influences such as salt-laden air, moisture, chemicals, etc., may not act on the implants.
- A careful inspection prior to surgery is recommended to ensure that the instruments or implants were not damaged during storage or previous procedures.

IMPORTANT NOTE TO THE OPERATING PHYSICIAN

- The activities / movement behaviour of patients has a significant influence on the usage life of the implants. The patient must be informed that any activity increases the risk of loss, bending or rupture of the implant components. It is essential to inform the patient about restrictions in his activities in the postoperative period and to monitor the patient postoperatively to assess the development of the fusion and the condition of the implant. Even if a solid bone fusion has occurred, implant components can still bend, rupture or loosen. This is why the patient must be informed that implant components can bend, rupture or loosen, even if restrictions in activities are observed. In the event of complications, the doctor must decide whether the implant should be examined, taking into account the well-being of the patient and the risks involved.
- The instructions in the operation instructions (Surgical Technique) are to be observed. You can download these from the internet from www.humantech-solutions.de or directly from your HumanTech representative. It is advised to use the instruments described in the surgical technique, which are intended to be used with the system. No liability can be accepted for the use of external instruments.
- In the area of the spinal cord / spinal marrow and nerve roots, extreme caution is required as any nerve damage can result in loss of neurological function.
- Rupture, slippage, or incorrect use of instruments or implants can cause injury to the patient or surgical staff.
- When carrying out PEEK implants particular care in preparing the intervertebral disc space is needed to keep the placement forces low,
- For a better ongrowing, the PEEK-Ti implants are coated with a titanium surface. Take care during insertion of these cages as the surface is structured and rough.
- Bone cement may not be used as this material makes removal of the components difficult or impossible. The heat generated by the curing process, can damage or deform the PEEK implants.
- Removed implants must be treated in such a way that reuse is not possible.

PRODUCT COMPLAINTS

Any persons working in healthcare (e.g. customer or user of this product system), who has any kind of complaint or is dissatisfied in dealing with the product, in terms of quality, identity, durability, consistency, resistance, safety, effectiveness, and / or its function should inform the respective HumanTech representative.

The HumanTech representative should be immediately notified in the case of an ADONIS® implant ever displaying a "malfunction" (i.e. it does not meet the performance specifications or does not function as intended), or if the suspicion arises that this could become the case.

If a malfunction of a HumanTech product ever occurs that caused the death of or serious injury to a patient, the representative must be notified immediately telephonically, per fax or in writing. Should you have a complaint, we would like to ask you to send us, in writing, the name, the item number and the lot number of the component, as well as your name and address along with a detailed error description.

FURTHER INFORMATION

In the case of complaints, suggestions or information on the content of this Instruction for Use or the use of the product, please contact us at the address mentioned above.

© HUMANTECH Germany GmbH. All rights reserved.

	Manufacturer
	Do not re-use
	Batch
	Reference number
	Read the documentation
	Keep dry
	Do not re-sterilize
	Don't use when packing damaged
	Date of expiry
	Sterilisation with ethylene oxide
	Attention

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© HUMANTECH Germany GmbH. Všechna práva vyhrazena.