 	
ADONIS® LUMBAR INTERBODY FUSION	
HumanTech Spine GmbH Gewerbestraße 5 D-71144 Steinborn Tel: +49 (0) 7157/5246-71 Fax: +49 (0) 7157/5246-33 e-mail: info@humantech-spine.de www.humantech-spine.de	 0297

Stand: 15.10.2020

DE

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Revision: 15.10.2020

EN

IMPORTANT NOTICE
 This instruction for use must be read carefully and its content must be adhered.

SCOPE
 The instructions for use is valid for sterile delivered implants:
 - ADONIS® PEEK ...
 - ADONIS® R-PEEK-Ti...
 - ADONIS® Ti... sterile
 - ADONIS® LLIF Ti 3D ...
 and unsterile delivered implants.
 - ADONIS® Ti...

BASIC STRUCTURE
 The ADONIS® Lumbar Interbody Fusion System is used as an implant system for the fusion of vertebrae bodies in spine column surgery (cages). The cages are available in different dimensions, so that adaptations taking into account the unique anatomy of individual patients are possible.

INDICATIONS
 ADONIS® is intended for the surgical treatment of diseases and injuries of the lumbar and lumbosacral spine in patients whose general skeletal growth is completed, especially for indications such as:

- Prolapsed intervertebral disc
- Hard prolapsed intervertebral disc
- Mechanical instabilities
- Calcification of the posterior longitudinal ligament
- Osteochondrosis
- Lumbal spinal stenosis

MATERIAL
 There are different implant product lines.
 Basic material PEEK Optima®:

ADONIS® PEEK
 For products of the ADONIS® PEEK line, the base body is entirely made of PEEK Optima®. To achieve a better radiographic contrast, radiographic markers made out of titanium alloy Ti6Al4V respectively tantalum have been incorporated into the implant.

ADONIS® R-PEEK-Ti
 For products of the ADONIS® R-PEEK-Ti line, the base body is also made of PEEK Optima®. In addition, the surfaces of the implant have a titanium coating. To achieve a better radiographic contrast, radiographic markers made out of titanium alloy Ti6Al4V respectively tantalum have been incorporated into the implant.

Basic material Titanium alloy:
ADONIS® Ti und Ti...sterile und ADONIS® LLIF Ti 3D
 Products of the ADONIS® Ti line consist of the titanium alloy Ti6Al4V.

AREA OF APPLICATION
 The area of application is the lumbal and lumbosacral spine.

GENERAL CONDITIONS FOR USE

- Only surgeons with the required professional training in the field of spinal surgery may insert the implants. The use of the implants must be decided in keeping with the surgical and medical indications, the potential risks and the limitations relating to this type of surgical intervention, as well as with the knowledge of contraindications, side effects and known precautionary measures and the knowledge of the condition and the physical, chemical and biological properties of the implant.
- This instruction for use provides essential information, but is not sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a clearly

defined post-operative treatment regimen, and to conduct scheduled follow-up examinations.

- The activities and physical activity of the patient have a significant influence on the implant's useful life. Patients must be informed, that every activity increases the risk of loss, loosening, dislocation, migration, bending or breakage of implant or implant components. Informing patients about limitations to their activities in the postoperative phase and postoperatively monitoring patients are crucial factors in assessing the development of the surgical result and the condition of the implant.
- The above-mentioned effects may occur even if the implant is well integrated and the activity restrictions are complied with.
- In case of complications, the surgeon must decide whether a revision of the implant or other measures should be taken, taking into account the well-being of the patient and the risks involved.
- Components of the ADONIS® system may not be replaced by components/products from other systems, from another source or from a different manufacturer. Furthermore, no direct connection of components/ products of the systems to components of other systems may be established. If this recommendation is not complied with or if the products are otherwise used or used improperly, HumanTech Spine GmbH assumes no liability responsibility.
- Following the notes in the operating instructions (surgical technique) is essential. You will find these in the Internet at www.humantech-spine.de or you can obtain them directly from your HumanTech representative. It is advisable to use the instruments described there and provided by the manufacturer. There may be complications if the implant is used with or without the dedicated instrument. For the use of foreign instruments, no liability can be accepted.
- Proceed with extreme caution in the region of the spinal cord and the roots of the nerves, since damage to the nerves can lead to the impairment of neurological functions.
- Breakage, slippage or incorrect use of the instruments or implants can injure the patient or the operating staff.
- Do not use bone cement, as this material makes the removal of the components difficult or impossible. The heat produced by the hardening process can damage or deform the PEEK implants.
- Residues consisting of implant material and/ or not from implant material should be removed.
- Damaged implants must not be implanted.
- Never reuse the implants. Even if the implant appears to be intact after revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break.
- Implants that have already come into contact with a patient's body fluids or tissues or have been soiled must not be reused.
- Handle removed implants in such a way that their reuse is not possible.
- Safety and compatibility of the device in the setting of magnetic resonance (imaging) have not been evaluated. No normal test or test of migration has been performed on the device in this setting.
- When using PEEK implants particular care in preparing the intervertebral disc space is needed to keep the placement forces low
- For a better ongrowing, the R-PEEK-Ti and 3D printed implants are coated with a titanium surface. Take care during insertion of these cages as the surface is structured and rough.

CONTRAINDICATIONS
 Contraindications may be either relative or absolute. The selection of a particular implant must be weighed carefully against the overall assessment of the patient. The following conditions can have an adverse impact on the chances of successful surgery:

- all destructive and inflammatory diseases of the vertebral body or movement segment such as tumours, spondylitis and fractures with a high degree of instability, with a strong static and structural change of the vertebral body and segment structure, which require vertebral body replacement and other stabilising measures for reconstruction and stabilisation of the segment acute infections or significant risks of infections (weakened immune system)
- indications of local inflammation
- fever or leucocytosis
- pathological obesity
- pregnancy
- psychological disorders
- any other medical or surgical condition hindering possible improvement through the use of the implant, such as the existence of congenital abnormalities, fracture in the vicinity of the operating point, elevated deposition rate not indicated by other clinical pictures, elevated leucocyte count (white blood cells), or distinct leftward shift in the white blood cell differential blood count.
- Diseases of the joints, bone resorption, osteopenia, osteomalacia and/or osteoporosis are relative contraindications, since these can limit the degree of correction and stabilization that can be achieved.
- Any neuromuscular disease that would place excess strain on the implant during the healing period.
- known hereditary or acquired brittleness of the bones or calcification problems
- spondylolisthesis that cannot be reduced to the first degree

- Suspicion of a metal allergy or intolerance, and documented metal allergy or intolerance. Appropriate tests should be carried out
- all cases in which the use of components of different metals or alloys is necessary
- already existing fusion in the segments treated
- all cases in which fusion is not necessary
- all cases in which the implant component selected for use is too large or too small to achieve a satisfactory result
- all patients with inadequate tissue structure on the operative side or inadequate bone or bone quality
- all patients in which the use of the implant would interfere with anatomical structures or restrict physiological performance
- all patients not willing to follow postoperative instructions
- all cases not described in the indications

SIDE EFFECTS / COMPLICATIONS
 The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

- loosening of the implant
- shifting of the implant
- breakage of the implant, respectively detached parts of the implant and wear
- foreign body reactions on the implant, including possible tumor formation, autoimmune disease and/or scar formation
- scar formation with possible neurological impairments, nerve compression or pain
- neurological complications
- paralysis
- soft tissue lesions
- neurological or spinal lesions of the dura mater as a result of surgical trauma
- The migration of the implant in the posterior direction can cause a loss of neurological function, the occurrence of radiculopathy, dural tears and/or pain.
- Neurovascular impairment, including paralysis
- temporary or permanent retrograde ejaculation in men, or other severe injuries
- Discharge of cerebrospinal fluid
- inability to perform daily tasks
- urinary retention or loss of control over the bladder, or other types of urinary system impairment
- alteration of the curvature and stiffness of the spine
- Failure to achieve the desired result of surgery, malpositioning of the implant
- partial loss of the degree of correction achieved during the operation
- superficial or deep-seated infections and inflammations, such as discitis or arachnoiditis
- bone fracture above or below the segments treated
- late osseointegration or no visible fusion and pseudarthrosis
- fracture, micro fracture, resorption, damage to or penetration of vertebral bone and/or the bone transplant and/or the bone transplant donor site above and below the segment treated
- adjacent segment disease
- bleeding and/or hematomas
- deep-seated venous thrombosis, thrombophlebitis, pulmonary embolism
- complications at the bone transplant donor site
- herniated disk, disk destruction or disk degeneration at, above or below the segment treated
- loss or enhancement of spinal mobility or functions
- damage to the reproductive system, sterility and sexual functional disturbances
- development of respiratory problems, for example pulmonary embolism, atelectasis, bronchitis, or pneumonia
- alteration of the mental state
- alteration in the growth of the fused spine
- allergic reaction to the materials used
- decrease in bone density
- bursitis
- death

PACKAGING, LABELING, TRANSPORT AND STORAGE

- The implant components must be handled, transported and stored with due care. Damages to the implant packaging or damages to the implant itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increase internal loads, which can result in a fracture of the implant.
- The implants and instruments should be stored at room temperature. Environmental influences such as salty air, humidity, chemicals etc. may not be allowed to act on the implants.
- The implants and instruments should be stored at room temperature. Environmental influences such as salty air, humidity, chemicals etc. may not be allowed to act on the implants.

The implants are delivered in both NON STERILE and STERILE packaged condition:
Unsterile delivered implants:
 UNSTERILE delivered implants are labeled as UNSTERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILISATION).
 The implants are supplied as an implant system in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

Sterile delivered implants:

STERILE delivered implants are sterilized by a validated gamma sterilization or ethylene oxide procedure and labeled as STERILE. Cleaning, disinfection and sterilization prior to use need not be performed.
 The implants are supplied individually packaged in a blister and protected by a covering box. The implants may be used only when the label on the outer packaging and also the inner packaging are intact. If the packaging is damaged or already open, the sterility of the implant is not guaranteed and the implant may not be used.
 The implants may not be used when the shelf life indicated has been exceeded.
 Processing, reprocessing, sterilization or re-sterilization of the products after opening the sterile packaging or with damaged sterile packaging is not intended.
 HumanTech Spine GmbH assumes no responsibility for the use of re-sterilized implants regardless of the person who performed the re-sterilization or the method used.
 The rules of asepsis must be observed when removing the implant from the sterile packaging. The sterile packaging may only be opened immediately before the implant is inserted. It is recommended to always have a replacement implant available. The implant must be removed from the packaging using appropriate aseptic provisos.

CLEANING, DISINFECTION AND STERILISATION
 UNSTERILE packaged delivered products of the ADONIS® Lumbar Interbody Fusion System must be cleaned, disinfected and sterilized prior to use. All necessary steps for cleaning, disinfection, maintenance and sterilization are described in the instruction "Processing unsterile delivered Spine". The latest version can always be found on our homepage <https://www.humantech-spine.de/378-EN-FUS.html>, as well as on request from HumanTech Spine GmbH.

DISPOSAL
 The products must be disposed of in accordance with the applicable local regulations, whereby the respective degree of contamination must be taken into account.

PRODUCT COMPLAINTS
 Any person involved in healthcare (e.g. customer or user of this product system) who has complaints of any kind or is dissatisfied using the product with respect to quality, identity, strength, durability, safety, effectiveness or function should notify the appropriate HumanTech representative. If a ADONIS® implant should ever malfunction (i.e. fails to fulfill the performance specifications or does not function in the foreseen manner), or if the above is suspected, the HumanTech representative should be informed immediately. If a HumanTech product should ever malfunction in a way that results in the death or serious injury of the patient or if it has contributed to the above, the representative is to be informed immediately by telephone, fax or in written form.

If you have any complaints then we ask you to state the name, article number and batch number of the components as well as your name and address along with a detailed description of the fault and forward this information to us in writing.

FURTHER INFORMATION
 If you have any complaints, suggestions or other points about the contents of this user manual or the use of the product, please contact the above address.

© HUMANTECH Spine GmbH. All rights reserved.

	Manufacturer
	Do not re-use
	Batch number
	Reference number
	Read the documentation
	Keep dry
	Unsterile
	sterilized by irradiation
	Sterilisation with ethylene oxide
	Do not re-sterilize
	Date of expiry

	Do not use when packaging is damaged
	Attention

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