	
SAMSON® VERTEBRAL BODY REPLACEMENT	
 <p>HumanTech Spine GmbH Gewerbestraße 5 D-71144 Steinbronn</p> <p>Tel: +49 (0) 7157/5246-71 Fax: +49 (0) 7157/5246-33 e-mail: info@humantech-spine.de www.humantech-spine.de</p>	

Stand: 25.10.2019

DE

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SAMSON® VERTEBRAL BODY REPLACEMENT

HumanTech Spine GmbH
Gewerbestraße 5
D-71144 Steinbronn
Tel: +49 (0) 7157/5246-71
Fax: +49 (0) 7157/5246-33
e-mail: info@humantech-spine.de
www.humantech-spine.de

Last updated: 25.10.2019

EN

IMPORTANT INFORMATION

This instruction for use must be read carefully and its content must be adhered to.

SCOPE

The instruction for use is applicable for all sterile and unsterile delivered implants of the SAMSON®-System.

BASIC STRUCTURE

The SAMSON® Vertebral Body Replacement System (VBR) is used as an implant system for the replacement of one or more vertebral bodies in spinal surgery. The VBR system offers a range of different dimensions, making it possible to perform adaptations which take into account the unique anatomy of each individual patient.

MATERIAL

All components of the SAMSON® vertebral body replacement system are made of titanium alloy Ti6Al4V according to ISO 5832-3.

INDICATIONS FOR USE

The system is in cases of cancerous, inflammatory and traumatic diseases which lead to instabilities around the anterior support, or compression of the neural structures, or for diseases which require treatment of an infection. The SAMSON® system has been designed for use with an additional dorsal fixation system. Implants of SAMSON® system are not intended for stand-alone use. SAMSON® is intended for patients whose general skeletal growth is finished. When using a SAMSON implant base bodies Size 5 or Size 6, or when using 12° angled endplates mounted at the caudal and cranial position on the implant base body, ventral support (such as the VENUStano®-Screw Rod System) is strongly advised.

GENERAL CONDITIONS/ INFORMATIONS FOR USE

- The implants must be implanted only by surgeons having undergone the necessary training in spinal surgery. The insertion of implants must be carried out in line with the surgical and medical indications, the potential dangers, and with the restrictions associated with this type of surgery. It should also be performed with an awareness of the contraindications, side effects, defined precautionary measures, and with regard for the structure and the physical, chemical and biological characteristics of the implant.
- The activities and physical activity of the patient have a significant influence on the implant's useful life. Patients must be informed that every activity increases the risk of loss, loosening, dislocation, migration, bending or breakage of implant components. Informing patients about limitations to their activities in the postoperative phase and postoperatively monitoring patients are crucial factors in assessing the development of the surgical result and the condition of the implant. The above-mentioned effects may occur even if the implant is well integrated and the activity restrictions are complied with. The patient must be informed of this. If such an effect occurs, the doctor must decide whether a revision of the implant or other measures should be taken, taking into account the well-being of the patient and the risks involved.
- This instruction for use provides essential information, but is not sufficient for the use of the system. This information is not a substitute for the following: the professional judgement and/or clinical skills and experience of the physician with regard to careful patient selection; preoperative planning and implant selection; knowledge of the anatomy and biomechanics of the spine; an understanding of the material and the mechanical properties of the used implants; training and skills in spinal surgery and the use of the instruments required for inserting the implants; the surgeon's ability to gain the patient's consent, adhere to a clearly

defined post-operative treatment regimen, and to conduct scheduled follow-up examinations.

- It is recommended not to replace components of SAMSON® system with components or products of other systems, or make a direct connection to components of other systems. HumanTech Spine GmbH shall accept no liability if this recommendation is not observed.
- Never reuse the implants. Even if the implant appears to be intact following revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break. Handle removed implants in such a way that their reuse is not possible.
- Implants must also not be re-used if they have come into contact with the body fluids or tissues of a third person. Cleanliness when processing can otherwise not be guaranteed.
- Handle removed implants in such a way that their reuse is not possible.
- Care should be taken to maintain the integrity of the endplates of the vertebral bodies. Damage to the endplates of the vertebral bodies or excessive partial abrasion of the endplates of the vertebral bodies can lead to sintering of the implant and to a loss of segmental stability.
- The endplates of the implant should be selected in a sufficiently large size in order to make use of the maximum contact surface of the vertebral body. This will ensure the greatest possible stability and counteract implant sintering. However, in order to avoid injury to the adjacent structures, the endplates of the implant should not protrude beyond the vertebral body.
- The implant must be implanted under prestress in order to avoid dislocation. Expansion of the implant must be performed under continuous x-ray control in order to avoid over-distraction, damage to the implant endplates or other implant components, e.g. of the posterior fixation.
- The instructions in the operation instructions (Surgical Technique) are to be observed. You can download these from the internet from www.humantech-spine.de or directly from your HumanTech representative. It is advised to use the instruments described in the surgical technique, which are intended to be used with the system. Complications may occur if the implant is inserted with or without the intended instruments. No liability can be accepted for the use of external instruments.
- Proceed with extreme caution in the region of the spinal cord and the roots of the nerves, since damage to the nerves can lead to the impairment of neurological functions.
- Breakage, slippage or incorrect use of the instruments or implants can injure the patient or the operating staff or result in more time being required for surgery.
- Bone cement must not be used, as the removal of the components becomes difficult or impossible.
- Residues consisting of implant material and/ or not from implant material should be removed.
- Damaged implants must not be implanted.
- Safety and compatibility of the device in the setting of magnetic resonance (imaging) have not been evaluated. No thermal test or test of migration has been performed on the device in this setting.
- The endplates of the implant must be assembled correctly to the preassembled implant base body.
- Careful inspection of the assembled SAMSON implants is recommended before and during surgery to ensure that the implants are undamaged during assembly or previous procedures and have been correctly assembled according to Surgical Technique.

CONTRAINDICATIONS

Contraindications may be either relative or absolute. The selection of a particular implant must be weighed carefully against the overall assessment of the patient. The following conditions can have an adverse impact on the chances of successful surgery:

- Extreme malalignments, impairing the stability of the instrumentation
- Suspicion of a metal allergy or intolerance, and documented metal allergy or intolerance. Appropriate tests should be carried out.
- All cases in which the selected implant components are too large or too small to achieve a satisfactory result
- Every patient in which the use of the implant would interfere with anatomical structures
- Any neuromuscular disease that would place excess strain on the implant during the healing period.
- Any patient who is not willing or able to follow post-operative instructions
- Morbid obesity
- Pregnancy
- Use in combination with implants made of different metals or alloys
- All cases which are not described in the indications
- Any other medical or surgical condition which prevents possible improvement through the use of the implant

SIDE EFFECTS / COMPLICATIONS

The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

- Missing intact connection segments, bone resorption, osteopenia, osteomalacia and/or osteoporosis, as well as a diffuse spinal tumours and infections, can, in addition to other causes, lead to the reduction of instrumentation stability, limit the degree of obtainable correction and stabilisation or increase the risk of sintering as well as implant loosening or migration.
 - Failure to achieve the desired result of surgery, malpositioning of the implant
 - Fracture, micro fracture, resorption, damage to or penetration of a vertebral body above or below the treated segment/s
 - Reduction of bone density
 - Change in the curvature and stiffness of the spine
 - Loss or increase in spinal mobility or functions
 - Physiological limitations, such as joint degeneration
 - Bending, fracture or loosening of the implant components, as well as abrasion
 - Delayed bone growth or no visible regeneration and pseudarthrosis, delayed healing or lack of healing
 - Neurological deterioration with paralysis, functional failures, loss of sensation and pain, e.g. caused by scarring
 - Soft-tissue lesions
 - Pressure on surrounding tissue or organs
 - The occurrence of laceration of the dura with risk of cerebrospinal fluid loss or a CSF fistula
 - Radiculopathy
 - Superficial or deep infection and inflammation as discitis, arachnoiditis, etc.
 - Inability to perform one's daily tasks
 - Dysfunctions of the bladder and rectum
 - Temporary or permanent retrograde ejaculation in men
 - Damage to the reproductive system, sterility and sexual dysfunctions
 - Deep vein thrombosis, thrombophlebitis, pulmonary embolus
 - Development of respiratory problems, including pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.
 - Lung injuries
 - Injuries of the large abdominal vessels
 - Bleeding and/or haematomas
 - Allergic reaction to the materials used
 - Foreign body reaction to the implants, including possible tumour formation, autoimmune disease and/or scarring
 - Metallosis
 - Bursitis
 - Death
- as well as all common surgical risks.

PACKAGING, LABELLING AND STORAGE

- Handle and store the implant components carefully. Damage to the implant can significantly reduce the strength and long-term stability of the implant system. Cracking and/or higher internal stresses can result, possibly resulting in the breakage of the implant.
- The implants and instruments should be stored at room temperature. Environmental factors, such as salt-laden air, humidity, chemicals, etc., must not be allowed to act on the implants.
- Thorough inspection is recommended before operating in order to ensure that the instruments or implants have not been damaged during storage or previous procedures.

The implants are delivered in both NON STERILE and STERILE packaged condition:

Unsterile delivered implants:

UNSTERILE delivered implants are labeled as UNSTERILE and must be cleaned and sterilized before use (see below REPROCESSING and STERILISATION).

The implants are supplied as a complete set or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilisation in the original pouch packaging is not allowed.

Sterile delivered implants:

STERILE delivered implants are sterilized by a validated gamma sterilization procedure and labeled as STERILE. Cleaning, processing and sterilization prior to use need not be performed.

The implants are supplied individually packaged in a blister-blister-combination and protected by a covering box. The implants may be used only when the label on the outer packaging and also the inner packaging are intact. If the packaging is damaged or already open, the sterility of the implant is not guaranteed and the implant may not be used.

The implants may not be used when the shelf life indicated has been exceeded. Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or with damaged sterile packaging is not intended.

REPROCESSING

Unsterile delivered products of the SAMSON® Vertebral Body Replacement System must be processed prior to sterilisation. The necessary steps are described in the reprocessing instructions for Ti implants.

STERILISATION

Unsterile delivered products of the SAMSON® Vertebral Body Replacement System must be sterilised or re-sterilised prior to use by means of validated steam sterilisation in an autoclave. The minimum recommended sterilisation parameters are as follows:

Steam	Pre-vacuum	134°C	5 minutes
-------	------------	-------	-----------

PRODUCT COMPLAINTS

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a SAMSON® implant ever exhibits "malfunctioning" (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, the representative should be informed immediately by phone, fax, or in writing.

Should you have a complaint, please provide the name, item number and batch number of the component along with your name, address and detailed description of the defect in writing.

FURTHER INFORMATION

In the event of complaints, suggestions or indications regarding the content of these instructions for use or the use of the product, please refer to the above-mentioned address.

© HUMANTECH Spine GmbH. All rights reserved.

	Manufacturer
	Do not re-use
	Batch number
	Reference number
	Read the documentation
	Keep dry
	Unsterile
	sterilized by irradiation
	Do not re-sterilize
	Date of expiry
	Do not use when packaging is damaged.
	Attention

	
SAMSON® VERTEBRAL BODY REPLACEMENT	
 <p>HumanTech Spine GmbH Gewerbestraße 5 D-71144 Steinbronn Tel: +49 (0) 7157/5246-71 Fax: +49 (0) 7157/5246-33 e-mail: info@humantech-spine.de www.humantech-spine.de</p>	

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© HUMANTECH Spine GmbH. Todos los derechos reservados.