# SAMSON®

## VERTEBRAL BODY REPLACEMENT

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DE

2 - 7

Instruction for use

EN

8 - 12

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## IMPORTANT INFORMATION

This instruction for use must be read carefully and its content must be adhered to.

## SCOPE

The instruction for use is applicable for all sterile and unsterile delivered implants of the SAMSON®-System.

## INTENDED USE

The SAMSON® Vertebral Body Replacement System (VBR) is an implant system for long-term use for anterior stabilization of the thoracic and lumbar spine as a replacement for one or more vertebral bodies.

Devices of the SAMSON® system are intended for use with an additional dorsal fixation system (VENUS®). Devices of the SAMSON® system are not intended for stand-alone use.

When using SAMSON implant base bodies Size 5 or Size 6, or when using 12° angled endplates mounted at the caudal and cranial position on the implant base body, ventral support (such as the VENUSnano®-Screw Rod System) is intended to be used.

## GENERAL PRODUCT DESCRIPTION

The SAMSON® Vertebral Body Replacement System consists of implant base bodies and implant endplates in various dimensions. The implant base bodies are available in different heights for different defect heights. The height of the implant base bodies can intraoperatively be adapted by an expansion mechanism to take into account the unique anatomy of the individual patient.

The implant endplates are available in different widths and angles for adaptation to different vertebral body structures and for mapping the existing or achievable lordotic or kyphotic curvature of the spine. The unique anatomy of the individual patient can be taken into account by the possibility of connecting the different implant components with each other via plate holding screws.

## INDICATIONS FOR USE

Devices of the SAMSON® system are intended for the surgical treatment of tumorous, inflammatory, and traumatic diseases and injuries of the upper thoracic to the lower lumbar spine, which lead to instabilities in the area of the anterior column or compression of the neural structures or for diseases requiring remediation as a result of an inflammation.

## CONTRAINDICATIONS

Contraindications may be either relative or absolute. The selection of a particular implant must be weighed carefully against the overall assessment of the patient.

Contraindications are:

- Acute infections or significant risks of infections (weakened immune system)
- Fever or leucocytosis
- Missing intact adjacent segments, joint diseases, bone resorption, osteopenia, osteomalacia and/or osteoporosis are relative contraindications, as these can limit the degree of achievable correction/stabilization
- Extreme malalignments, impairing the stability of the instrumentation
- Suspected allergy or intolerance, as well as documented allergy or intolerance to the material used (Ti6Al4V). Appropriate test must be carried out.
- All cases where a combined use of implants made of different metals or alloys or from different sources is required.
- Any patient who is not willing or able to follow postoperative instruction
- Morbid obesity
- Pregnancy

The following conditions can have an adverse impact on the chances of successful surgery:

- Wound healing disorders
- Any patient with inadequate tissue structure at the surgical site
- All cases in which the selected implant components are too large or too small to achieve a satisfactory result
- Every patient in whom the use of the implant would interfere with anatomical structures
- Any neuromuscular disease that would place an exceptionally high load on the implant during the healing period.

## INTENDED PATIENT GROUP

All patients that are at least 18 years old and skeletally matured at the time of surgery.

## INTENDED USERS

The devices of the SAMSON® system may only be implanted by surgeons who have completed the necessary training in spinal surgery.

The

- professional judgement and/or clinical skills and experience of the surgeon with regard to careful patient selection
- preoperative planning and implant selection
- knowledge of the anatomy and biomechanics of the spine

- understanding of the material and mechanical properties of the used implants
- training and skills in spinal surgery and the use of the instruments required for inserting the implants
- surgeon's ability to gain the patient's consent to adhere to a clearly defined postoperative treatment regimen and to conduct scheduled follow-up examinations

is/are necessary for the use of the devices of the system. The surgeons are supported in the operating theater environment by surgical assistants. General knowledge in surgical procedures, the surgical technique and special training and skills regarding the correct use of the devices of the system and the use of the instruments required for inserting the implants is/are necessary for the use of the devices.

### INTENDED ENVIRONMENT

The devices have to be used in an operating theater. A sterile environment and compliance with the general hygiene regulations that apply in operating theaters are basic requirements for using the devices.

It must be assured that the surgical field is fully illuminated with a view of the situs.

The implantation of the product in the patient takes place under general anesthesia. In order to be able to perform the operation under general anesthesia, anesthetic gas or other anesthetics are required.

It has to be assured that C-Arm (X-ray) and the system-specific instruments are available. Furthermore, general surgical instruments are required for the performance of the surgery.

### SIDE EFFECTS / ADVERSE EVENTS

The side effects and complications listed are not only due to the implants, but often also to the surgical procedure and include, but are not limited to:

- Implant loosening
- Implant migration
- Implant subsidence
- Implant breakage
- Malpositioning of the implant
- Instability of the instrumentation
- Bending, fracture or loosening of the implant components
- Expansion of the implant not possible
- Signs of abrasion, for example occurrence of microparticles in the area of the implant (metallosis)
- Revision surgery or requirement of further surgeries
- Missing intact connection segments, bone resorption, osteopenia, osteomalacia and/or osteoporosis, as well as diffuse spinal tumours and infections, can, in addition to other causes, lead to the reduction of instrumentation stability, limit the degree of obtainable correction and stabilisation or increase the risk of sintering as well as implant loosening or migration.
- Changes in the curvature and stiffness of the spine

- Loss or increase of spinal mobility or functions
- Partial loss of the degree of correction achieved during the surgery
- Physiological restrictions, such as negative impairment or degeneration of adjacent segment(s)
- Delayed bone growth or no visible fusion and pseudarthrosis, delayed healing or lack of healing
- Herniated disc, disc destruction or degeneration in, above or below the treated segment(s)
- Changed growth of the fused spine
- Fracture, micro fracture, resorption, damage to, necrosis or penetration of a vertebral body above or below the treated segment(s)
- Reduction of bone density
- Neurological deterioration with numbness, paralysis, functional failures, loss of sensitivity, reflex changes and pain, e.g. caused by scarring or durotomy
- Neurological complications like injury or irritation of the spinal cord, the root nerves (radiculopathy) or peripheral nerves
- Soft-tissue complications like lesions or heterotopic ossification
- Pressure on or injury of surrounding tissue, nerves or organs
- Neurological or spinal lesion of the dura mater due to surgical trauma with risk of loss of cerebrospinal fluid (CSF) or CSF fistula
- Injury of the lymphatic vessels/ lymph leakage
- Bleeding and/or haematomas
- Injuries of the large abdominal and/or other blood vessels
- Development of respiratory problems, including atelectasis, bronchitis, pneumonia, pleural fluid collections, pneumothorax, hemothorax, etc.
- Lung injuries
- Urinary retention or loss of bladder control, or other types of impairment of the urological system; dysfunction of the rectum
- Temporary or permanent retrograde ejaculation in men
- Damage to the reproductive system, sterility and sexual dysfunctions
- Superficial or deep infections and inflammation as discitis, sepsis, arachnoiditis, etc.
- Inability to perform daily tasks
- Deep vein thrombosis, thrombophlebitis, pulmonary embolus
- Change of mental state
- Allergic reaction to the materials used
- Foreign body reaction to the implants, including possible tumour formation, autoimmune disease and or scar formation
- Failure to achieve the desired result of surgery
- Back or leg cramps
- Headache
- Insomnia
- Myocardial infarction
- Bursitis
- Death

as well as all general surgical risks, like high blood pressure, cardiopulmonary complications, delirium, coagulation disorder, coagulopathy, incisional hernia, blunt trauma and fever.

## WARNINGS AND PRECAUTIONS

### General

The user has to ensure that the most current versions of the complete product materials are on hand and considered. These are also available at: [www.humantech-spine.de](http://www.humantech-spine.de)

### MR safety information

SAMSON® Vertebral Body Replacement Implants have been examined in non-clinical-testings in accordance with ASTM norms and determined as “MR conditional”. Patients with SAMSON® implants can be safely examined via MRI on following conditions:

- Static magnetic field of 1.5 T to 3 T
- Spatial gradient magnetic field up to max. 1100 Gauss/cm or 11 T/m
- Whole body specific absorption rate (SAR) up to max. 2 W/kg in normal operating mode and at max. 15 min continuous scan duration

SAMSON® Vertebral Body Replacement implant can heat up to 4,7 °C with continuous scanning time of 15 min. and under the scanning conditions defined above. Under these conditions, a patient can be examined at low risk. To keep the risks of warming low, the scan duration should be kept as short as possible and the SAR as low as possible.

The product has not been examined for artifacts in MR environment. MR image quality may be reduced if the examined area is in relative proximity to a SAMSON® Vertebral Body Replacement implant.

A thorough risk-benefit assessment by the attending doctor must be carried out.

### Preoperative

- The use of implants for implantation must be decided upon in accordance with the surgical and medical indications, the potential risks and restrictions associated with this type of surgical procedure, as well as in knowledge of the contraindications, side effects and defined precautions, and in awareness of both the nature and the physical, chemical and biological characteristics of the implant
- The correct selection and placement of implants has to be ensured or to be checked based on suitable pre-, intra- and postoperative diagnostic procedures.
- The activities / movement behaviour of the patient has a significant influence on the service life of the implant. The patient must be informed that any activity increases the risk of loss, bending or fracture of the implant components. It is crucial to inform the patient about restrictions in activities in the post-operative phase and to monitor the patient post-

operatively in order to assess the development of the fusion and the condition of the implant. Even when a solid bone fusion has occurred, implant components may still bend, break or loosen. Therefore, the patient must be informed that implant components can also bend, break or loosen if the restrictions in activities are properly followed

- The patient has to be informed in detail about the risks of the procedure and the use of the implants, including any revisions that may be necessary.
- Errors in implant selection can lead to premature clinical implant failure. The number of segments to be served is to be determined carefully

### Intraoperative

- The instructions in the surgical technique must be adhered to. You can download these at [www.humantech-spine.de](http://www.humantech-spine.de) or obtain them directly from your HumanTech representative. It is advised to use the instruments described in the surgical technique, which are intended to be used with the system. No liability can be assumed for the use of external/third-party instruments
- Damaged implants must not be used
- Components of the SAMSON® system may not be replaced by components / products from other systems from another source or from a different manufacturer. Furthermore, no direct connection of components / products of the system to components of other systems may be established. If this is not complied with or if the products are otherwise used or used improperly, HumanTech Spine GmbH assumes no responsibility
- Breakage, slippage or incorrect use of the instruments or implants can result in an injury of the patient or the surgical personell or result in an extension of the duration of the surgery
- Care should be taken to maintain the integrity of the endplates of the vertebral bodies. Damage to the endplates of the vertebral bodies or excessive partial abrasion of the endplates of the vertebral bodies can lead to sintering of the implant and to a loss of segmental stability
- The endplates of the implant should be selected in a sufficiently large size in order to make use of the maximum contact surface of the vertebral body. This will ensure the greatest possible stability and counteract implant sintering. However, in order to avoid injury to the adjacent structures, the endplates of the implant should not protrude beyond the vertebral body
- The endplates of the implant must be assembled correctly to the preassembled implant base body
- Careful inspection of the assembled SAMSON implants is recommended before and during surgery to ensure that the implants are undamaged during assembly or previous procedures and have been correctly assembled according to the instructions in

- the surgical technique
- The implant must be implanted under pre-stress in order to avoid dislocation. Expansion of the implant must be performed under continuous x-ray control in order to avoid over-distraction, damage to the implant endplates or other implant components, e.g. the posterior fixation
- Proceed with extreme caution in the region of the spinal cord and the nerve roots as any damage to the nerves can result in the negative impairment or loss of neurological functions
- Should an implant component or parts of an instrument loosen or break in situ, these residues must be completely removed
- Implants that have come into contact with a patient's body fluids or tissues or have been soiled must not be reused
- Bone cement must not be used, as the removal of the components becomes difficult or impossible

### Postoperative

- The patient should be instructed to inform his treating physician immediately regarding any unusual changes in the surgical area
- The patient should be monitored if a change in the area where the implant(s) is/are located is noticed. The treating physician should assess the possibility of clinical implant failure and discuss with the patient the actions required to achieve further healing
- In case of complications, the surgeon has to decide whether a revision of the implant should be carried out, taking into account the patient's condition and the possible risks involved
- Never reuse the implants. Even if the implant appears to be intact after the revision, alterations within the implant or minute defects resulting from the loading and stressing to which the implant has been subjected can cause the implant to break
- Handle removed implants in such a way that their reuse is not possible
- An implant pass must be given to the patient

### MATERIAL

All components of the SAMSON® vertebral body replacement system including all components of preassembled assemblies are made of titanium alloy Ti6Al4V according to ISO 5832-3.

### PACKAGING, LABELING, TRANSPORT AND STORAGE

- The devices must be handled, transported and stored with care. Damages to the device packaging or damages to the device itself can significantly reduce the performance, stability and durability of the implant system. It can lead to cracks and/or increased internal loads, which can result in a fracture of the device

- The implants and corresponding instruments should be stored at room temperature. Environmental factors such as salty air, humidity, chemicals, etc., must not be allowed to act on the devices
- Before the surgery, a careful inspection of the components of the SAMSON® Vertebral Body Replacement system to be used must be carried out in order to avoid damages caused by storage, transport or previous procedures. The sterilization containers, trays and the associated covers must not be damaged

The devices are delivered in both NON STERILE and STERILE packaged condition:

#### Unsterile delivered devices:

UNSTERILE delivered devices are labelled as UNSTERILE and must be cleaned, disinfected and sterilized before use (see CLEANING, DISINFECTION and STERILISATION).

The devices are supplied as an implant system in equipped trays in sterilization containers or individually packaged. Individual packaging must be unopened at the time of delivery. Sterilization in the original packaging is not allowed. The sterilization containers, trays and the associated covers must not be damaged.

#### Sterile delivered devices:

STERILE delivered devices are sterilized by a validated gamma sterilization procedure and labelled as STERILE. Cleaning, disinfection and sterilization prior to use need not be performed.

The devices are supplied packaged in a blister-blister-combination and protected by an outer carton. The devices may only be used when the label on the outer packaging and also the inner packaging are intact. If the packaging is damaged or already open, the sterility of the device is not guaranteed, and the device must not be used.

The devices must not be used if the indicated expiry date has been exceeded. Processing, reprocessing, sterilization or resterilization of the products after opening the sterile packaging or if the sterile packaging is damaged is not intended.

HumanTech Spine GmbH assumes no responsibility for the use of re-sterilized devices regardless of the person who performed the resterilization or the method used.

The rules of asepsis must be observed when removing the devices from the sterile packaging. The sterile packaging may only be opened immediately prior to the insertion of the device. It is recommended that a replacement device is always kept available. The implant must be removed from the packaging using appropriate aseptic precautions.

### CLEANING, DISINFECTION AND STERILISATION

UNSTERILE packaged delivered products of the SAMSON® Vertebral Body Replacement System must be cleaned, disinfected and sterilized prior to use. All

necessary steps for cleaning, disinfection, maintenance and sterilization are described in the instruction “Processing instructions for unsterile delivered medical devices Spine”. The latest version can always be found on our homepage <https://www.humantech-spine.de/378-en-IFUs.html>, as well as on request from HumanTech Spine GmbH.

The unsterile packaged delivered preassembled implant basebodies of the SAMSON® Vertebral Body Replacement System are not intended to be disassembled prior to processing.

## DISPOSAL

Any HumanTech implant that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol.

Devices must be disposed of as a healthcare medical device in accordance with hospital procedures.

## PRODUCT COMPLAINTS

Any person operating within the healthcare system (e.g. customer or user of this product system) who has complaints of any kind or who is dissatisfied when using the product in terms of quality, identity, durability, stability, safety, effectiveness and/or function, should inform the relevant HumanTech representative.

If a SAMSON® implant ever exhibits “malfunctioning” (i.e. does not correspond to the performance specifications or does not work as intended), or something occurs which would suggest this, the HumanTech representative should be informed immediately.

If a HumanTech product ever exhibits a defect which has caused or contributed to the death or serious injury of a patient, the representative should be informed immediately by phone, fax, or in writing.



Should you have a complaint, please provide the name, item number and batch number of the component along with your name, address, and detailed description of the defect in writing.

## FURTHER INFORMATION

In the event of complaints, suggestions, or indications regarding the content of these instructions for use or the use of the product, please refer to the above-mentioned address.

© HUMANTECH Spine GmbH. All rights reserved.

## EXPLANATION OF SYMBOLS

	Manufacturer
	Do not re-use
	Batch number
	Reference number
	Quantity
	Medical device
	Unique Device Identification
	Read the documentation
	Keep dry
	Unsterile
	sterilized by irradiation
	Do not re-sterilize
	Date of expiry
	Do not use when packaging is damaged.
	Attention
	MR conditional

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