HERO®

CERVICAL PLATE SYSTEM

HumanTech Germany GmbH
Gewerbestraße 5
D-71144 Steinenbronn
Tel: +49 (0) 7157/5246-71
Fax: +49 (0) 7157/5246-33
e-mail: info@humantech-solutions.com
www.humantech-solutions.com

Stand: 20.07.2011

DE

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EN

IMPORTANT NOTICE

The users confirm that they have read and accepted the conditions of this product insert as a contractual document.

BASIC DESIGN

The HERO® Cervical Plate System, which is used as an implant system for fixing in spinal surgery, consists of plates and screws of various sizes making adjustments possible, taking into account the unique pathology of the individual patient.

MATERIAL

All components are made of Ti6Al4V ELI, a titanium alloy in compliance with ASTM F136-98 and ISO 5832-3.

INDICATIONS FOR USE

The Hero® Cervical Plate System is designed for the surgical treatment of diseases and injuries of the cervical spine from C2 to C7 in patients in whom general skeletal growth has ceased, in particular for indications such as:

- degenerative diseases of the discs
- fractures
- tumors
- pseudarthrosis
- **spinal canal stenosis**
- cervical myelopathy
- deformities (i.e. cyphosis, lordosis and/or scoliosis)
- surgical revisions

AREAS OF APPLICATION

The area of application is the cervical spine from C2 to C7.

GENERAL CONDITIONS FOR USE

- The implants may only be implanted by surgeons who have successfully completed the necessary training in the field of spinal surgery. The decision to use implants must be made in accordance with the surgical and medical indications, taking into account the potential dangers and the limitations associated with this type of surgical intervention, as well in knowledge of the contraindications, side effects and the defined precautionary measures and in awareness of the properties and the physical, chemical and biological characteristics of the implant.
- It is recommended that the HERO® Cervical Plate System is not used in conjunction with implants from another supplier or another manufacturer. If this recommendation is not respected then HumanTech Germany GmbH can accept no liability.
- Under no circumstance may the implants be re-used. Even if the implant appears intact after the revision, there may be alterations within the implant or small defects, which are attributable to the action of loads and strain and which may result in fracture of the implant.

CONTRAINDICATIONS

The contraindications include, but are not limited to:

- All destructive and inflammatory diseases of the cervical vertebrae and of the motor segments such as tumors, spondylitis and fractures with a high degree of instability, with pronounced static and structural alteration of the spine and the segment structure, which require a vertebral body replacement and other measures to reconstruct and stabilize the segment.
- Acute infections or significant risks of infection (weakened immune system).
- Signs of local inflammations

- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Physical diseases
- Severe anatomical deformities caused by congenital abnormalities
- Any other medical or surgical condition, which prevents a possible improvement by the use of the implant, such as the presence of congenital abnormalities, a fracture close to the operation site, increase of the deposition rate which is not explained by other disease patterns, elevated white blood cell count (WBC) or a pronounced leftward shift in the differential WBC.
- Articular diseases, bone resorption, osteopenia, osteomalacia and/or osteoporosis are relative contraindications as these may limit the degree of achievable correction or stabilization.
- Known hereditary or acquired bone brittleness or calcification problems
- Suspected allergy or intolerance as well as documented allergy or intolerance of the material used.
- All cases in which the use of components made of various metals or alloys is required.
- All cases in which no fusion is required
- All cases in which the implant components selected for use are too big or too small to achieve a satisfactory result.
- Any patient with inadequate tissue structure at the surgical site or with inadequate bone stock or bone quality
- Any patient in whom the use of implant is hindered by anatomical structures or where these structures would restrict the physiological performance.
- Any patient who is unwilling to follow postoperative instructions.
- All cases not described under the indications.

SIDE EFFECTS

Complications may arise if the implant is implanted with or without the required instruments.

Side effects include, but are not limited to:

- Loosening of the implant
- Fracture of the implant
- Foreign body reactions to the implants including possible formation of tumors, autoimmune disease or scar formation.
- Neurological complications, paralysis, lesions of soft tissue and/or migration of the implant
- Neurological or spinal lesion of the dura mater due to surgical trauma
- Swallowing difficulties (dysphagia), hoarseness.
- Damage to the lymph vessels / lymph leakage.
- Changes to the curvature and rigidity of the spinal column
- Partial loss of degree of correction achieved in operation
- Superficial or internal infections and inflammations such as discitis, arachnoiditis
- Delayed engraftment of bone or no visible fusion and pseudarthrosis
- Fracture, microfracture, resorption, damage or penetration of a vertebra and/or bone transplant and/or bone transplant extraction site above or below the treated segment.
- Bleeding and/or hematomas
- Deep vein thrombosis, thrombophlebitis, lundenembolus
- Complications at bone transplant donor site
- Herniated disc, spinal disc destruction or degeneration, above or below the treated segment
- loss or increase of spinal mobility or functions
- Damage to the reproductive system, sterility and sexual dysfunction.
- Development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia etc.
- Changes in the mental condition
- Changed growth of the fused spinal column
- Allergic reaction to the Ti6Al4V alloy
- Appearance of microparticles in the vicinity of the implant
- Reduction of bone density
- Bursitis
- Death

PACKAGING, LABELING AND STORAGE

- The implants are NON STERILE. Before use they must be cleaned and sterilized (see below)
- The implants are delivered in packaged form. The packaging must be intact at the time of delivery.
- The implants are delivered as a complete implant set or as individual components.
- These implants and instruments are kept in specially designed trays/containers, which can be directly sterilized in the validated sterilization processes.
- The implant components must be handled and stored with due care. Damage to the implant can significantly reduce the robustness and

durability of the implant system. It can lead to cracks and/or increase internal loads, which can in turn result in a fracture of the implant.

- The implants and instruments should be stored at room temperature. Environmental influences such as salty air, humidity, chemicals etc. may not be allowed to act on the implants.
- It is to be recommended that the instruments and implants are carefully inspected before the operation to ensure that they have not been damaged during storage or preliminary procedures.

STERILIZATION

Implants of the HERO® Cervical Plate System must be sterilized or resterilized before use by means of validated steam sterilization processes in autoclaves. The correct process parameters in relation to sterilization temperature and cycle durations should be obtained from the appropriate instructions from the manufacturer. Sterilization should be performed by steaming according to the specifications of EN ISO 17665-1. The use of flash sterilization is not recommended. The minimum sterilization parameters are as follows:

Steam	Pre-vacuum	132	≥ 5 minutes
Steam	Pre-vacuum	121	20 minutes
Steam	Gravitation	121	30 minutes

Non-recommended processes are possible; however they must be completely validated by the user.

IMPORTANT NOTICE FOR THE ATTENDING SURGEON

- Implants of the HERO® Cervical Plate System have a restricted useful life. The activities / movement patterns of the patient have a significant impact on the useful life. The patient must be informed that every activity increases the risk of a loss, deformation or fracture of the implant components. It is of decisive importance to inform the patient about the restrictions to his or her activities in the post-operative phase and to monitor the patient in this phase so that the development of the fusion and the condition of the implant can be assessed. Even if a solid bone fusion occurs, the implant components may still bend, break or become loose. Therefore the patient must be informed that implant components may bend, fracture or become loose, even if the restrictions in activities are observed. If an implant breakage occurs, the physician must decide whether a revision of the implant should be performed, taking into account the condition of the patient and the possible associated risks.
- Informing the patient of how he or she can reduce loading on the implants is important so that the risk of re-occurrence of the clinical problem due to a failed fixing can be reduced.
- The instructions in the operation instructions (Surgical Technique) are to be observed. You can download these from the internet from www.humantech-solutions.de or directly from your HumanTech representative.
- Bending of the plates can have an adverse effect on the biomechanical properties of the implant. Deformations in the vicinity of locating holes can adversely affect the anchoring of the system. Any bending in this area is to be avoided.
- Proceed with extreme care in the vicinity of the spinal cord and the nerve roots, as any damage to nerves can result in the loss of neurological functions.
- Fracture, slipping or incorrect use of instruments or implants can result in injury to the patient or surgical personnel.
- Bone cement may not be used as it makes the removal of components difficult or even impossible.
- Once removed implants are to be treated so as to render any re-use impossible.

PRODUCT COMPLAINTS

Any person involved in healthcare (e.g. customer or user of this product system) who has complaints of any kind or is dissatisfied using the product with respect to quality, identity, strength, durability, safety, effectiveness or function should notify the appropriate HumanTech representative. If a HERO® implant should ever malfunction (i.e. fails to fulfill the performance specifications or does not function in the foreseen manner), or if the above is suspected, the HumanTech representative should be informed immediately. If a HumanTech product should ever malfunction in a way that results in the death or serious injury of the patient or if it has contributed to the above, the representative is to be informed immediately by telephone, fax or in written form. If you have any complaints then we ask you to state the name, article number and batch number of the components as well as your name and address along with a detailed description of the fault and forward this information to us in writing.

FURTHER INFORMATION

If you have any complaints, suggestions or other points about the contents of this user manual or the use of the product, please contact the above address.

© HUMANTECH Germany GmbH. All rights reserved.

	Manufacturer
	For single use only
	Lot number
	Order number
	Read the documentation
	Protect from moisture
	Non sterile
	Do not use if packaging has been damaged

 	
HERO® CERVICAL PLATE SYSTEM	
HumanTech Germany GmbH Gewerbestraße 5 D-71144 Steinenbronn Tel.: +49 (0) 7157/5246-71 Fax.: +49 (0) 7157/5246-33 e-mail: info@humantech-solutions.com www.humantech-solutions.com	 

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© HUMANTECH Germany GmbH. Tüm hakları saklıdır.